
Heraeus

PALACOS[®] MV

Medium viscosity, radiopaque bone cement

PALACOS® MV

EN	Instructions for use Working times (Revision status: 2022-02)	5 98
AR	تعليمات الاستخدام أوقات العمل (تاريخ المراجعة 2022-02)	8 98
BG	Инструкции за употреба Време за обработка при (Дата на актуализиране на текста: 2022-02)	11 98
CS	Návod k použití Doby zpracování (Poslední revise: 2022-02)	14 98
DA	Brugsanvisning Arbejdstider (Aktualiseret: 2022-02)	17 98
DE	Gebrauchsanweisung Verarbeitungszeiten (Revisionsstand: 2022-02)	20 98
EL	Οδηγίες χρήσης Χρόνοι εργασίας (Ημερομηνία αναθεώρησης: 2022-02)	23 98
ES	Instrucciones de uso Tiempos de manipulación (Última revisión: 2022-02)	26 98
ET	Kasutusjuhend Tööjad (Viimati muudetud: 2022-02)	29 98
FI	Käyttöohje Työstöajat (Revisiostatus: 2022-02)	32 98
FR	Notice d'utilisation Temps de travail (Date de révision: 2022-02)	35 98
HR	Upute za uporabu Vremena obrade (Stanje revizije: 2022-02)	38 98
HU	Használati utasítás Megmunkálási (Aktualizálás időpontja: 2022-02)	41 98

IT	Istruzioni per l'uso	44
	Tempi di lavorazione	98
	(Ultima revisione: 2022-02)	
KO	사용 설명서	47
	작업 시간	98
	(2022년 02월 수정)	
LT	Naudojimo instrukcija	50
	Darbo laikas	98
	(Tikrinta: 2022-02)	
LV	Lietošanas instrukcija	53
	Apstrādes laiks	98
	(Pārbaudes statuss: 2022-02)	
NL	Gebruiksaanwijzing	56
	Verwerkingstijden	98
	(Revisiestatus: 2022-02)	
NO	Brugsanvisning	59
	Virketider	98
	(Revisjonsstatus: 2022-02)	
PL	Instrukcja użytkownika	62
	Czasy obróbki	98
	(Ostatnia aktualizacja: 2022-02)	
PT	Instruções de utilização	65
	Fases de aplicação	98
	(Data da última actualização: 2022-02)	
PT (BR)	Nota: Apenas para fins brasileiros!	
	Instruções de utilização	68
	Fases de aplicação	98
	(Data da última actualização: 2022-02)	
RO	Instrucțiuni de utilizare	71
	Timpi de prelucrare	98
	(Stadiu revizie: 2022-02)	
RU	Инструкция по применению	74
	Временные фазы для смешивания	98
	(Ревизионный статус: 2022-02)	
SK	Návod na použitie	77
	Časy spracovateľnosti	98
	(Posledná revízia textu: 2022-02)	
SL	Navodilo za uporabo	80
	Časovne opredelitve	98
	(Datum revizije: 2022-02)	

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Properties

PALACOS® MV is a fast-curing, radiopaque, poly-(methyl methacrylate)-based bone cement.

PALACOS® MV contains the X-ray contrast medium zirconium dioxide. To improve visibility in the surgical field

PALACOS® MV has been coloured with chlorophyll (E141).

The bone cement is prepared directly before use by mixing a polymer powder component with a liquid monomer component. A ductile dough forms which cures within a few minutes.

Composition

A pack of **PALACOS® MV** contains one/or two bag/s of cement powder (polymer powder) and one/two amber glass ampoule/s (monomer liquid).

Polymer powder ingredients:

Poly-(methylacrylate, methyl methacrylate), zirconium dioxide, benzoyl peroxide and colorant E141.

Monomer liquid ingredients:

Methyl methacrylate, N,N-dimethyl-p-toluidine, hydroquinone and colorant E141.

The cement powder is triple packaged. The outer, unsterile protective aluminium bag contains a polyethylene paper bag (peel-off) which is unsterile on the outside and sterile on the inside. Inside this there is another sterile polyethylene paper bag containing the cement powder. The brown glass ampoule with the sterile-filtered monomer liquid is also sterile packed in an ethylene oxide sterilised individual blister pack.

Indication

PALACOS® MV is indicated for

- Stable attachment of total or partial joint endoprostheses in bone
- Filling and stabilising bone defects within the scope of internal fixation treatment or for endoprosthesis revision surgery
- Primary and secondary coverage of skull bone defects

Contraindication

PALACOS® MV must not be used

- If it is known or supposed that the patient is hypersensitive to ingredients of the bone cement
- During pregnancy or breast-feeding

Target population

There is very little evidence with children and adolescents, therefore it is not recommended to use **PALACOS® MV**. If no other option is available, for example with present surgical trauma, the decision whether to use **PALACOS® MV** lies with the attending surgeon.

Target user group

Healthcare professionals in a clinical context.

Side effects

Rarely there is a temporary drop in blood pressure after preparation of the prosthesis bed or directly after the implantation of PMMA bone cements and the endoprosthesis.

In single cases serious complications can occur such as severe allergic reactions which are associated with cardiac arrest, anaphylactic shock or even sudden death.

In order to avoid pulmocardiovascular complications such as pulmonary embolism and cardiac arrest it is recommended that the implantation site be irrigated thoroughly with an isotonic solution (pulse lavage application) before the bone cement is introduced.

In the event of pulmonary or cardiovascular events it is necessary to monitor blood volume and possibly increase it.

In the case of acute respiratory insufficiency anaesthesiological measures should be taken.

The following adverse effects have been observed when using poly(methyl methacrylate) bone cements: thrombophlebitis, haemorrhage, trochanteric bursitis.

Other side effects observed: myocardial infarction, brief cardiac arrhythmia, cerebrovascular incident. Moreover, complications are possible that can occur during any surgical procedure.

Precautions**Use by operating staff**

Before using **PALACOS® MV** the user should be well-acquainted with its properties, handling and application.

The user is advised to practise the entire procedure of mixing, handling and introducing **PALACOS® MV** before using it for the first time. Detailed knowledge is necessary even if mixing systems and syringes are being used for application of the cement.

The liquid monomer is highly volatile and flammable. Ignition of monomer fumes caused by use of electrocautery devices in surgical sites near freshly implanted bone cements has been reported. The monomer is also a strong lipid solvent and should not come into direct contact with the body. When handling the monomer or the prepared cement **PALACOS® MV** gloves must be worn which provide the necessary protection against penetration of the skin by the monomer methyl methacrylate.

Gloves made of PVP (three-ply polyethylene, ethylene vinyl alcohol copolymer, polyethylene) and Viton®/butyl gloves have proved to provide good protection over a lengthy period. For safety's sake it is recommended that two pairs of gloves be worn over one another, e.g. one polyethylene surgical glove over an inner pair of standard latex surgical gloves.

The use of latex or polystyrene-butadiene gloves on their own is inadequate. Please make enquiries with your supplier to establish which gloves are suitable for such an application.

The monomer vapours can irritate respiratory passages and eyes and possibly cause damage to the liver.

Skin irritations have been described which are due to contact with the monomer.

Manufacturers of soft contact lenses recommend removing the lenses in the presence of damaging or irritant vapours. Since soft contact lenses are permeable to liquids and gases, they should not be worn in the operating theatre if methyl methacrylate is being used.

Use on the patient

Blood pressure, pulse and breathing must be monitored carefully during and directly after introduction of the bone cement. Any significant change in these vital signs must be eliminated without delay by taking appropriate action. When using **PALACOS® MV** the prepared bone should be carefully cleaned, aspirated and dried just before the bone cement is placed.

Incompatibilities

Aqueous solutions (e.g. ones containing antibiotics) must not be added to the bone cement because they have a considerable detrimental effect on the physical and mechanical properties of the cement.

Amount required

After mixing the cement powder with the monomer liquid, a fast-curing ductile dough develops which is introduced to bone cavities for anchorage and/or filling purposes.

PALACOS® MV is coloured green in order to make the cement clearly visible in the surgical field.

A dose is prepared by mixing the entire contents of the bag of cement powder with all the monomer liquid of an ampoule. The quantity of cement dough required depends on the specific surgical intervention and on the technique being used.

At least one additional dose of **PALACOS® MV** should be available before commencing the operation.

Preparation:

Preparation:

Before opening the non-sterile aluminium protective pouch (compare above) move the contents down (1 polyethylene-paper bag) by shaking or tapping in order to ensure that when the bag is cut open at the top the content is not damaged. The polyethylene-paper bag and the ampoule may only be opened under sterile conditions. For this purpose the sterile components (inner PE-paper bag and glass ampoule) are sterile when they are delivered.

Opening under sterile conditions:

The opening flaps from the bag top assist to detach the foil from the paper.

In order to grip as much of the opening flaps as possible, the side of the paper/foil should be kept between thumb, index and middle finger.

Please use the whole thumb surface to grip foil and paper side and take off each side evenly.

Open the outer PE-paper bag at the special point under sterile conditions so that the inner PE-paper bag remains sterile when it is removed. Also open the blister pack at the special point under sterile conditions so that the glass ampoule remains sterile when it is removed.

Before opening the inner PE-paper bag, moved the contents down by shaking or tapping in order to ensure that when the bag is cut open at the top no powder is lost. To make it easier to open the glass ampoule the latter is provided with a pre-determined breaking point at the transition to the head of the ampoule.

Do not open the ampoule over the mixing device to prevent contamination of the cement with glass fragments.

The 10ml- and 20ml- ampoule is provided with a snapping-off device (tube) in order to facilitate the opening procedure. If this is the case, take hold of the fitted snapping-off device instead of the ampoule head and break off the ampoule head over it. When the ampoule head has snapped off, it remains inside the tube.

Mixing the components:

For mixing of **PALACOS® MV** all Heraeus Medical mixing systems are suitable. Should the user want to use

PALACOS® MV in combination with other mixing systems, he is advised to check their suitability prior to first use.

It is advisable to first measure out the liquid and then add the powder. If this order is reversed, powder nests are more likely to form as a result of polymerisation commencing immediately at the surface. Both components, i.e. the relative proportions of powder and monomer, are precisely matched. The pouch and ampoule must therefore be emptied completely if an optimal mix is to be achieved.

The components can be mixed with a vacuum mixing system or by hand.

The mixing, waiting, working and curing times of **PALACOS® MV** are shown in the diagrams at the end of the instructions for use. Please note that these are stated for guidance only because working time and curing time depend on temperature, mixing and humidity, whereby the direct ambient temperatures are important, e.g. of cement powder, mixing system, bench and hands. A high temperature accelerates the waiting, working and curing times.

Preparation with a vacuum mixing system

In order to obtain a bone cement with a reduced number of air inclusions, the liquid and the powder are mixed in a vacuum. For this purpose an airtight mixing system has to be used, ensuring that a sufficient vacuum builds up in the mixing vessel quickly (pressure approx. 200 mbar absolute). Perform filling and mixing under sterile conditions. The mixing time is 30 seconds unless otherwise recommended. For details of mixing technique please refer to the instructions for the mixing system being used. The result is a homogeneous, green, doughy compound which can be worked when it no longer adheres to the rubber gloves. Always mix the entire contents of a bag with the entire contents of an ampoule of monomer liquid.

Preparation mixing by hand

The cement components should only be filled into the mixing vessel just before mixing. Perform filling and mixing under sterile conditions – mixing time is 30 seconds. During this time the two components are mixed with one another while stirring evenly. The result is a homogeneous, green, doughy compound which can be worked when it no longer adheres to the rubber gloves. Always mix the entire contents of a bag with the entire contents of an ampoule of monomer liquid.

Use of the bone cement

To ensure adequate fixation the prosthesis must be introduced and held within the time window allowed for working, until the bone cement has set completely. Remove any surplus cement as long as it is still soft.

Cranial vault reconstruction

In the management of major skull defects, after careful preparation of the bone gap the dura mater is first of all covered with moist cotton wool or cellulose.

A thin plastic or aluminium foil is laid over it for further protection. The mixed pasty dough is introduced to the prepared bone gap and moulded at the bone margins to the required thickness of 4–5 mm. While it is setting, irrigate with normal saline solution in order to dissipate the heat caused by polymerisation. When the reconstruction has almost hardened it is then removed, corrected at the margins and provided with perforations, through which epidural fluid can drain away and connective tissue can grow. After removal of the cotton wool/cellulose and plastic foil the prosthesis is fixed in place at three or four points with non-absorbable suture.

Storage

Do not store above 25°C (77°F).

Shelf life/sterility

The shelf life is printed on the folding box, the protective aluminium bag and the inner bag.

Do not use **PALACOS® MV** if this date indicated has expired.

The information on the blister pack may be different for production reasons.

The contents of opened or damaged aluminium sachets or ampoule blisters must not be re-sterilized and thus must be discarded. If the cement powder has turned yellow, do not use **PALACOS® MV**. The cement powder and packagings are sterilised by gassing with ethylene oxide. The monomer liquid has been sterilised by filtration.

Disposal

Single components of the bone cement, cured solid material as well as (uncleaned) packaging material must be disposed by following the regulations of the local authorities.

Dispose of the polymer component in an authorized waste facility. The liquid component should be evaporated under a well-ventilated hood or absorbed by an inert material and transferred in a suitable container for disposal.

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Heraeus Medical GmbH