
Heraeus

# PALACOS<sup>®</sup> LV+G

Low viscosity, radiopaque bone cement  
containing Gentamicin

# PALACOS® LV+G

EN	<b>Instructions for use</b> <b>Working times</b> (Revision status: 2019-07)	<b>5</b> <b>136</b>
EN (CA)	<b>Note: Applicable in Canada!</b> <b>Instructions for use</b> <b>Working times</b> (Revision status: 2019-11)	<b>9</b> <b>136</b>
AR	<b>تعليمات الاستخدام</b> <b>أوقات العمل</b> (تاريخ المراجعة: 2019-07)	<b>13</b> <b>136</b>
BG	<b>Инструкции за употреба</b> <b>Време за обработка при</b> (Дата на актуализиране на текста: 2019-07)	<b>17</b> <b>136</b>
CS	<b>Návod k použití</b> <b>Doby zpracování</b> (Poslední revise: 2019-07)	<b>21</b> <b>136</b>
DA	<b>Brugsanvisning</b> <b>Arbejdstider</b> (Aktualiseret: 2019-07)	<b>25</b> <b>136</b>
DE	<b>Gebrauchsanweisung</b> <b>Verarbeitungszeiten</b> (Revisionsstand: 2019-07)	<b>29</b> <b>136</b>
EL	<b>Οδηγίες χρήσης</b> <b>Χρόνοι εργασίας</b> (Ημερομηνία αναθεώρησης: 2019-07)	<b>33</b> <b>136</b>
ES	<b>Instrucciones de uso</b> <b>Tiempos de manipulación</b> (Última revisión: 2019-07)	<b>37</b> <b>136</b>
ET	<b>Kasutusjuhend</b> <b>Tööjad</b> (Viimati muudetud: 2019-07)	<b>41</b> <b>136</b>
FI	<b>Käyttöohje</b> <b>Työstöajat</b> (Revisiostatus: 2019-07)	<b>45</b> <b>136</b>
FR	<b>Notice d'utilisation</b> <b>Temps de travail</b> (Date de révision: 2019-07)	<b>49</b> <b>136</b>
FR (CA)	<b>Nota: Applicable au Canada!</b> <b>Notice d'utilisation</b> <b>Temps de travail</b> (Date de révision: 2019-11)	<b>53</b> <b>136</b>

HR	<b>Upute za uporabu Vremena obrade</b> (Stanje revizije: 2019-07)	<b>57</b> <b>136</b>
HU	<b>Használati utasítás Mégmunkálási</b> (Aktualizálás időpontja: 2019-07)	<b>61</b> <b>136</b>
IT	<b>Istruzioni per l'uso Tempi di lavorazione</b> (Ultima revisione: 2019-07)	<b>65</b> <b>136</b>
KO	<b>사용 설명서 작업 시간</b> (2019년 07월 수정)	<b>69</b> <b>136</b>
LT	<b>Naudojimo instrukcija Darbo laikas</b> (Tikrinta: 2019-07)	<b>73</b> <b>136</b>
LV	<b>Lietošanas instrukcija Apstrādes laiks</b> (Pārbaudes statuss: 2019-07)	<b>77</b> <b>136</b>
NL	<b>Gebruiksaanwijzing Verwerkingstijden</b> (Revisiestatus: 2019-07)	<b>81</b> <b>136</b>
NO	<b>Brugsanvisning Virketider</b> (Revisjonsstatus: 2019-07)	<b>85</b> <b>136</b>
PL	<b>Instrukcja użytkowania Czasy obróbki</b> (Ostatnia aktualizacja: 2019-07)	<b>89</b> <b>136</b>
PT	<b>Instruções de utilização Fases de aplicação</b> (Data da última actualização: 2019-07)	<b>93</b> <b>136</b>
PT (BR)	<b>Nota: Apenas para fins brasileiros! Instruções de utilização Fases de aplicação</b> (Data da última actualização: 2019-07)	<b>97</b> <b>136</b>
RO	<b>Instrucțiuni de utilizare Timpi de prelucrare</b> (Stadiu revizie: 2019-07)	<b>100</b> <b>136</b>
RU	<b>Инструкция по применению Временные фазы для смешивания</b> (Ревизионный статус: 2019-07)	<b>104</b> <b>136</b>

SK	<b>Návod na použitie</b>	<b>108</b>
	<b>Časy spracovateľnosti</b> (Posledná revízia textu: 2019-07)	<b>136</b>
SL	<b>Navodilo za uporabo</b>	<b>112</b>
	<b>Časovne opredelitve</b> (Datum revizije: 2019-07)	<b>136</b>
SV	<b>Bruksanvisning</b>	<b>116</b>
	<b>Hanteringstider</b> (Revision status: 2019-07)	<b>136</b>
TH	<b>คำแนะนำวิธีใช้</b>	<b>120</b>
	<b>เวลาในการใช้งาน</b> (แก้ไขปรับปรุงครั้งสุดท้ายเมื่อ: 2019-07)	<b>136</b>
TR	<b>Kullanma Talimatı</b>	<b>124</b>
	<b>İşleme süreleri</b> (Bilgilerin son güncellenme tarihi: 2019-07)	<b>136</b>
UK	<b>Інструкція для застосування</b>	<b>128</b>
	<b>Робочий час</b> (Редакція: 2019-07)	<b>136</b>
ZH	<b>使用说明书</b>	<b>第132页</b>
	<b>加工时间</b> (修订状态: 2019 年 07 月)	<b>第136页</b>

**Properties**

**PALACOS® LV+G** is a slow-curing, radiopaque, poly(methyl methacrylate)-based bone cement. It contains the aminoglycoside antibiotic gentamicin to protect the cured cement and contiguous tissue against contamination by microbes sensitive to gentamicin.

**PALACOS® LV+G** contains the X-ray contrast medium zirconium dioxide. To improve visibility in the surgical field **PALACOS® LV+G** has been coloured with chlorophyll (E141). The bone cement is prepared directly before use by mixing a polymer powder component with a liquid monomer component. A ductile dough forms which cures within a few minutes.

**Composition**

A pack of **PALACOS® LV+G** contains one or two bag/s of gentamicin-containing cement powder (polymer powder) and one/two amber glass/es ampoule (monomer liquid).

Cement powder ingredients:

Gentamicin (in the form of gentamicin sulphate), poly(methylacrylate, methyl methacrylate), zirconium dioxide, benzoyl peroxide and colorant E141.

**Ingredients**

Pack size	40
Gentamicin base	1.0g
Cement powder	41.7 g
Monomer liquid	20ml

The mass ratio is about 31 percent monomer liquid to 69 percent cement powder.

Monomer liquid ingredients:

Methyl methacrylate, N,N-dimethyl-p-toluidine, hydroquinone and colorant E141.

The cement powder is triple packaged. The outer, unsterile protective aluminium bag contains a polyethylene paper bag (peel-off) which is unsterile on the outside and sterile on the inside. Inside this there is another sterile polyethylene paper bag containing the cement powder. The brown glass ampoule with the sterile-filtered monomer liquid is also sterile packed in an ethylene oxide sterilised individual blister pack.

**Composition****The cement powder contains:**

poly(methyl acrylate, methyl methacrylate)	80 %
zirconium dioxide	15 %
hydrous benzoyl peroxide	1 %
gentamicin sulphate	4 %

**The monomer liquid contains:**

methyl methacrylate	98 %
N,N-dimethyl-p-toluidine	2 %

**Other constituents:**

In the cement powder: chlorophyll VIII (colorant E141)  
In the monomer liquid: chlorophyll VIII (colorant E141)  
in an oily solution, hydroquinone

**Indication**

The use of **PALACOS® LV+G** with gentamicin is indicated for

- stable anchoring of total or partial joint endoprostheses in bone or filling and stabilising bone defects within the scope of internal fixation treatment or for endoprostheses revision surgery
- primary and secondary plastic reconstruction of bone defects.

**Contraindication**

**PALACOS® LV+G** with gentamicin must not be used in the following cases

- In case of suspected or proven hypersensitivity to components of the bone cement
- During pregnancy or breast-feeding
- In cases of severe renal insufficiency
- In the presence of an active or incompletely treated infection at the bone site caused by gentamicin non-sensitive strains
- **PALACOS® LV+G** is not suitable for primary or secondary plastic reconstruction of skull bone defects.

**Target population**

There is very little evidence with children and adolescents, therefore it is not recommended to use **PALACOS® LV+G**.

If no other option is available, for example with present surgical trauma, the decision whether to use **PALACOS® LV+G** lies with the attending surgeon.

**Target user group**

Healthcare professionals in a clinical context.

## Side effects

### General

Rarely there is a temporary drop in blood pressure after preparation of the prosthesis bed or directly after the implantation of bone cements and the endoprosthesis. In individual cases severe complications can occur such as severe allergic reactions which are associated with cardiac arrest, anaphylactic shock or even sudden death. In order to avoid pulmocardiocardiovascular complications such as pulmonary embolism and cardiac arrest it is recommended that the implantation site be irrigated thoroughly with an isotonic solution (pulse lavage application) before the bone cement is introduced. In the event of pulmonary or cardiovascular events it is necessary to monitor blood volume and possibly increase it. In the case of acute respiratory insufficiency anaesthesiological measures should be taken. The following adverse effects have been observed when using poly(methyl methacrylate) bone cements: thrombophlebitis, haemorrhage, trochanteric bursitis. Other side effects that were observed are heterotopic bone formation, osteolysis due to cement fragments, myocardial infarction, brief cardiac arrhythmia, cerebrovascular incident.

### Gentamicin

Due to the gentamicin contained in **PALACOS® LV+G**, side effects that are typical for gentamicin can occur:

- impairment to auditory and vestibular nerves
- renal toxicity: Adverse renal effects have been reported, as demonstrated by the presence of cells or protein in the urine or by rising of serum creatinine or oliguria. They occur more frequently in patients with a history of renal impairment. Application of clinical doses of gentamicin has led to the occasional appearance of Fanconi's syndrome or a Bartter's-like syndrome
- neuromuscular blockade/ neurotoxicity: Serious adverse effects on both vestibular and auditory branches of the 8th cranial nerves have been reported, primarily in patients with renal impairment. Symptoms include dizziness, vertigo, tinnitus, roaring in the ears and hearing loss, which, as with the other aminoglycosides, may be irreversible. Other factors which may increase the risk of toxicity include dehydration and previous exposure to other ototoxic drugs
- rare cases of paresthesia, tetany and amyasthenia
- rare cases of allergic reactions (exanthema, urticaria, anaphylactic reactions).

Like all aminoglycosides gentamicin is also potentially nephrotoxic and/or ototoxic. In most cases, damage to hearing caused by gentamicin is irreversible, whereas renal damage is reversible.

Other reported adverse reactions possibly related to gentamicin include: respiratory depression, lethargy, confusion, depression, visual disturbances, decreased appetite, weight loss, and hypotension or hypertension; rash, itching, urticaria, generalised burning, laryngeal oedema, anaphylactoid reactions, fever, and headache, nausea, vomiting, increased salivation, and stomatitis; purpura, pseudotumor cerebri, acute organic brain syndrome, pulmonary fibrosis, alopecia, joint pain, transient hepatomegaly, and splenomegaly.

The use of **PALACOS® LV+G**, results in a desired high bacte-

ricidal but locally limited concentration of gentamicin. Therefore, the occurrence of the afore-mentioned side effects is highly unlikely. Care should be taken and the level of gentamicin should be monitored in patients with impaired renal function, especially if treated with parenteral aminoglycosides or muscle relaxants simultaneously. The same applies to patients with previous neuromuscular disorders (e.g. myasthenia gravis, Parkinson's disease). Allergic reactions can occur dosage-independent.

## Interactions

### Gentamicin

Owing to the administration of muscle relaxants and ether, the neuromuscular blocking properties of gentamicin may be intensified. However, in view of the very low serum levels this is relatively unlikely in patients with healthy kidneys. The probability of interactions occurring increases in proportion to the serum levels of gentamicin, e.g. in patients with impaired renal function.

Interactions of gentamicin have been reported to occur with the following groups of drugs:

Concurrent or sequential use of neurotoxic and / or nephrotoxic pharmaceuticals, including other aminoglycosides, polymyxin B, colistin, cisplatin, ciclosporines, foscarnet, vancomycin, amphotericin B, clindamycin and cephalosporins.

Potent diuretics: Ethacrynic acid, furosemide, or other potent diuretics which may themselves cause ototoxicity or enhance aminoglycoside toxicity by altering antibiotic concentrations in serum and tissue.

Muscle relaxants: in particular succinyl choline and tubocurarine, decamethonium, halogenated hydro-carbon inhalation anaesthetics or opioid analgesics. Gentamicin may interact with these medications and result in skeletal muscle weakness and respiratory depression (apnea).

Concurrent application of these medications and gentamicin during surgery or in the post-operative period should be monitored with caution, especially if there is a possibility of incomplete reversal of neuromuscular blockade post-operatively.

Moreover, complications are possible that can occur during any surgical procedure.

## Precautions

### Use by operating staff

Before using **PALACOS® LV+G** the user should be well-acquainted with its properties, handling and application. The user is advised to practise the entire procedure of mixing, handling and introducing **PALACOS® LV+G** before using it for the first time. Detailed knowledge is necessary even if mixing systems and syringes are being used for application of the cement.

The liquid monomer is highly volatile and flammable. Ignition of monomer fumes caused by use of electrocautery devices in surgical sites near freshly implanted bone cements has been reported. The monomer is also a strong lipid solvent and should not come into direct contact with the body.

When handling the monomer or the prepared cement **PALACOS® LV+G** gloves must be worn which provide the necessary protection against penetration of the skin by the monomer methyl methacrylate.

Gloves made of PVP (three-ply polyethylene, ethylene vinyl alcohol copolymer, polyethylene) and Viton®/butyl gloves have proved to provide good protection over a lengthy period. For safety's sake it is recommended that two pairs of gloves be worn over one another, e.g. one polyethylene surgical glove over an inner pair of standard latex surgical gloves.

The use of latex or polystyrene-butadiene gloves on their own is inadequate. Please make enquiries with your supplier to establish which gloves are suitable for such an application. The monomer vapours can irritate respiratory passages and eyes and possibly cause damage to the liver. Skin irritations have been described which are due to contact with the monomer.

Manufacturers of soft contact lenses recommend removing the lenses in the presence of damaging or irritant vapours. Since soft contact lenses are permeable to liquids and gases, they should not be worn in the operating theatre if methyl methacrylate is being used.

#### Use on the patient

Blood pressure, pulse and breathing must be monitored carefully during and directly after introduction of the bone cement. Any significant change in these vital signs must be eliminated without delay by taking appropriate action.

When using **PALACOS® LV+G** the prepared bone should be carefully cleaned, aspirated and dried just before the bone cement is placed.

#### **Incompatibilities**

Aqueous solutions (e.g. ones containing antibiotics) must not be added to the bone cement because they have a considerable detrimental effect on the physical and mechanical properties of the cement.

#### **Amount required**

After mixing the cement powder with the monomer liquid, a slow-curing ductile dough develops which is introduced to bone cavities for anchorage and/or filling purposes.

**PALACOS® LV+G** is coloured green in order to make the cement clearly visible in the surgical field.

A dose is prepared by mixing the entire contents of the bag of cement powder with all the monomer liquid of an ampoule. The quantity of cement dough required depends on the specific surgical intervention and on the technique being used.

At least one additional dose of **PALACOS® LV+G** should be available before commencing the operation.


#### **Preparation**

##### Preparation:

Before opening the non-sterile aluminium protective pouch (compare above) move the contents down (1 polyethylene-paper bag) by shaking or tapping in order to ensure that when the bag is cut open at the top the content is not damaged.

The polyethylene-paper bag and the ampoule may only be opened under sterile conditions. For this purpose the sterile components (inner PE-paper bag and glass ampoule) are sterile when they are delivered.

#### Opening under sterile conditions:

The opening flaps from the bag top assist to detach the foil from the paper.

In order to grip as much of the opening flaps as possible, the side of the paper/foil should be kept between thumb, index and middle finger.

Please use the whole thumb surface to grip foil and paper side and take off each side evenly.

Open the outer PE-paper bag at the special point under sterile conditions so that the inner PE-paper bag remains sterile when it is removed. Also open the blister pack at the special point under sterile conditions so that the glass ampoule remains sterile when it is removed.

Before opening the inner PE-paper bag, moved the contents down by shaking or tapping in order to ensure that when the bag is cut open at the top no powder is lost. To make it easier to open the glass ampoule the latter is provided with a predetermined breaking point at the transition to the head of the ampoule.

Do not open the ampoule over the mixing device to prevent contamination of the cement with glass fragments.

The 10 ml- and 20 ml- ampoule is provided with a snapping-off device (tube) in order to facilitate the opening procedure. If this is the case, take hold of the fitted snapping-off device instead of the ampoule head and break off the ampoule head over it. When the ampoule head has snapped off, it remains inside the tube.

#### Mixing the components:

It is advisable to first measure out the liquid and then add the powder. If this order is reversed, powder nests are more likely to form as a result of polymerisation commencing immediately at the surface. Both components, i.e. the relative proportions of powder and monomer, are precisely matched. The pouch and ampoule must therefore be emptied completely if an optimal mix is to be achieved.

The components can be mixed with a vacuum mixing system or by hand.

The mixing, waiting, working and curing times of **PALACOS® LV+G** are shown in the diagrams at the end of the instructions for use. Please note that these are stated for guidance only because working time and curing time depend on temperature, mixing and humidity, whereby the direct ambient temperatures are important, e.g. of cement powder, mixing system, bench and hands. A high temperature accelerates the waiting, working and curing times.

#### **Preparation with a vacuum mixing system**

In order to obtain a bone cement with a reduced number of air inclusions, the liquid and the powder are mixed in a vacuum. For this purpose an airtight mixing system has to be used, ensuring that a sufficient vacuum builds up in the mixing vessel quickly (pressure approx. 200 mbar absolute).

Perform filling and mixing under sterile conditions.

The mixing time is 30 seconds unless otherwise recommended. For details of mixing technique please refer to the instructions for the mixing system being used.

Always mix the entire contents of a bag with the entire contents of an ampoule of monomer liquid.

**Preparation mixing by hand**

The cement components should only be filled into the mixing vessel just before mixing. Perform filling and mixing under sterile conditions – mixing time is 30 seconds. During this time the two components are mixed with one another while stirring evenly. Always mix the entire contents of a bag with the entire contents of an ampoule of monomer liquid.

**Use of the bone cement**

The bone cement can be applied as soon as the doughy bone cement no longer adheres to the gloves. The application period depends on the temperature of the material and the room temperature. To ensure adequate fixation the prosthesis must be introduced and held within the time window allowed for application, until the bone cement has set completely. Remove any surplus cement as long as it is still soft.

**Storage**

Do not store above 25 °C (77 °F).

**Shelf life/sterility**

The shelf life is printed on the folding box, the protective aluminium bag and the inner bag. Do not use **PALACOS® LV+G** if this date indicated has expired. The information on the blister pack may be different for production reasons. The contents of opened or damaged aluminium sachets or ampoule blisters must not be re-sterilized and thus must be discarded. If the cement powder has turned yellow, do not use **PALACOS® LV+G**. The cement powder and packagings are sterilised by gassing with ethylene oxide. The monomer liquid has been sterilised by filtration.

**Disposal**

Single components of the bone cement, cured solid material as well as (uncleaned) packaging material must be disposed by following the regulations of the local authorities. Dispose of the polymer component in an authorized waste facility. The liquid component should be evaporated under a well-ventilated hood or absorbed by an inert material and transferred in a suitable container for disposal.

EN (CA)

# PALACOS® LV+G

**Note: Applicable in Canada!**

## Properties

**PALACOS® LV+G** is a slow-curing, radiopaque, poly(methyl methacrylate)-based bone cement. It contains the aminoglycoside antibiotic gentamicin to protect the cured cement and contiguous tissue against contamination by microbes sensitive to gentamicin.

**PALACOS® LV+G** contains the X-ray contrast medium zirconium dioxide. To improve visibility in the surgical field **PALACOS® LV+G** has been coloured with chlorophyll (E141). The bone cement is prepared directly before use by mixing a polymer powder component with a liquid monomer component. A ductile dough forms which cures within a few minutes.

## Composition

A pack of **PALACOS® LV+G** contains one or two bag/s of gentamicin-containing cement powder (polymer powder) and one/two amber glass/es ampoule (monomer liquid).

### Cement powder ingredients:

Gentamicin (in the form of gentamicin sulphate), poly(methylacrylate, methyl methacrylate), zirconium dioxide, benzoyl peroxide and colorant E141.

## Ingredients

Pack size	40
Gentamicin base	1.0g
Cement powder	41.7 g
Monomer liquid	20ml

The mass ratio is about 31 percent monomer liquid to 69 percent cement powder.

### Monomer liquid ingredients:

Methyl methacrylate, N,N-dimethyl-p-toluidine, hydroquinone and colorant E141.

The cement powder is triple packaged. The outer, sterile protective aluminium bag contains a polyethylene paper bag (peel-off) which is unsterile on the outside and sterile on the inside. Inside this there is another sterile polyethylene paper bag containing the cement powder. The brown glass ampoule with the sterile-filtered monomer liquid is also sterile packed in an ethylene oxide sterilised individual blister pack.

## Composition

### The cement powder contains:

poly(methyl acrylate, methyl methacrylate)	80 %
zirconium dioxide	15 %
hydrous benzoyl peroxide	1 %
gentamicin sulphate	4 %

### The monomer liquid contains:

methyl methacrylate	98 %
N,N-dimethyl-p-toluidine	2 %

## Other constituents:

In the cement powder: chlorophyll VIII (colorant E141)  
In the monomer liquid: chlorophyll VIII (colorant E141) in an oily solution, hydroquinone

## Indication

The use of **PALACOS® LV+G** with gentamicin is indicated for

- stable anchoring of total or partial joint endoprostheses in bone or filling and stabilising bone defects within the scope of internal fixation treatment or for endoprostheses revision surgery.

## Contraindication

**PALACOS® LV+G** with gentamicin must not be used in the following cases

- In case of suspected or proven hypersensitivity to components of the bone cement
- During pregnancy or breast-feeding
- In cases of severe renal insufficiency
- In the presence of an active or incompletely treated infection at the bone site caused by gentamicin non-sensitive strains
- **PALACOS® LV+G** is not suitable for primary or secondary plastic reconstruction of skull bone defects.

## Target population

There is very little evidence with children and adolescents, therefore it is not recommended to use **PALACOS® LV+G**. If no other option is available, for example with present surgical trauma, the decision whether to use **PALACOS® LV+G** lies with the attending surgeon.

## Target user group

Healthcare professionals in a clinical context.

## Side effects

### General

Rarely there is a temporary drop in blood pressure after preparation of the prosthesis bed or directly after the implantation of bone cements and the endoprosthesis.

In individual cases severe complications can occur such as severe allergic reactions which are associated with cardiac arrest, anaphylactic shock or even sudden death.

In order to avoid pulmocardiovascular complications such as pulmonary embolism and cardiac arrest it is recommended that the implantation site be irrigated thoroughly with an isotonic solution (pulse lavage application) before the bone cement is introduced. In the event of pulmonary or cardiovascular events it is necessary to monitor blood volume and possibly increase it.

In the case of acute respiratory insufficiency anaesthesiological measures should be taken.

The following adverse effects have been observed when using poly(methyl methacrylate) bone cements: thrombophlebitis, haemorrhage, trochanteric bursitis.

Other side effects that were observed are heterotopic bone formation, osteolysis due to cement fragments, myocardial infarction, brief cardiac arrhythmia, cerebrovascular incident.

### Gentamicin

Due to the gentamicin contained in **PALACOS® LV+G**, side effects that are typical for gentamicin can occur:

- impairment to auditory and vestibular nerves
- renal toxicity: Adverse renal effects have been reported, as demonstrated by the presence of cells or protein in the urine or by rising of serum creatinine or oliguria. They occur more frequently in patients with a history of renal impairment. Application of clinical doses of gentamicin has led to the occasional appearance of Fanconi's syndrome or a Bartter's-like syndrome
- neuromuscular blockade/ neurotoxicity: Serious adverse effects on both vestibular and auditory branches of the 8th cranial nerves have been reported, primarily in patients with renal impairment. Symptoms include dizziness, vertigo, tinnitus, roaring in the ears and hearing loss, which, as with the other aminoglycosides, may be irreversible. Other factors which may increase the risk of toxicity include dehydration and previous exposure to other ototoxic drugs
- rare cases of paresthesia, tetany and amyosthenia
- rare cases of allergic reactions (exanthema, urticaria, anaphylactic reactions).

Like all aminoglycosides gentamicin is also potentially nephrotoxic and/or ototoxic. In most cases, damage to hearing caused by gentamicin is irreversible, whereas renal damage is reversible.

Other reported adverse reactions possibly related to gentamicin include: respiratory depression, lethargy, confusion, depression, visual disturbances, decreased appetite, weight loss, and hypotension or hypertension; rash, itching, urticaria, generalised burning, laryngeal oedema, anaphylactoid reactions, fever, and headache, nausea, vomiting, increased salivation, and stomatitis; purpura, pseudotumor cerebri, acute organic brain syndrome, pulmonary fibrosis, alopecia, joint pain, transient hepatomegaly, and splenomegaly.

The use of **PALACOS® LV+G**, results in a desired high bactericidal but locally limited concentration of gentamicin.

Therefore, the occurrence of the afore-mentioned side effects is highly unlikely. Care should be taken and the level of gentamicin should be monitored in patients with impaired renal function, especially if treated with parenteral aminoglycosides or muscle relaxants simultaneously. The same applies to patients with previous neuromuscular disorders (e.g. myasthenia gravis, Parkinson's disease). Allergic reactions can occur dosage-independent.

## Interactions

### Gentamicin

Owing to the administration of muscle relaxants and ether, the neuromuscular blocking properties of gentamicin may be intensified. However, in view of the very low serum levels this is relatively unlikely in patients with healthy kidneys. The probability of interactions occurring increases in proportion to the serum levels of gentamicin, e.g. in patients with impaired renal function.

Interactions of gentamicin have been reported to occur with the following groups of drugs:

Concurrent or sequential use of neurotoxic and / or nephrotoxic pharmaceuticals, including other aminoglycosides, polymyxin B, colistin, cisplatin, ciclosporines, foscarnet, vancomycin, amphotericin B, clindamycin and cephalosporins.

Potent diuretics: Ethacrynic acid, furosemide, or other potent diuretics which may themselves cause ototoxicity or enhance aminoglycoside toxicity by altering antibiotic concentrations in serum and tissue.

Muscle relaxants: in particular succinyl choline and tubocurarine, decamethonium, halogenated hydro-carbon inhalation anaesthetics or opioid analgesics. Gentamicin may interact with these medications and result in skeletal muscle weakness and respiratory depression (apnea).

Concurrent application of these medications and gentamicin during surgery or in the post-operative period should be monitored with caution, especially if there is a possibility of incomplete reversal of neuromuscular blockade post-operatively.

Moreover, complications are possible that can occur during any surgical procedure.

## Precautions

### Use by operating staff

Before using **PALACOS® LV+G** the user should be well-acquainted with its properties, handling and application. The user is advised to practise the entire procedure of mixing, handling and introducing **PALACOS® LV+G** before using it for the first time. Detailed knowledge is necessary even if mixing systems and syringes are being used for application of the cement.

The liquid monomer is highly volatile and flammable. Ignition of monomer fumes caused by use of electrocautery devices in surgical sites near freshly implanted bone cements has been reported. The monomer is also a strong lipid solvent and should not come into direct contact with the body.

When handling the monomer or the prepared cement **PALACOS® LV+G** gloves must be worn which provide the necessary protection against penetration of the skin by the monomer methyl methacrylate.

Gloves made of PVP (three-ply polyethylene, ethylene vinyl alcohol copolymer, polyethylene) and Viton®/butyl gloves have proved to provide good protection over a lengthy period. For safety's sake it is recommended that two pairs of gloves be worn over one another, e.g. one polyethylene surgical glove over an inner pair of standard latex surgical gloves.

The use of latex or polystyrene-butadiene gloves on their own is inadequate. Please make enquiries with your supplier to establish which gloves are suitable for such an application. The monomer vapours can irritate respiratory passages and eyes and possibly cause damage to the liver. Skin irritations have been described which are due to contact with the monomer.

Manufacturers of soft contact lenses recommend removing the lenses in the presence of damaging or irritant vapours. Since soft contact lenses are permeable to liquids and gases, they should not be worn in the operating theatre if methyl methacrylate is being used.

#### Use on the patient

Blood pressure, pulse and breathing must be monitored carefully during and directly after introduction of the bone cement. Any significant change in these vital signs must be eliminated without delay by taking appropriate action.

When using **PALACOS® LV+G** the prepared bone should be carefully cleaned, aspirated and dried just before the bone cement is placed.

#### **Incompatibilities**

Aqueous solutions (e.g. ones containing antibiotics) must not be added to the bone cement because they have a considerable detrimental effect on the physical and mechanical properties of the cement.

#### **Amount required**

After mixing the cement powder with the monomer liquid, a slow-curing ductile dough develops which is introduced to bone cavities for anchorage and/or filling purposes.

**PALACOS® LV+G** is coloured green in order to make the cement clearly visible in the surgical field.

A dose is prepared by mixing the entire contents of the bag of cement powder with all the monomer liquid of an ampoule. The quantity of cement dough required depends on the specific surgical intervention and on the technique being used.

At least one additional dose of **PALACOS® LV+G** should be available before commencing the operation.


#### **Preparation**

##### Preparation:

Before opening the non-sterile aluminium protective pouch (compare above) move the contents down (1 polyethylene-paper bag) by shaking or tapping in order to ensure that when the bag is cut open at the top the content is not damaged.

The polyethylene-paper bag and the ampoule may only be opened under sterile conditions. For this purpose the sterile components (inner PE-paper bag and glass ampoule) are sterile when they are delivered.

#### Opening under sterile conditions:

The opening flaps from the bag top assist to detach the foil from the paper.

In order to grip as much of the opening flaps as possible, the side of the paper/foil should be kept between thumb, index and middle finger.

Please use the whole thumb surface to grip foil and paper side and take off each side evenly.

Open the outer PE-paper bag at the special point under sterile conditions so that the inner PE-paper bag remains sterile when it is removed. Also open the blister pack at the special point under sterile conditions so that the glass ampoule remains sterile when it is removed.

Before opening the inner PE-paper bag, moved the contents down by shaking or tapping in order to ensure that when the bag is cut open at the top no powder is lost. To make it easier to open the glass ampoule the latter is provided with a predetermined breaking point at the transition to the head of the ampoule.

Do not open the ampoule over the mixing device to prevent contamination of the cement with glass fragments.

The 10 ml- and 20 ml- ampoule is provided with a snapping-off device (tube) in order to facilitate the opening procedure. If this is the case, take hold of the fitted snapping-off device instead of the ampoule head and break off the ampoule head over it. When the ampoule head has snapped off, it remains inside the tube.

#### Mixing the components:

It is advisable to first measure out the liquid and then add the powder. If this order is reversed, powder nests are more likely to form as a result of polymerisation commencing immediately at the surface. Both components, i.e. the relative proportions of powder and monomer, are precisely matched. The pouch and ampoule must therefore be emptied completely if an optimal mix is to be achieved.

The components can be mixed with a vacuum mixing system or by hand.

The mixing, waiting, working and curing times of **PALACOS® LV+G** are shown in the diagrams at the end of the instructions for use. Please note that these are stated for guidance only because working time and curing time depend on temperature, mixing and humidity, whereby the direct ambient temperatures are important, e.g. of cement powder, mixing system, bench and hands. A high temperature accelerates the waiting, working and curing times.

#### **Preparation with a vacuum mixing system**

In order to obtain a bone cement with a reduced number of air inclusions, the liquid and the powder are mixed in a vacuum. For this purpose an airtight mixing system has to be used, ensuring that a sufficient vacuum builds up in the mixing vessel quickly (pressure approx. 200 mbar absolute).

Perform filling and mixing under sterile conditions.

The mixing time is 30 seconds unless otherwise recommended. For details of mixing technique please refer to the instructions for the mixing system being used.

Always mix the entire contents of a bag with the entire contents of an ampoule of monomer liquid.

**Preparation mixing by hand**

The cement components should only be filled into the mixing vessel just before mixing. Perform filling and mixing under sterile conditions – mixing time is 30 seconds. During this time the two components are mixed with one another while stirring evenly. Always mix the entire contents of a bag with the entire contents of an ampoule of monomer liquid.

**Use of the bone cement**

The bone cement can be applied as soon as the doughy bone cement no longer adheres to the gloves. The application period depends on the temperature of the material and the room temperature. To ensure adequate fixation the prosthesis must be introduced and held within the time window allowed for application, until the bone cement has set completely. Remove any surplus cement as long as it is still soft.

**Storage**

Do not store above 25 °C (77 °F).

**Shelf life/sterility**

The shelf life is printed on the folding box, the protective aluminium bag and the inner bag. Do not use **PALACOS® LV+G** if this date indicated has expired. The information on the blister pack may be different for production reasons. The contents of opened or damaged aluminium sachets or ampoule blisters must not be re-sterilized and thus must be discarded. If the cement powder has turned yellow, do not use **PALACOS® LV+G**. The cement powder and packagings are sterilised by gassing with ethylene oxide. The monomer liquid has been sterilised by filtration.

**Disposal**

Single components of the bone cement, cured solid material as well as (uncleaned) packaging material must be disposed by following the regulations of the local authorities. Dispose of the polymer component in an authorized waste facility. The liquid component should be evaporated under a well-ventilated hood or absorbed by an inert material and transferred in a suitable container for disposal.

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Heraeus Medical GmbH