
Heraeus

OSTEOPAL[®] plus

Radiopaque bone cement
for filling and stabilising vertebral bodies
with an extended application time

OSTEOPAL® plus

| | | |
|---------|---|------------------------|
| EN | Instructions for use Working times (Revision status: 2019-12) | 5 90 |
| EN (US) | Note: Instructions for use for sale within the US Instructions for use Working times (Revision status: 2019-12) | 8 90 |
| BG | Инструкция за употреба Времена за обработка (Версия 2019-12) | 12 90 |
| CS | Návod k použití Doby zpracování (Poslední revise: 2019-12) | 15 90 |
| DA | Brugsanvisning Bearbejdningstider (Aktualiseret: 2019-12) | 18 90 |
| DE | Gebrauchsanweisung Verarbeitungszeiten (Revisionsstand: 2019-12) | 21 90 |
| EL | Οδηγίες χρήσης Χρόνοι εργασίας (Ημερομηνία αναθεώρησης: 2019-12) | 24 90 |
| ES | Instrucciones de uso Tiempos de manipulación (Última revisión: 2019-12) | 27 90 |
| ET | Kasutusjuhend Tööajad (Viimati muudetud: 2019-12) | 30 90 |
| FI | Käyttöohje Työstöajat (Revisiostatus: 2019-12) | 33 90 |
| FR | Notice d'utilisation Temps de travail (Date de révision: 2019-12) | 36 90 |
| HR | Upute za uporabu Vremena obrade (Stanje revizije: 2019-12) | 39 90 |
| HU | Használati utasítás Megmunkálási idők (Aktualizálás időpontja: 2019-12) | 42 90 |

| | | |
|----|--|------------------------|
| IT | Istruzioni per l'uso Tempi di lavorazione (Ultima revisione: 2019-12) | 45 90 |
| LT | Naudojimo instrukcija Darbo laikas (Tikrinta: 2019-12) | 48 90 |
| LV | Lietošanas instrukcija Apstrādes laiks (Pārbaudes statuss: 2019-12) | 51 90 |
| NL | Gebruiksaanwijzing Verwerkingstijden (Revisiestatus: 2019-12) | 54 90 |
| NO | Brugsanvisning Virketider (Revisjonsstatus: 2019-12) | 57 90 |
| PL | Instrukcja użytkowania Czasy obróbki (Ostatnia aktualizacja: 2019-12) | 60 90 |
| PT | Instruções de utilização Fases de aplicação (Data da última actualização: 2019-12) | 63 90 |
| RO | Instrucțiuni de utilizare Timpi de prelucrare (Stadiu revizie: 2019-12) | 66 90 |
| RU | Инструкция по применению Временные фазы для смешивания (Ревизионный статус: 2019-12) | 69 90 |
| SK | Návod na použitie Časy spracovateľnosti (Posledná revízia textu: 2019-12) | 72 90 |
| SL | Navodilo za uporabo Časovne opredelitve (Datum revizije: 2019-12) | 75 90 |
| SV | Bruksanvisning Hanteringstider (Revision status: 2019-12) | 78 90 |
| TH | คำแนะนำวิธีใช้ เวลาในการใช้งาน (แก้ไขปรับปรุงครั้งล่าสุดเมื่อ: 2019-12) | 81 90 |
| TR | Kullanma Talimatı Çalışma süreleri (Bilgilerin son güncellenme tarihi: 2019-12) | 84 90 |

ZH

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Properties

OSTEOPAL® plus is a radiopaque, low-viscosity bone cement, based on polymethylmethacrylate with an extended application phase, used to fill and stabilise vertebral bodies.

OSTEOPAL® plus contains zirconium dioxide as an X-ray contrast agent. The colouring agent chlorophyll (E141) is mixed with **OSTEOPAL® plus** to improve visibility of the surgical field.

The bone cement is prepared immediately prior to use by mixing the polymer powder component and the liquid monomer component. A low viscosity paste is created, and with the aid of an application system placed in the vertebral body, where it cures.

Composition

A pack of **OSTEOPAL® plus** 1x20 contains a pouch of cement polymer powder and a brown glass ampoule containing the cement monomer liquid.

Constituents of the polymer powder:

poly(methyl acrylate, methyl methacrylate), zirconium dioxide, benzoyl peroxide, colouring agent E141.

Constituents of the monomer liquid:

methylmethacrylate, N,N-dimethyl-p-toluidine, hydroquinone.

The packaging for the cement powder is sterile. The outer, non-sterile aluminium protective bag contains a polyethylene-paper pouch, which is sterile on the inside. This pouch contains an additional sterile PE-paper bag that contains the cement powder.

The brown glass ampoule containing the sterile-filtered monomer liquid is also sterile packed in an ethylene oxide sterilised individual blister pack.

Indications

OSTEOPAL® plus is used for augmenting and stabilising of vertebral bodies.

- for compression fractures of the vertebral body
- for vertebral body tumours (metastases and myelomas)
- for symptomatic vertebral haemangiomas

In every case, the vertebroplasty and kyphoplasty only represent stabilising forms of treatment. This does not constitute therapy of the underlying disease (osteoporosis, tumour disease).

Contraindications

Haemorrhagic diathesis and infections represent absolute contraindications. Lesions of the vertebral body with epidural extension are relative contraindications due to the risk of spinal cord compression.

OSTEOPAL® plus must not be used

- in case of suspected or proven hypersensitivity to components of the bone cement
- during pregnancy and breast-feeding

Target patient population

Adult population.

There is very little evidence with children and adolescents. Therefore, it is not recommended to use **OSTEOPAL® plus**. If no other option is available, the decision whether to use **OSTEOPAL® plus** lies with the attending surgeon.

Side effects

There are frequent reports of a temporary reduction in blood pressure immediately after implantation of PMMA bone cements. Less frequent cases of hypotension accompanied by anaphylaxis have been described; including anaphylactic shock, cardiac arrest and sudden death. These cardiovascular and respiratory side effects, also known as implantation syndrome or bone cement syndrome, result from infiltration of bone marrow components into the venous vascular system. In the case of pulmonary or cardiovascular complications, monitoring and possibly an increase in blood volume are necessary. In the event of acute respiratory insufficiency, anaesthesiological measures should be undertaken.

Moreover, complications which may occur with any surgical procedure are conceivable.

The requirements for a contemporary cementing technique must be taken into consideration when using

OSTEOPAL® plus to limit undesirable side effects and to be in a position to guarantee stable and long-lasting anchoring of the bone cement in the vertebral body.

Interactions

Are not yet known.

Precautions

Use by surgical personnel

Prior to using **OSTEOPAL® plus**, the user should be familiar with its properties, handling and application. The user is recommended to practice the mixing and application process in advance. Detailed knowledge of the mixing systems and syringes for applying the cement is a prerequisite.

The liquid monomer is highly volatile and flammable. Ignition of monomer fumes caused by use of electrocautery devices in surgical sites near freshly implanted bone cements has been reported. The monomer is a strong lipid solvent and should not come into direct contact with the body.

When handling the monomer or the cement, it is always necessary to wear protective gloves that ensure the necessary protection against ingress of the monomer (methylmethacrylate). PVP gloves (three-layer polyethylene, ethylene vinyl alcohol copolymer, polyethylene) and Vitonbutyl gloves are highly recommended over a long period of use. It is recommended to wear two pairs of gloves over each other (to be on the safe side), e.g. a polyethylene surgical glove over an inner latex standard surgical glove.

The use of latex or polystyrene-butadiene gloves alone is insufficient. In addition, the relevant information should be obtained in advance from the manufacturer/supplier.

The monomer vapours can irritate the respiratory tract and eyes, and possibly cause liver damage. Skin irritations have been described which are attributable to contact with the monomer.

The manufacturers of soft contact lenses recommend removing these lenses from the eyes prior to handling MMA and the irritant vapours that result. As soft contact lenses are permeable to liquids and gases, they should not be worn in the operating theatre if methylmethacrylate is used.

Use on the patient

Blood pressure, pulse and breathing must be carefully monitored during and immediately after implantation of the cement. Every significant change in these vital signs must be rectified with the appropriate measures.

The leakage of cement can cause damage to the paravertebral structures. Complications, such as spinal cord compression, intercostal neuralgia, escape of cement into the intervertebral space, erivertebral filling of veins and arteries (risk of embolism), infections and post-procedural pain are possible. In order to avoid the escape of cement and to recognise adverse events in good time, the application must take place using imaging methods (real time display).

The means of immediate operative procedure to surgically rectify the complications described must also be available. A thorough pre-operative radiological examination must be carried out to preclude possible risks (e.g. vertebral body lesions, vascularisation of the vertebral body or oedema). Insufficient filling of the vertebral body with bone cement can lead to insufficient acute pain reduction and reduced long-term stability of the vertebral body treated.

Incompatibilities

Aqueous solutions (e.g. containing antibiotics) must not be mixed with the bone cement, as they significantly impair the physical and mechanical properties of the cement.

Dosage

A dose is prepared by mixing the complete contents of a cement powder pouch with all the monomer liquid in an ampoule. The quantity to be applied depends on the respective anatomical conditions. There is no information available on the maximum amount of bone cement that can be inserted and the maximum number of vertebral bodies to be treated for a vertebroplasty or kyphoplasty.

Generally one or two portions are used. However, this depends on the surgical technique applied, as well as the size of the respective vertebral body and the defect.


At least one additional pack of **OSTEOPAL® plus** should be available before commencing the operation.

Procedure for use

Preparation:

Prior to opening the unsterile aluminium protective pouch, the contents are shaken or tapped down so that the contents are not damaged on opening. The inner polyethylene-paper pouch containing the powder and the ampoule must only be opened under sterile conditions. The sterile components (inner polyethylene-paper pouch and glass ampoule) are handed over under sterile conditions.

Opening under sterile conditions:

The opening flaps from the bag top assist to detach the foil from the paper.

In order to grip as much of the opening flaps as possible, the side of the paper/foil should be kept between thumb, index and middle finger.

Please use the whole thumb surface to grip foil and paper side and take off each side evenly.

The outer peel-off pouch is opened under sterile conditions and in the intended place so that the inner polyethylene-paper pouch can be removed sterile. The blister pack is also opened under sterile conditions and at the intended place so the glass ampoule can be removed sterile.

Before opening the inner polyethylene-paper pouch, the contents are moved down by shaking or tapping to ensure that there is no loss of powder when the top edge is cut open.

Do not open the ampoule over the mixing device to prevent contamination of the cement with glass fragments. To make it easier to open the glass ampoule, it has a predetermined breaking point at the transition to the head of the ampoule. A breaking aid (tube) is attached to the glass ampoule to make it easier to open the ampoule. The breaking aid is gripped rather than the ampoule head and the ampoule head is broken off beyond this. The broken off ampoule head remains in the breaking aid.

Mixing the components:

The liquid from the ampoule is put in a sterile mixing vessel. The powder from the open inner pouch is then added. The mixture is stirred with a sterile spatula or spoon until a homogeneous paste is produced. The mixture should be stirred for 30 seconds irrespective of the ambient temperature. Strict adherence with the instructions for mixing the powder and liquid components can minimise complications.

The cement components can be mixed in a mixing system with or without a vacuum. The mixing time for mixing with or without a vacuum is 30 seconds irrespective of the ambient temperature. Details on the mixing systems may be found in their instructions for use.

The processing time and polymerisation are strongly dependent on the temperature and environment. Higher temperatures accelerate the curing time, lower temperatures extend it.

Processing

The viscosity increases with advancement of polymerisation, i.e. the advancement of the processing phase. The cement paste should be put into an application system immediately after mixing, as it still has low viscosity and may be easy to aspirate. To avoid vascular escape of cement, the cement should be applied in a paste state.

OSTEOPAL® plus can be inserted into the vertebral body with the aid of an application system approved for percutaneous vertebroplasty or kyphoplasty, which allows constant and controlled injection. The handling of the system is described in the manufacturer's instructions for use. During intravertebral application, constant X-ray screening (latero-

lateral) in real time is necessary. In the event of paravertebral escape of cement, cement injection must be interrupted immediately and can be continued after the viscosity of the cement has increased. If the vertebral filling is not sufficient, a further contralateral approach is possible. After augmentation, a mandrin should be inserted in the injection needle so that no residues of cement remain in the soft tissue after removal of the injection needle.

The mixing, processing and curing times of **OSTEOPAL® plus** are shown in the diagram at the end of these instructions.

The values apply to using bone fillers with a diameter of 3.5 mm (MAXXSPINE Ltd, 65307 Bad Schwalbach, Germany). The cement can still be removed from the bone fillers at the end of the processing period, but no longer binds free from creases. Therefore, the end of the processing time refers to the state of the cement, not to the possible end time for removing cement from the bone fillers. There may be changes in the processing for other application systems. Cannulas with a diameter below 1.8 mm (1.3G) should not be used. The patient must remain immobilised until the cement is fully cured.

Storage

The cement must be stored unopened and protected from light at a maximum temperature of 25 °C (77 °F) in a dry, clean place in the original packaging.

Shelf life/sterility

The expiry date is specified on the folding box, the aluminium pouch and the inner pouch.

OSTEOPAL® plus must not be used after the date specified. For production reasons, the shelf lives of the individual components may differ from that specified on the folding box.

The content of non-used, opened or damaged packages must not be resterilised and must therefore be discarded. If the cement powder has turned yellow, **OSTEOPAL® plus** must no longer be used.

Disposal

Single components of the bone cement, cured solid material as well as (uncleaned) packaging material must be disposed by following the regulations of the local authorities. Dispose of the polymer component in an authorized waste facility. The liquid component should be evaporated under a well-ventilated hood or absorbed by an inert material and transferred in a suitable container for disposal.

EN (US)

OSTEOPAL® plus

Note: Instructions for use for sale within the US

Uses and properties

OSTEOPAL® plus is a setting polymer for use in augmentation and stabilization of vertebral bodies. Mixing of the two component system, consisting of a powder and a liquid, initially produces a paste, which is injected in the vertebral body for stabilization. Insoluble zirconium(IV)oxide is included in the cement powder as an X-ray contrast medium. The chlorophyll additive serves as optical marking of the bone cement at the site of the operation.

Indications

OSTEOPAL® plus bone cement is indicated for the treatment of pathological fractures of the vertebral body due to osteoporosis, cancer, or benign lesions using a vertebroplasty or balloon kyphoplasty procedure.

Contraindications

- The use in cases of haemorrhaging diathesis and infections are absolutely contraindicated. Lesions of the vertebral body with epidural extension are relative contraindications (danger of spinal cord compression).
- **OSTEOPAL® plus** must not be used in cases of known hypersensitivity to constituents of the bone cement.
- **OSTEOPAL® plus** must not be used during pregnancy or the nursing period.
- Allergies or sensitivity to any of its chemical components
- Cases of active or incompletely treated infections
- Coagulation disorders, or with severe cardiopulmonary disease
- Spinal stenosis (>20% by retropulsed fragments)
- Compromised of the vertebral fractures due to posterior involvement
- Patient clearly improving on medical therapy
- Prophylaxis in metastatic or osteoporotic patients with no evidence of acute fracture

- Non-pathological, acute traumatic fractures of the vertebra
- Vertebral plana (collapse >90%)
- Compromise of the vertebral body or the walls of the pedicles
- Unstable vertebral fractures due to posterior involvement

Target patient population

Adult population.

There is very little evidence with children and adolescents. Therefore, it is not recommended to use **OSTEOPAL® plus**. If no other option is available, the decision whether to use **OSTEOPAL® plus** lies with the attending surgeon.

Warnings

- Read and understand instructions
- Have specific training and experience to be thoroughly familiar with the properties, handling characteristics, and application of the product and percutaneous cement delivery.
- Heraeus Medical as a manufacturer does not recommend a surgical technique, it is the responsibility of the physician to determine the appropriateness of this product and specific technique for each patient.
- Hypotensive reactions have occurred between 10 and 165 seconds. They have lasted from 30 seconds to 5 minutes. Some patients have progressed to cardiac arrest. Monitor patients carefully for any change in blood pressure during and immediately following the application of bone cement. Adverse patient reactions affecting the cardiovascular system have been associated with the use of bone cements.
- MMA can be hypersensitivity in susceptible persons, which may result in anaphylactic response.
- Use in pregnancy and children has not been established.

| Composition | OSTEOPAL® plus 1x20 |
|---|--|
| | 1 pouch of 26.0g powder contains: |
| poly (methyl acrylate, methyl methacrylate) | 14.2g |
| zirconium dioxide | 11.7g |
| benzoyl peroxide | 0.1g |
| | 1 ampoule of 10ml liquid contains: |
| methyl methacrylate | 9.35g |
| N,N-dimethyl-p-toluidine | 0.05g |
| Other constituents: | In the powder: chlorophyll VIII In the liquid: hydroquinone |

- Not recommended in patients that do not exhibit a pathological condition, such as primary or secondary osteoporosis or a tumor, which would impair the ability of the patient to heal using conservative treatment methods.
- Strict adherence to good surgical principles and techniques. Deep wound infection is a serious postoperative complication and may require total removal of the embedded cement. Deep wound infection may be latent and not manifest itself even for several years postoperatively.
- Follow handling, mixing and preparation instructions carefully.
- Do not resterilize. Single patient use only. Sterile only if package is unopened and undamaged.
- Caution should be used to prevent excessive exposure to the concentrated vapors of the monomer, which may produce irritation of the respiratory tract, eyes and possible liver.
- Do not use monomer if the liquid component shows any sign of thickening or premature polymerization.
- Do not allow the liquid component to contact rubber or latex gloves. The liquid component is a powerful lipid solvent. Should contact occur, the gloves may dissolve and tissue damage may occur. Wearing a second pair of gloves may diminish the possibility of hypersensitivity reactions.
- Do not allow personnel wearing contact lenses to be near or involved in mixing of the bone cement.
- Use appropriate imaging techniques to verify correct needle placement, absence of damage to surrounding structures, and appropriate location of the injected bone cement. Use imaging, such as fluoroscopy, to assess the ability of the vertebra to contain the injected bone cement.
- Avoid over pressurization of the bone cement because this may lead to leakage of the bone cement beyond the site of its intended application.
- Cement leakage may cause tissue damage, nerve or circulatory problems, and other serious adverse events.
- Leaks can also occur when injecting if the needle is in a vein or if unseen microfracture are prevalent.
- If bone cement is seen outside of the vertebral body or in the circulatory system during the procedure, immediately stop the injection.
- Inadequate fixation or unanticipated postoperative events may affect the cement-bone interface and lead to micro-motion of the cement against bone surface. A fibrous tissue layer may develop between the cement and the bone. Long-term follow-up is advised for all patients on a regular scheduled basis.
- The completion of polymerization occurs in the patient and is an exothermic reaction with considerable liberation of heat. The released heat may damage bone or other tissues surrounding the implant. According to ISO 5833 standard, the temperature can be as high as 95 °C. The long term effect of the heat produced in situ, positioning of the patient should maintained securely to obtain proper fixation. For proper fixation, 1 – 2 hours or longer may be required, as determined by the patient's medical condition and the attending physician.
- The long term effects of the bone cement in the spine have not yet been established.

Precautions

- Do not use this product after expiration date printed on the package. This device may not be safe or effective beyond its expiration date.
- Follow the handling and mixing instructions carefully to avoid contact dermatitis. Strict adherence to the instructions for mixing the powder and liquid components may reduce the incidence of this complication.
- Adequately ventilate the operating room to eliminate as much monomer vapor as possible. The liquid monomer is highly volatile and flammable. Ignition of monomer fumes caused by the use of electrocautery devices in surgical sites near freshly implanted bone cements has been reported.
- When working with the monomer or the cement, gloves must be worn to assure adequate protection against the penetration of the monomer (methyl methacrylate) into the skin. PVP (three-layer polyethylene, ethylene-vinyl alcohol-copolymer, polyethylene) and Viton/Butyl gloves have proved to provide good protection over an extended period. Putting on two pairs of gloves – one pair of polyethylene surgical gloves over a pair of standard latex surgical gloves has also proved to offer adequate protection. The use of latex or polystyrene-butadiene gloves, however, must be avoided. Please request the confirmation of your glove supplier whether the respective gloves are suitable for the use of this cement.

Potential adverse events associated with vertebroplasty procedure include

- Pneumonia
- Intercostal neuralgia
- Collapse of a vertebra adjacent to the one injected, due to an osteoporotic disease
- Pneumothorax
- Cement extrusion into soft tissue
- Fracture of a pedicle
- Rib fracture in patients with diffuse osteopenia, especially during thoracic vertebroplasty procedures, due to the significant downward force exerted during needle insertion
- Compression of the spinal cord with paralysis or loss of feeling
- Cement leakage between intervertebral discs.
- Serious adverse events, some with fatal outcome, associated with the use of acrylic bone cements for vertebroplasty or kyphoplasty include myocardial infarction, cardiac arrest, cerebrovascular accident, pulmonary embolism, and cardiac embolism. Although the majority of these adverse events present early within the post-operative period, there have been some reports of diagnoses beyond a year or more after the procedure. Other reported adverse events for acrylic bone cements intended for vertebroplasty or kyphoplasty include: Leakage of the bone cement beyond the site of its intended application with introduction into the vascular system resulting in embolism of the lung and/or heart or other clinical sequelae.
- The most frequent adverse reactions reported with acrylic bone cements are a transitory fall in blood pressure, thrombophlebitis, hemorrhage and hematoma, loosening or displacement of the prosthesis, superficial or deep wound infection, trochanteric bursitis, and short-term

cardiac conduction irregularities. Other reported adverse events for acrylic bone cements include pyrexia due to an allergy to the bone cement, hematuria, dysuria, bladder fistula, and delayed sciatic nerve entrapment due to extrusion of the bone cement beyond the region of its intended application.

Important Physician Information

- Percutaneous vertebroplasty procedures should only be performed in medical settings in which emergency decompressive surgery is available.
- Adverse reactions affecting the cardiovascular system have been attributed to recent data that indicate the monomer undergoes rapid hydrolysis to methacrylic acid, and that a significant fraction of the circulating methacrylate is in the form of circulating concentrations of methyl methacrylate/methacrylic acid and changes in blood pressure has not been established.
- The physician is responsible for any complication or harmful consequences, which may result from an erroneous indication or operating technique, from inappropriate use of the material, or non-observation of the safety instructions that appear in the directions for use.
- Additives (such as antibiotics) are not to be mixed with the bone cement, as this will alter cement properties.

Interactions

Not known.

Incompatibilities

Aqueous (e.g. antibiotic-containing) solutions must not be mixed with the bone cement, as this reduces the strength considerably.

Dosage

The amount of bone cement required depends on the patient's anatomy. There is no study about the maximal amount of bone cement or the maximal amount of treatable vertebral bodies by vertebroplasty. This is subject to the surgeon's responsibility. One or more complete units (the contents of one sachet and one ampule) must always be mixed together.

Notes on use


Before using **OSTEOPAL® plus** the surgeon must be thoroughly familiar with its properties, handling and use and must have read the relevant literature. In particular, the surgeon must be proficient in the special application techniques required for low viscosity cements. An appropriate up to date cementing technique must be used with **OSTEOPAL® plus** to limit the side effects and to achieve a stable and long-lasting anchorage of the cement in the vertebral body.

The application must be performed with an application instrument certified for vertebroplasty. For special application techniques, the relevant instructions of the application system must be consulted by the surgeon. In order to recognize immediately cement leakage, an imaging procedure (real time) during the application of the cement is recommended.

For recognition of a vertebral extension a phlebography before injection of the cement is recommended. The protective aluminium packaging and outer non-sterile

polyethylene sachet („peel-off“ package) and the ampule packaging should be opened maintaining sterility by the circulating nurse. The ampule and powder sachet are then taken out under sterile conditions and placed on a sterile working surface. The ampule is opened by breaking the neck and the inner sachet is cut open with sterile scissors.

Opening under sterile conditions:

The opening flaps from the bag top assist to detach the foil from the paper.

In order to grip as much of the opening flaps as possible, the side of the paper/foil should be kept between thumb, index and middle finger.

Please use the whole thumb surface to grip foil and paper side and take off each side evenly.

Mixing

1. Empty the entire content of the liquid ampoule into an appropriate sterile, inert mixing device. Add the entire contents of the powder to the liquid. Don't add the liquid to the powder.
2. Mix the liquid and powder components with a sterile mixing spoon or porcelain or stainless-steel spatula for a minimum of 30 sec and until the powder is completely saturated with the liquid. The liquid monomer and the powder component should be vigorously mixed to create a homogenous mixture, which is essential to obtain consistent cement. Follow the instructions provided by the mixing device manufacturer.
3. Care must be taken to minimize air entrapment within the cement when using a cement delivery system.
4. The physician should, by specific training and experience, be thoroughly familiar with the properties, handling characteristics, application of the bone cements, and amount of cement required per procedure. Because the handling and curing characteristics of this bone cement may vary with temperature, humidity, and mixing technique, they are best determined by the physician's actual experience.
5. Do not apply the cement until two minutes after completely mixing the liquid component to the powder component.
6. Select an appropriate method to apply the cement to bone. If applicable, follow the device manufacturer's instructions.

Working procedure

The working time and the rate of polymerization are strongly dependent on ambient temperature and on temperature of the components. The hardening time is shorter at high temperatures and longer at low temperatures. Viscosity increases as polymerization progresses, i.e. as the cement hardens. After mixing, the cement is poured directly into a cement syringe, since at this time **OSTEOPAL® plus** still has a low viscosity and is easily pourable.

The application of the cement should be performed during a pasty condition phase of the cement to avoid vascular migration. **OSTEOPAL® plus** can be applied with a percutaneous vertebroplasty certified cement syringe into the vertebral body which allows a controlled constant injection.

The method of application is given in the instructions for the system used by the manufacturer and must be consulted.

During intravertebral application a strict latero-lateral real time radioscopy is necessary. In case of paravertebral leakage the injection must be stopped immediately and can be proceeded when the viscosity of the cement is increased. In cases of insufficient vertebral filling a contralateral approach can be performed. After vertebral augmentation a mandrin should be placed into the vertebroplasty injection needle to avoid cement residues being deposited in the soft tissues on withdrawal of the application needle. At a room and a material temperature of 17°C, 20°C and 23°C (63°F, 68°F and 73°F) the schedule of the temperature/time diagram applies:

The reference data belong to bonefiller (MAXXSPINE Ltd, 65307 Bad Schwalbach, Germany) with a canula diameter of 3.5mm. This time may vary depending on the syringe system used. No canula with a diameter less than 1.8mm should be used. The patients should be kept immobilised until the end of the setting time as illustrated in the graph.

Other uses

OSTEOPAL® plus is not certified for other indications than percutaneous vertebroplasty and kyphoplasty. The application in another indication (i.e. endoprosthesis) is not investigated.

Storage

Do not store above 25°C (77°F).

Shelf life/sterility

The expiry date is printed on the outer carton, on the protective aluminium packaging and on the inner sachet.

OSTEOPAL® plus must not be used after the expiry date. Once removed from the protective packaging, the powder sachet cannot be resterilized and must therefore be discarded.

OSTEOPAL® plus is sterilized with ethylene oxide gas and must not be resterilized. If the polymer powder shows a yellow discoloration, it must not be used.

Disposal

Single components of the bone cement, cured solid material as well as (uncleaned) packaging material must be disposed by following the regulations of the local authorities.

Dispose of the polymer component in an authorized waste facility. The liquid component should be evaporated under a well-ventilated hood or absorbed by an inert material and transferred in a suitable container for disposal.

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Heraeus Medical GmbH