

PALAMIX[®]

QUALITY AND FLEXIBILITY IN MIXING



HIGH CEMENT QUALITY AND SAFETY


PALAMIX®: FOR SUCCESSFUL SURGICAL OUTCOME

Using modern cementing technique in arthroplasty is key to ensure long survival rates for endoprostheses and a low risk of revision.^{1,2} One of its crucial success factors is a homogeneously mixed bone cement.³

Vacuum mixing improves cement homogeneity by reducing porosity and strengthens the cement prosthesis interface.⁴ The disposable vacuum mixing system PALAMIX® enables standardised mixing of homogeneous bone cement of reproducibly high quality.


VACUUM MIXING

- porosity of 0.1–1%⁵
- substantially decreases cement implant interface voids⁴
- increases force required to separate cement and implant and improves fatigue strength⁴
- reduces the longer term risk of revision in THA¹
- substantially reduces MMA* fumes (<10 ppm)^{6,7}



NON-VACUUM MIXING

- porosity of 5–16%⁵
- weaker cement implant bonds⁴
- decreases the fatigue life of the cement³
- higher MMA fumes exposure⁶



SAFETY FOR STAFF AND PATIENTS



To ensure a high level of safety, PALAMIX® has a filling funnel with two separate sections for the bone cement components. A particle filter in the liquid chamber protects against small glass particles in the bone cement. It avoids injuries and preserves the cement's mechanical quality.

Additionally, the PALAMIX® vacuum hose comes with an active charcoal filter that minimises MMA fumes during mixing.⁷

* MMA=methyl methacrylate


(1) Malchau H, Herberts P, Ahnfelt L. Prognosis of total hip replacement in Sweden. Follow up of 92,675 operations performed 1978–1990. Acta Orthop Scand. 1993; 64(5): 497–506.
(2) Breusch S, Malchau H. The Well Cemented Total Hip Arthroplasty. Theory and Practice. Springer Verlag 2005; 147–148.
(3) Dunne NJ et al. The relationship between porosity and fatigue characteristics of bone cements. Biomaterials 2003; 24(2): 239–245.
(4) Geiger MH et al. The clinical significance of vacuum mixing bone cement. Clin Orthop Relat Res 2001; 382: 258–266.
(5) Wang JS. The Benefit of Vacuum Mixing. In: The Well Cemented Total Hip Arthroplasty 2005. Springer Verlag; 107–112.

MORE CONVENIENCE AND FLEXIBILITY

MIXING AND COLLECTING

With PALAMIX®, the mixing of the components and the collection of the bone cement is carried out under vacuum, which

- facilitates mixing
- saves preparation time before application
- reduces porosity⁵
- results in a homogeneously mixed bone cement that builds the basis for long term success in arthroplasty.⁸




APPLYING AND PRESSURISING

PALAMIX® provides various options for application and pressurisation of bone cement. Different nozzles are available for individual procedures and allow flexibility in surgery. With PALAMIX®, it is even possible to expel the mixed bone cement remaining in the nozzle, thereby minimising cement waste. A spatula clip can be attached to the short nozzle to apply a thin cement layer on the bone or prosthesis.

Pressurisation results in greater penetration of the bone, improved bone cement interface and increased fatigue strength of the cement.⁹ The PALAMIX® system contains two pressurisers for both, hip and knee arthroplasty. The soft material and shape of the femoral pressuriser enhance the sealing of the femur during cement application.¹⁰

In knee surgery, the included knee pressuriser helps to increase the cement's interdigitation into the tibial plateau.¹⁰



(6) Kuehn KD. PMMA Cements. Springer Verlag 2014; 262.

(7) Jelecevic J et al. Methyl methacrylate levels in orthopedic surgery: comparison of two conventional vacuum mixing systems. Ann Occup Hyg. 2014; 58(4): 493–500.

(8) Breusch SJ, Kühn KD. Bone cements based on polymethylmethacrylate. Der Orthopäde 2003; (32): 4–50.

(9) Wang JS. Femoral Pressurisation. In: The Well Cemented Total Hip Arthroplasty 2005. Springer Verlag; 160–163.

(10) Data on file at Heraeus Medical GmbH.

PALAMIX®

ADVANTAGES AT A GLANCE

- homogeneous bone cement of reproducibly high quality
- convenient and intuitive handling
- time saving thanks to collection under vacuum
- various application options for flexibility in surgery
- safety for staff and patients



PALAMIX®	Description	Content	REF
PALAMIX® uno	Vacuum mixing system with collection under vacuum for up to two pouches (2x40)	10	66057893
PALAMIX® duo	Vacuum mixing system with collection under vacuum, with two cartridges for up to four pouches (4x40)	10	66057897
PALAMIX® medium nozzle	Flexible, conical nozzle; Ø 8.7–12.6 mm	10	66043960
PALAMIX® slim nozzle	For use with low-viscosity bone cements; Ø 7 mm	10	66036747
PALAMIX® cement gun	Reusable cement gun	1	66036163
PALAMIX® vacuum pump	Vacuum pump	1	66036748

Simply order from Heraeus.

www.heraeus-medical.com