
COPAL[®] spacem

High viscosity, radiopaque bone cement
for creating spacers in two-stage revision

COPAL® spacem

EN	Instructions for use Working times (Revision status: 2020-03)	4 80
CS	Návod k použití Doby zpracování (Poslední revise: 2020-03)	7 80
DA	Brugsanvisning Arbejdstider (Aktualiseret: 2020-03)	10 80
DE	Gebrauchsanweisung Verarbeitungszeiten (Revisionsstand: 2020-03)	13 80
EL	Οδηγίες χρήσης Χρόνοι εργασίας (Ημερομηνία αναθεώρησης: 2020-03)	17 80
ES	Instrucciones de uso Tiempos de manipulación (Última revisión: 2020-03)	21 80
ET	Kasutusjuhend Tööjad (Viimati muudetud: 2020-03)	24 80
FI	Käyttöohje Työstöajat (Revisiostatus: 2020-03)	27 80
FR	Notice d'utilisation Temps de travail (Date de révision: 2020-03)	30 80
HU	Használati utasítás Megmunkálási (Aktualizálás időpontja: 2020-03)	34 80
IT	Istruzioni per l'uso Tempi di lavorazione (Ultima revisione: 2020-03)	37 80
LT	Naudojimo instrukcija Darbo laikas (Tikrinta: 2020-03)	40 80
LV	Lietošanas instrukcija Apstrādes laiks (Pārbaudes statuss: 2020-03)	43 80

NL	Gebruiksaanwijzing Verwerkingstijden (Revisiestatus: 2020-03)	46 80
NO	Brugsanvisning Virketider (Revisjonsstatus: 2020-03)	50 80
PL	Instrukcja użytkowania Czasy obróbki (Ostatnia aktualizacja: 2020-03)	53 80
PT	Instruções de utilização Fases de aplicação (Data da última actualização: 2020-03)	57 80
PT (BR)	Nota: Apenas para fins brasileiros! Instruções de utilização Fases de aplicação (Data da última actualização: 2020-03)	60 80
RU	Инструкция по применению Временные фазы для смешивания (Ревизионный статус: 2020-03)	64 80
SK	Návod na použitie Časy spracovateľnosti (Posledná revízia textu: 2020-03)	68 80
SL	Navodilo za uporabo Časovne opredelitve (Datum revizije: 2020-03)	71 80
SV	Bruksanvisning Hanteringstider (Revision status: 2020-03)	74 80
ZH	使用说明书 加工时间 (修订状态: 2020年3月)	第77页 第80页

Properties

COPAL® spacem is a fast-curing cement-like polymer based on poly(methylacrylate, methyl methacrylate) with no added antibiotics.

COPAL® spacem contains calcium carbonate as an x-ray contrast agent. The colouring agent chlorophyll (E141) is added to improve the visibility of spacers made of **COPAL® spacem** in the surgical area.

COPAL® spacem is prepared immediately before use by mixing a polymer powder component and a liquid monomer component. This forms a pre-polymerised dough-like mass, which sets in a few minutes.

Composition

A package of **COPAL® spacem** contains two bags with powder (polymer powder) and two brown glass ampoules (monomer liquid).

Constituents of the polymer powder

Poly(methylacrylate, methyl methacrylate), calcium carbonate, benzoyl peroxide, colour E141.

Constituents of the monomer liquid

Methyl methacrylate, N,N-dimethyl-p-toluidine, hydrochinone, colour E141.

The cement powder is triple packaged. The outer, unsterile protective aluminium bag contains a polyethylene paper bag (peel-off) which is unsterile on the outside and sterile on the inside. Inside this there is another sterile polyethylene paper bag containing the cement powder. The brown glass ampoule with the sterile-filtered monomer liquid is also sterile packed in an ethylene oxide sterilised individual blister pack.

Indication

COPAL® spacem is indicated for extracorporeal manufacture of non-articulating and articulating temporary spacers for two-stage revisions with the following goals:

- deadspace management
- stabilisation of the resected area
- possibility of mobilisation

Contraindication

COPAL® spacem must not be used on the patient

- for permanent, stable anchorage of total or partial joint endoprostheses in the bone
- for suspected or confirmed hypersensitivity to constituents of the polymer
- during pregnancy and lactation

Target population

There is very little evidence with children and adolescents, therefore it is not recommended to use **COPAL® spacem**. If no other option is available, for example with present surgical trauma, the decision whether to use **COPAL® spacem** lies with the attending surgeon.

Target user group

Healthcare professionals in a clinical context.

Side effects

When using self-curing polymers (such as poly(methylacrylate, methyl methacrylate) bone cements) directly in the organism, a temporary drop in blood pressure after preparation of the implant site or immediately after implantation of the temporary spacer has been observed in rare cases. In isolated cases serious complications may be encountered, including cardiac arrest, anaphylactic shock and sudden death. This side effect is also possible with **COPAL® spacem**, but is very unlikely because the temporary spacer is manufactured extracorporally.

We recommend irrigating the implant site thoroughly with an isotonic solution before placement of the temporary spacer (pulse lavage irrigation) to prevent pulmocardiovascular complications such as lung embolism and cardiac arrest. In the event of pulmonary or cardiovascular events, the blood volume must be monitored and increased if necessary. Anaesthetic measures may be required in the case of acute respiratory problems.

The following undesirable effects have been observed in association with the use of self-curing poly(methyl methacrylate) polymers (poly(methyl methacrylate) bone cements): thrombophlebitis, haemorrhage. Other side effects that have been observed are: myocardial infarction, short-term cardiac dysrhythmia, cerebral stroke. The complications that may be encountered in any surgical procedure are also possible.

The following spacer-related effects have been observed: spacer migration or dislocation, spacer fracture, bone fracture.

Precautions for use

Use by operating theatre staff

Before working with **COPAL® spacem** users should be thoroughly familiar with its properties, handling and application. Users are recommended to practice the complete mixing and handling procedure of **COPAL® spacem** before using it for the first time. Thorough knowledge is also necessary if mixing systems and syringes are used for processing the self-curing polymer.

The liquid monomer is highly volatile and flammable. Ignition of monomer fumes caused by use of electrocautery devices in surgical sites near freshly implanted bone cements has been reported. The monomer is also a strong lipid solvent and should not come into direct contact with the body. Gloves must be worn when working with the monomer or the mixed, self-curing **COPAL® spacem** polymer to ensure the required protection against penetration of the monomer methyl methacrylate into the skin. PVP gloves (three layers of polyethylene, ethylene vinyl alcohol copolymer, polyethylene) and Viton®/butyl gloves have been proven over an extended period to give adequate protection. For the sake of safety, we recommend using two pairs of gloves at a time, e.g. polyethylene surgical gloves over an inner pair of standard latex surgical gloves. The use of latex or polystyrolbutadiene gloves alone is not sufficient. Please ask your suppliers which gloves are suitable for this application.

The monomer vapours may irritate the eyes and respiratory passages and may also cause liver damage. Skin irritations caused by contact with the monomer have been described. Manufacturers of soft contact lenses recommend removing lenses of this type in the presence of harmful or irritating vapours. Soft contact lenses are permeable to liquids and gases and therefore should not be worn in the operating theatre if methyl methacrylate is being used.

Mechanical stability

The size and shape of the temporary spacers must be planned to ensure that the spacers are strong enough for the purpose. It is important to avoid thin sections, notches, sharp corners and sharp edges on the spacers. We recommend integrated metal reinforcements in the spacers. The metal reinforcements must be large enough to provide strength. Make sure that if metal reinforcements are used, the metal surface is clean and grease-free to ensure optimum bonding of the pre-polymerised mass to the metal reinforcement to provide adequate strength. It is also important to ensure that only standard surgical metals and alloys that are biocompatible are used. **COPAL® spacem** is ready for the manufacture of spacers. When manufacturing temporary spacers, it is important to ensure when mixing the two components, polymer powder and monomer liquid, that the mixing ratio of the two constituents is not changed (see the instructions below in the second paragraph under "Dosage"). If additional solid substances must be mixed with the **COPAL® spacem** powder on medical orders, a quantity of 4 g for the complete temporary spacer must not be exceeded with reference to the specific requirements of the patient. If this requirement is ignored, the physical and mechanical properties of the temporary spacer made of **COPAL® spacem** may be significantly compromised. The chemical composition of **COPAL® spacem** also means that compared to spacers made of conventional bone cement an increased volume of water-soluble substances may be released from the polymer matrix into the organism.

Application on patients

Blood pressure, pulse and respiration must be carefully monitored during and immediately after implantation of temporary spacers made of **COPAL® spacem**. Any significant changes in these vital signs must be corrected by appropriate means without delay. When using temporary spacers made of **COPAL® spacem** the prepared bone must be thoroughly irrigated and aspirated shortly before placement of the spacers. The curing process of spacers made of **COPAL® spacem** generates heat, the temperature of which may increase above the value for tissue compatibility. Therefore, temporary spacers made of **COPAL® spacem** should not be implanted into the bone until they are fully cured and have cooled down to room temperature. The use of **COPAL® spacem** is not recommended for children and adolescents, because there is no clinical experience with these groups.

Incompatibilities

No additional liquid solutions apart from the included monomer constituent may be mixed with the self-curing **COPAL® spacem** polymer. Any substances that could impair the curing process of the self-curing polymer must not be added.

Dosage

After mixing the polymer powder with the monomer liquid, a fast-curing pre-polymerised doughy mass is formed. The temporary spacer can be manufactured by manual forming or by placing the pre-polymerised mass into suitable moulds (such as plastic or metal moulds). **COPAL® spacem** has a green colour to ensure that the resulting temporary spacer remains clearly visible in the surgical site.

One dose is prepared by mixing the complete contents of a bag of polymer powder with the complete contents of an ampoule of monomer liquid, even if the complete amount of the resulting mass is not required for the manufacture of a temporary spacer. The mixing ratio of polymer powder and monomer liquid must not be changed when mixing the two constituents. The amount of pre-polymerised mass used depends on the size and volume of the required temporary spacer.

At least one additional dose of **COPAL® spacem** should be available before starting to manufacture the temporary spacer to ensure that the whole spacer can be manufactured from one pre-polymerised mass. One or more additional mixing systems should also be available. If larger quantities of **COPAL® spacem** are mixed (e.g. 5 packages for one hip spacer), it may be necessary to mix the polymer in several mixing processes, depending on the maximum capacity of the mixing system. In this case it is important to ensure that the additional pre-polymerised mass is always processed before the setting phase of the initially applied pre-polymerised mass begins.


Preparation and Mixing

Preparation

We recommend chilling the constituents to 4–7 °C for at least 24 hours beforehand. The constituents should not be removed from the refrigerator and added to the mixing container until the mixing process is about to start. Chilling extends the forming and setting times, because the initial viscosity for **COPAL® spacem** is slightly reduced compared to an unchilled **COPAL® spacem**.

Before opening the non-sterile aluminium bag (see above), shake or tap the contents down to prevent damage to the contents when opening the bag. The polyethylene paper bag and the ampoule must be opened only under sterile conditions. The sterile constituents (inner polyethylene paper bag and glass ampoule) are passed into the sterile area unopened.

Opening under sterile conditions

The opening flaps from the bag top assist to detach the foil from the paper.

In order to grip as much of the opening flaps as possible, the side of the paper/foil should be kept between thumb, index and middle finger.

Please use the whole thumb surface to grip foil and paper side and take off each side evenly.

The outer polyethylene paper bag is opened at the place of use to allow the sterile inner polyethylene paper bag to be removed under sterile conditions. The blister package is also opened under sterile conditions and at the place of use to allow the glass ampoule to be removed while sterile.

Before opening the inner polyethylene paper bag it is shaken or tapped to shake the contents down to ensure that no powder is lost when the top edge of the bag is cut open. Do not open the ampoule over the mixing device to prevent contamination of the cement with glass fragments. The glass ampoule has a break point at the transition to the head of the ampoule and a tube to break it to make it easier to open. Grip the tube instead of the head of the ampoule and break it off together with the ampoule head. The broken off head of the ampoule remains in the tube.

Mixing the constituents

For mixing of **COPAL® spacem** all Heraeus Medical mixing systems are suitable. Should the user want to use **COPAL® spacem** in combination with other mixing systems, he is advised to check their suitability prior to first use. The constituents are mixed under sterile conditions. The mixing time is 30 seconds unless otherwise recommended. The complete contents of one bag must always be mixed with the complete contents of one ampoule. The constituents can be mixed with a vacuum mixer or manually. The mixing, modelling and setting times of **COPAL® spacem** are shown in the diagrams at the end of the instructions for use. Note that the times are only estimates, because the modelling and setting times depend on the temperature, the mixing process and the humidity. The ambient temperatures of items such as polymer powder, mixing system, table and hands all play a role in deciding the times. Higher temperatures reduce the modelling and setting times.

Vacuum mixing: The pre-polymerised mass is stirred under vacuum to reduce the number of bubbles in the temporary spacer to a minimum and to obtain optimum mechanical strength. This process requires the use of a closed airtight mixing system and fast establishment of an adequate vacuum in the mixing container (approx. 200 mbar absolute pressure). For details of the mixing technique please see the instructions for use of the mixing system. The mixing phase (30 seconds) prepares a homogenous, green, doughy mass that can be modelled by hand directly after mixing or by means of suitable moulds to produce a temporary spacer.

Hand mixing: The liquid monomer constituent is placed in the container first and then the polymer powder. During the mixing phase (30 seconds) the two constituents are mixed together by regular, even stirring. This produces a homogenous, green, doughy mass that can be manually modelled to a temporary spacer after the end of the waiting phase. If the temporary spacer is produced by means of suitable moulds, the waiting phase can be dispensed with.

Use of the polymer

To manufacture a temporary spacer with **COPAL® spacem**, the pre-polymerised mass must be modelled manually or placed in suitable moulds during the modelling phase. After the modelling phase the pre-polymerised mass cures and when cooled to room temperature it can be placed in the prepared bone cavity.

For optimum mechanical strength of the temporary spacer, the use of metal reinforcements is recommended. Note that only biocompatible metals and alloys approved for surgery can be used. The pre-polymerised mass is modelled around the metal reinforcements during the modelling phase. The pre-polymerised mass can be modelled around the metal by hand or placed in moulds that contain the metal reinforcement. If a mould is used, the metal reinforcements must be positioned to allow the pre-polymerised mass to surround them completely. The pre-polymerised mass cures after the modelling phase to form a metal-reinforced spacer.

Storage

Do not store above 25 °C.

Shelf life / sterility

The shelf life of the self-curing polymer is printed on the box, the aluminium bag and the inner bag. **COPAL® spacem** must not be used after the date printed on the box. The shelf life of the individual constituents may vary from that printed on the box because of production conditions. The contents of opened or damaged aluminium sachets or ampoule blisters must not be re-sterilized and thus must be discarded.

Disposal

Single components of the bone cement, cured solid material as well as (uncleaned) packaging material must be disposed by following the regulations of the local authorities. Dispose of the polymer component in an authorized waste facility. The liquid component should be evaporated under a well-ventilated hood or absorbed by an inert material and transferred in a suitable container for disposal.

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Heraeus Medical GmbH