

Summary of Safety and Clinical Performance

PALACOS[®] R+G

- Patient Section -

Document number : 61339
Effective date : 10.12.2024

English

Titel: SSCP PALACOS R+G

Doc.-No.: 61339

Page 2 of 17

1 Table of contents

	Page
1 Table of contents.....	2
2 Abbreviations / Explanations	3
3 General Information.....	3
3.1 Relevant information for Users/Healthcare Professionals.....	3
3.2 Relevant Information for patients	4
3.2.1 Background information	4
3.2.2 Device identification and general information.....	5
3.2.2.1 Products (device trade names) covered by this document.....	5
3.2.2.2 Manufacturer name and address.....	5
3.2.2.3 Basic UDI-DI number of the concerned product.....	5
3.2.2.4 Year of first CE-mark	5
3.2.3 Intended use of the device	5
3.2.3.1 Intended purpose.....	5
3.2.3.2 Indications and intended patient groups	5
3.2.3.3 Contraindications/ advice against treatment.....	6
3.2.3.4 Lifetime of the device.....	6
3.2.4 Device description.....	6
3.2.5 Risks and warning.....	7
3.2.6 Summary of clinical evaluation and post-market clinical follow-up.....	8
3.2.7 Possible diagnostic or therapeutic alternatives	10
References	12

List of Tables

	Page
Table 1: Composition of PALACOS® R+G.....	6
Table 2: Frequencies of side effects.....	7
Table 3: Frequencies of residual risks	8

Titel: SSCP PALACOS R+G

Doc.-No.: 61339

Page 3 of 17

2 Abbreviations / Explanations

ALBC	Antibiotic-loaded bone cement
AUS	Australia
BCIS	Bone cement implantation syndrome
BfArM	Federal Institute for Drugs and Medical Devices (Germany) [<i>Bundesinstitut für Arzneimittel und Medizinprodukte</i>]
CAN	Canada
CE	Conformité Européenne
CER	Clinical Evaluation Report
CH	Switzerland
CND	<i>Classificazione Nazionale dei Dispositivi medici</i> [National Classification of Medical Devices]
CT	Computed tomography
CS	common specifications as defined in the MDR
DIN	German standard [<i>Deutsches Institut für Normung</i>]
EMDN	European Medical Device Nomenclature
EN	European Standard [<i>Europäische Norm</i>]
EU	European Union
FDA	Food and Drug Organization (USA)
GER	Germany
IFU	Instructions for Use
ISO	International Organization for Standardization
MDD	Medical Device Directive
MDR	Medical Device Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices, amending Directive 2001/83/EC, Regulation (EC) No 178/2002 and Regulation (EC) No 1223/2009 and repealing Council Directives 90/385/EEC and 93/42/EEC and the amendments of 2017/45 (2020/561 and 2023/607)
MHRA	Medicines and Healthcare Products Regulatory Agency in the UK
MRI	Magnetic resonance imaging
N/A	Not applicable
NB	Notified Body
PMCF	Post-Market Clinical Follow-Up
PMMA	poly (methyl methacrylate)
PMS	Post-Market Surveillance
SRN	single registration number for an economic operator
SSCP	Summary of Safety and Clinical Performance
Swissmedic	Swiss Agency for Therapeutic Products (Switzerland)
TGA	Therapeutic Goods Administration (Australia)
TPLC	FDA Total Product Life Cycle (TPLC)
UDI-DI	Unique Device Identification - device identifier
UK	United Kingdom
US/USA	United States of America

3 General Information

This document applies to implantable class IIb, and class III medical devices developed by Heraeus Medical GmbH and is established to comply with the Medical Device Regulation (MDR) 2017/745 (EU) of 5th April 2017, valid from May 2021.

The Summary of Safety and Clinical Performance (SSCP) is intended to provide a summary of clinical data pertinent to the safety and clinical performance of the medical device. The SSCP is an important source of information for intended users – both healthcare professionals and if relevant for patients. It is one of several means intended to fulfil the MDR objectives, to enhance transparency and provide adequate access to information.

3.1 Relevant information for Users/Healthcare Professionals

Please refer to the English language version of this document for the contents of this section.

Titel: SSCP PALACOS R+G

Doc.-No.: 61339

Page 4 of 17

3.2 Relevant Information for patients

The following chapters provide a summary of the safety and clinical performance of the device intended for patients.

This Summary of Safety and Clinical Performance (SSCP) provides public access to an updated summary of the main aspects of the safety and clinical performance of the device. The information presented below addresses patients or lay persons. The first part of the document shows a more extensive summary of safety and clinical performance prepared for healthcare professionals.

The SSCP does not provide general advice on the treatment of a medical condition. Please contact your doctor/surgeon in case you have questions about your medical condition or about the use of the device in your situation. This SSCP does not replace an Implant Card or the Instruction for Use (IFU) to provide information on the safe use of the device.

3.2.1 Background information

PALACOS® R+G is a bone cement. It is based on a biologically safe material called poly (methyl methacrylate) (PMMA). This material has a long history of safe use in humans. PALACOS® R+G was the basis for development of the bone cements PALACOS® MV+G, PALACOS® LV+G, and PALACOS® fast R+G. All four bone cements together comprise the product family PALACOS® +G bone cements.

PALACOS® +G bone cements are used in adults such as elderly patients with degenerative joint disease. Osteoarthritis is an example for such a joint disease. Osteoarthritis is the most common form of arthritis and affects millions of people worldwide. It occurs when the protective cartilage that cushions the ends of the bones wears down over time. Patients with trauma after severe accidents with several fractures in a bone can also be considered for treatment with bone cements. The bone cement is used to anchor total or partial joint endoprostheses. It attaches endoprostheses firmly and stably to the bone. Endoprostheses are medical devices used to replace parts of the inside of your body. Hip, knee, or shoulder joints can be replaced by an endoprostheses, for example.

Arthroplasty is a surgical procedure to restore the function of a joint. Primary arthroplasty refers to the first joint replacement. Revision arthroplasty refers to follow-up surgery on the same joint. In total joint replacement parts of a joint are removed and replaced by an implant, the endoprosthesis. In partial joint replacement artificial surfaces replace only the movable surfaces of a joint. The healthy parts of the joint stay intact.

Bone cements can also treat cases of bone loss. For example, after severe accidents with multiple fractures in a bone. The name of this surgical technique is reconstruction of bone. It restores bone continuity mainly in patients suffering from tumor of the bone or in trauma.

Your doctor/surgeon applies the bone cement during surgery. The instructions for use give directions.

Your doctor/surgeon takes care of the following aspects during your surgery:

- The bone cement is applied to your carefully cleaned, aspirated, and dried bone.
- Your prosthesis is put in place and held until the bone cement has set completely.
- During and immediately after the bone cement is applied, your doctor/surgeon will monitor your blood pressure, pulse, and breathing carefully. This ensures early detection and treatment of adverse events such as low blood pressure and cardiac arrest. Drops in blood pressure have occurred remotely and shortly after application of bone cement. However, consequences such as cardiac arrest are only reported in very few cases.

It is safe to have magnetic resonance tests (MRI) with PALACOS® +G bone cements. But the composition of the prosthesis you receive together with the bone cement may affect your ability to have magnetic resonance tests. You will receive an implant card for the bone cement that was used. Additionally, you will receive an implant card for the prosthesis. Please keep these documents and provide them in future examinations (e.g., X-ray, CT scan, MRI).

Titel: SSCP PALACOS R+G

Doc.-No.: 61339

Page 5 of 17

3.2.2 Device identification and general information

3.2.2.1 Products (device trade names) covered by this document

- PALACOS® R+G

3.2.2.2 Manufacturer name and address

Heraeus Medical GmbH
Philipp-Reis-Str. 8/13
61273 Wehrheim
Germany

3.2.2.3 Basic UDI-DI number of the concerned product

The unique device identification (UDI) consists of a series of numbers with letters. It allows the unmistakable identification of a specific medical device on the market. An UDI device identifier (UDI-DI) is specific to a device, connecting the product to the information on the EUDAMED database.

The following UDI-DI numbers are assigned to the different products:

Product	Basic UDI-DI
PALACOS® R+G	4260102130101010001AJ

3.2.2.4 Year of first CE-mark

Before a medical device is introduced on the market in the European Union, it needs to show that the product fulfills the requirements. The so-called CE-certification documents the fulfilment, and the CE-mark is placed on the product. The legal requirements for medical devices have changed in May 2021. Then, the Medical Device Regulation (MDR) replaced the Medical Device Directive (MDD).

The following table contains the detailed information about the different products. The table lists the year of the first CE-mark under MDR and under MDD.

Product	Year of first CE-mark under MDR	Year of first CE-mark under MDD
PALACOS® R+G	2022	2005

3.2.3 Intended use of the device

3.2.3.1 Intended purpose

PALACOS® +G bone cements are intended for stable anchoring of total or partial joint replacements (endoprostheses) in living bone as well as for reconstruction of bone.

3.2.3.2 Indications and intended patient groups

PALACOS® +G bone cements are indicated for surgical treatment such as

- anchoring of endoprosthesis in primary and revision arthroplasty procedures of
 - hip
 - knee
 - ankle
 - shoulder
 - elbow
- reconstruction of bone via induced membrane technique after tumor surgery and / or trauma

Titel: SSCP PALACOS R+G

Doc.-No.: 61339

Page 6 of 17

These treatments are typically conducted in adults, predominantly elderly patients with osteoarthritis and patients with trauma.

3.2.3.3 Contraindications/ advice against treatment

PALACOS® +G bone cements must not be used in the following cases:

- known or suspected intolerance against parts of the bone cement or against the antibiotic gentamicin
- infection at the site of the body where the surgery is planned
- patients with impaired kidney function
- reconstruction of skull bone defects
- spinal surgery
- children

3.2.3.4 Lifetime of the device

There is no general factor influencing the lifetime of the PALACOS® +G bone cements. The general provisions for the prosthesis they anchor also apply to the bone cements. The actual lifetime of these bone cements can be influenced by factors such as your medical situation and your lifestyle.

3.2.4 Device description

PALACOS® +G bone cements are based on a biologically safe material called poly (methyl methacrylate) (PMMA) which has a long history of safe use in humans.

Composition

The cement consists of two main components, a powder, and a liquid. The table below shows the composition of the components. Mixing of the components starts a chemical reaction. This so-called polymerization forms a soft dough. The dough becomes more and more solid over time. Your surgeon determines the right time for the application of the dough to the bone. There it hardens completely. In addition, the cement contains an antibiotic (gentamicin). Your treating surgeon chose the antibiotic to prevent an infection.

PALACOS® R+G contains:

Table 1: Composition of PALACOS® R+G

Constituents	PALACOS® R+G
Powder:	
PMMA copolymer <i>Polymer (powder component)</i>	82 %
zirconium dioxide <i>X-ray contrast medium (enabling visualization with X-ray, CT or MRI)</i>	15%
benzoyl peroxide <i>Chemical component initiating the polymerization reaction</i>	1 %
gentamicin sulfate <i>(Antibiotic)</i>	2 %
Liquid:	
methyl methacrylate <i>Monomer (liquid component)</i>	98 %
N, N-dimethyl-p-toluidine <i>Chemical component accelerating the polymerization reaction</i>	2 %

The data is rounded

Other constituents:

- Powder: chlorophyll-copper-complex (E141) (Food colorant. Improving visibility of the bone cement in the surgical field)
- Liquid: chlorophyll-copper-complex (E141), hydroquinone (chemical component stabilizing the chemical reaction)

Titel: SSCP PALACOS R+G

Doc.-No.: 61339

Page 7 of 17

Traces of histamine may be present in these bone cements. But no manufacturing residuals that could pose a risk to you have been found. Be aware that the composition table shows the constituents before mixture of the bone cement components. The methyl methacrylate is completely used up during setting and forms the hardened bone cement. PALACOS® R+G bone cement is intended for single-use and are supplied sterile.

3.2.5 Risks and warning

Contact your doctor/surgeon if you believe that you are experiencing side effects. This applies for side effects related to the device or its use and if you are concerned about risks. This document does not replace a consultation with your doctor/surgeon if needed.

Side effects are events that are known when using the device. They can be caused by the device.

Residual risks are risks which cannot be controlled by the device manufacturer. They are mostly related to the surgical procedure in general or to the handling of the user.

Adverse events are events that can occur in a clinical investigation. They have a negative impact mostly on the patient. No causal relationship with the device must be present.

Heraeus Medical GmbH has a risk management process that complies with harmonized risk management guidelines. It ensures that the benefits of using the medical device are greater than any potential risks.

Side effects and residual risks of the device can occur with different frequencies. As an example, if a side effect occurs in less than 1% of cases (< 1%), the side effect will occur in less than 1 in 100 surgeries.

Side effects

Frequencies are taken from literature, < 0.0001% can include potential risks not reported in literature so far.

Table 2: Frequencies of side effects

Frequency	Side effect
Immune system	
4.79%*	• hypersensitivity / allergic reaction and local reaction which may include inflammation, induration, erythema, pruritus, or pain
< 0.0001%	• anaphylactic shock
Kidney and Urinary Tract	
< 0.0001%	• renal impairment
Musculoskeletal System	
36.65%*	• ossification
22.78%*	• osteolysis due to bone cement fragments
Skin and Subcutaneous Tissue	
< 0.0001%	• rash
< 0.0001%	• urticaria

*cases reported to Heraeus Medical GmbH showed frequencies lower than in literature

Residual risks

Residual risks listed below are procedure related risks which are beyond the control of the manufacturer because they are procedure or user related. Frequencies are taken from literature, < 0.0001% can include potential risks not reported in literature so far.

Titel: SSCP PALACOS R+G

Doc.-No.: 61339

Page 8 of 17

Table 3: Frequencies of residual risks

Frequency	Residual Risk
Vascular System, Heart, Respiratory System, Blood and Lymphatic System, Nervous System	
Bone Cement Implantation Syndrome (BCIS): Insertion of bone cement may produce a high medullary pressure that forces bone marrow constituents into the venous vascular system, resulting in fat and marrow emboli. To avoid BCIS, it is recommended that the implantation site is cleaned thoroughly with pulsatile, high pressure, high-volume lavage using an isotonic solution and dried before the bone cement is introduced. The bone cement should be applied retrogradely under sustained low pressure into the medullary canal. Subsequently, the prosthesis should be introduced slowly into the cemented medullary canal. In case of pulmonary or cardiovascular events, it is necessary to monitor blood volume and possibly increase it. In case of acute respiratory failure, anesthesiologic measures should be taken. In general, adverse reactions of BCIS might include low blood pressure/hypotension, hypoxia, bradycardia, tachycardia, pulmonary hypertension, thrombosis, embolism, pulmonary embolism, myocardial infarction, cerebrovascular accident, respiratory arrest, and cardiac arrest.	
25.48%*	BCIS grade 1 moderate hypoxia (SpO2 < 94%) or hypotension [fall in systolic blood pressure > 20%]
8.17%*	BCIS grade 2 severe hypoxia (SpO2 < 88%) or hypotension [fall in systolic blood pressure > 40%] or unexpected loss of consciousness.
3.37%*	BCIS grade 3 cardiovascular collapse, requiring CPR
Nervous System	
< 0.0001%	• numbness
Blood and Lymphatic System	
< 0.0001%	• hypovolemia
Musculoskeletal System	
15.22%*	• aseptic loosening
< 0.0001%	• unequal limb length
0.27%*	• loss of range of motion
< 0.0001%	• ambulation difficulties
Infection	
3.58%*	• bacterial infection including cellulitis, and / or osteomyelitis
Generalized Disorders	
2,56%	• inflammation
< 0.0001%	• swelling / edema
0.94%*	• fibrosis
< 0.0001%	• heat necrosis

*cases reported to Heraeus Medical GmbH showed frequencies lower than in literature

Please contact your healthcare professional if you have any questions.

Reporting of side effects, residual risks, or adverse events

If you experience any of these side effects or residual risks, or if you notice any adverse events not listed in this document, contact your doctor/surgeon immediately. You can also contact Heraeus Medical GmbH directly using the following email-address: hm.vigilance.medical@heraeus.com

3.2.6 Summary of clinical evaluation and post-market clinical follow-up

PALACOS® R+G was the first bone cement with an antibiotic introduced in 1972. All further products of the PALACOS® +G product family like PALACOS® MV+G, PALACOS® LV+G and PALACOS® fast R+G are based on PALACOS® R+G. They own few modifications in terms of product characteristics. These products were developed and placed on the market between 1975 (PALACOS® LV+G), 1998 (PALACOS® MV+G) and 2013 (PALACOS® fast R+G). PALACOS® +G bone cements treated a total of about 30 million patients worldwide so far.

Titel: SSCP PALACOS R+G

Doc.-No.: 61339

Page 9 of 17

The product range of PALACOS® +G bone cements can be considered as state-of-the-art in the field of stable anchoring of joint endoprostheses as well as for reconstruction of bone.

The manufacturer performs the analysis of any clinical data regularly. Sources can be endoprosthesis registries and scientific publications, for example. These activities are called post-market clinical follow-up measures. They allow the continuous proof of the benefit/risk ratio of the medical device. Registries are databases which collect long-term results after application of products in patients. These databases can be initiated by governmental authorities, medical societies, or manufacturers. In most cases they collect data from hospitals or private practices on a regional or national level.

The following clinical benefits and outcome parameters relate to the use of the bone cements:

- Stable fixation of the endoprosthesis with a low risk of revision surgery. This is evaluated based on long-term data from regional or national registries.
- Improvement of impaired body function with a high patient satisfaction. This is evaluated based on quality-of-life data from registries.
- Relief of symptoms related to the surgical procedure with high patient success. This is evaluated based on quality-of-life data from registries.
- Application of bone cements in combination with an antibiotic with a low risk of infection. This is evaluated based on revisions that are caused by infections, compared to the overall number of revisions (based on data from registries).
- Local use of an antibiotic within the bone cement can result in a reduced risk for side effects compared to oral or intravenous administration of the antibiotic. This is evaluated based on complaints reported to manufacturer, evaluation of databases and data regarding the development of the medical device.
- The reconstruction of bone via induced membrane technique can result in the preservation of function of the limb or of the limb itself. This is evaluated by determination of the union of bone defects after tumor surgery and / or trauma.

The above-mentioned clinical benefits and clinical outcome parameters are important to decide on the benefit/risk ratio of PALACOS® R+G. The manufacturer evaluates the achievement of these clinical benefits.

The analysis revealed that PALACOS® R+G performed as expected in all aspects of the above-listed outcome parameters:

- Stable fixation was analyzed by the rate at which operations needed to be repeated (revision rate). The rate was in a range comparable with the current state-of-the-art. For example, the revision rate of PALACOS® R+G was reported to be 6.2% for primary hip and 5.6% for primary knee, which is comparable to benchmark standards (range for hip: 6.2% - 9.5% after 15 years; range for knee: 5.6% - 16.2% after 15 years).
- Impaired body function was evaluated through questionnaires. In these, patients have reported on how much they are impacted in their daily activities. In all cases, PALACOS® R+G were comparable to [the](#) current state-of-the-art.
- Relief of symptoms was evaluated through questionnaires. In these, patients have reported on how much better their joint was after the surgery. In all cases, PALACOS® R+G were comparable to [the](#) current state-of-the-art.
- The number of re-operations because of an infection at the site of surgery was comparable to the current state-of-the-art in patients who underwent their first operation with PALACOS® R+G and for revision surgeries.

Titel: SSCP PALACOS R+G

Doc.-No.: 61339

Page 10 of 17

- PALACOS® R+G contains an antibiotic that can also be given directly into the veins. From this it is known that too high amounts can cause severe side effects. In a clinical study, it was measured how high up the blood concentrations of antibiotics released from the bone cement would go after an operation with PALACOS® R+G. The result was that the values remained far below the levels which can lead to severe side effects.
- Reconstruction of bone via induced membrane technique was analyzed by the rate of successful bone union after two-stage reconstruction. The performance of PALACOS® R+G was comparable to the state-of-the-art.

Additionally, the scientific literature for PALACOS® R+G was thoroughly evaluated, and 80 scientific publications were identified and analyzed. It can be summarized that all data show favorable clinical results for PALACOS® R+G.

In conclusion, the success rates of the clinical benefits were comparable to or better than the current state-of-the-art.

Therefore, the manufacturer confirms that the benefits outweigh the risks for the indications of PALACOS® R+G:

- anchoring of endoprosthesis in primary and revision arthroplasty procedures of
 - hip
 - knee
 - ankle
 - shoulder
 - elbow
- reconstruction of bone via induced membrane technique after tumor surgery and / or trauma

The following activities are planned to ensure safety and performance of PALACOS® +G bone cements:

- Device Registry Analysis, to monitor the safety and performance of PALACOS® +G bone cements
- Screening of Scientific Literature, to monitor the safety and performance of PALACOS® +G bone cements
- Authority Databases (adverse events and recalls), to monitor the safety of PALACOS® +G bone cements

The same activities are performed for similar products, in order to detect potential safety or performance issues early. The results will be summarized in reports. These activities will be conducted on an annual basis in connection with the continuous updates of the clinical evaluations.

3.2.7 Possible diagnostic or therapeutic alternatives

General information

Contact your doctor/surgeon when you consider alternative treatments. Depending on your individual situation, two treatment approaches are possible. On the one hand conservative treatment such as physiotherapy or pain medication without surgery is possible. On the other hand, surgical treatment such as joint surgery like hip replacement surgery could be reasonable. Choice of treatment depends on your specific condition and your doctor's opinion.

Joint surgery

If possible, your doctor/surgeon will try to treat defective joints by other means. If all other treatment options fail, a reconstructive joint surgery may be necessary. This means, the complete joint or only parts of the joint are replaced by an endoprosthesis. Joint surgeries and endoprosthesis revision operation as well as the use of PMMA bone cements are very well-established procedures in joint replacement surgery.

PMMA is widely and successfully used for the fixation of various endoprostheses since decades. At present, PMMA is still the most commonly used fixation material in primary joint surgeries. Uncemented procedures have also been used in primary joint surgeries. However, current data do not allow to determine if cementless or cemented generally perform better in joint surgeries. The advantage of the cemented procedures using PMMA is the long-term experience with this material. Also, the majority of orthopedic surgeons is familiar with the use of PMMA. Furthermore, bone cement can apply local antibiotics. This allows for infection prevention in patients at risk for infection. In addition, bone cements generally spread the force of movement evenly into the bone. Especially in

Titel: SSCP PALACOS R+G

Doc.-No.: 61339

Page 11 of 17

patients with poor bone substance this is an advantage. Your doctor/surgeon will decide on the procedure that fits to your specific clinical condition best.

There is no other treatment option than a surgery in patients with suspected or confirmed infection of the implanted device (so-called prosthetic joint infections). Such a revision surgery can be either a one-stage or a two-stage surgery. A so-called one-stage surgery takes place in a single surgical step. The surgeon removes the infected prosthesis and bone cement, cleans the surgical site thoroughly, and places a new prosthesis. A so-called two-stage approach consists of two separate surgeries. During the first surgery, the surgeon removes the infected prosthesis and bone cement, cleans the surgical site thoroughly, and places a provisional spacer. This ensures proper treatment of the infection. The spacer also provides a limited range of motion during the time until the second operation. After the infection is cured, the second surgery takes place. The surgeon removes the provisional spacer and places a new permanent prosthesis. The attending surgeon will choose the appropriate surgical approach according to the patient's situation.

Reconstruction of bone

Oncological treatment or trauma may lead to bone loss. PMMA bone cement is able to fill certain bone defects depending on the depth and surface of the defect. The method "Induced Membrane Technique" can support new bone growth in an area where part of a bone had to be removed due to cancer or was lost due to trauma. For this approach bone cement is only placed between the ends of a defect for a short period of time. The bone cement is not fixed to the bone.

For larger defects further therapy options have to be considered. Therapy options like human tissues from donors, metal implants or custom-made prostheses are available. The attending surgeon will choose the appropriate surgical approach according to the patient's situation.

Titel: SSCP PALACOS R+G

Doc.-No.: 61339

Page 12 of 17

References

1. Wahlig H, Dingeldein E, Buchholz HW, Buchholz M, Bachmann F. Pharmacokinetic study of gentamicin-loaded cement in total hip replacements. Comparative effects of varying dosage. *J Bone Joint Surg Br.* 1984;66(2):175-9.
2. Trela-Larsen L, Sayers A, Blom AW, Webb JCJ, Whitehouse MR. The association between cement type and the subsequent risk of revision surgery in primary total hip replacement. *Acta Orthop.* 2018;89(1):40-6.
3. van der Voort P, Valstar ER, Kaptein BL, Fiocco M, van der Heide HJL, Nelissen RGHH. Comparison of femoral component migration between Refobacin bone cement R and Palacos R + G in cemented total hip arthroplasty: A randomised controlled roentgen stereophotogrammetric analysis and clinical study. *Bone Joint J.* 2016;98-B(10):1333 – 41.
4. Jørgensen PB, Lamm M, Søballe K, Stilling M. Equivalent hip stem fixation by Hi-Fatigue G and Palacos R + G bone cement: a randomized radiostereometric controlled trial of 52 patients with 2 years' follow-up. *Acta Orthop.* 2019;90(3):237 – 42.
5. Tabori-Jensen S, Mosegaard SB, Hansen TB, Stilling M. Inferior stabilization of cementless compared with cemented dual-mobility cups in elderly osteoarthritis patients: a randomized controlled radiostereometry study on 60 patients with 2 years' follow-up. *Acta Orthop.* 2020;91(3):246-53.
6. Turnbull GS, Marshall C, Nicholson JA, MacDonald DJ, Clement ND, Breusch SJ. The Olympia anatomic polished cemented stem is associated with a high survivorship, excellent hip-specific functional outcome, and high satisfaction levels: follow-up of 239 consecutive patients beyond 15 years. *Arch Orthop Trauma Surg.* 2021.
7. Siepen W, Zwicky L, Stoffel KK, Ilchmann T, Clauss M. Prospective two-year subsidence analysis of 100 cemented polished straight stems - a short-term clinical and radiological observation. *BMC Musculoskeletal Disord.* 2016;17(1):395.
8. Jain S, Magra M, Dube B, Veysi VT, Whitwell GS, Aderinto JB, Emerton ME, Stone MH, Pandit HG. Reverse hybrid total hip arthroplasty. *Bone Joint J.* 2018;100-B(8):1010 – 1017.
9. Lindalen E, Dahl J, Nordsletten L, Snorrason F, Høvik, Röhrli S. Reverse hybrid and cemented hip replacement compared using radiostereometry and dual-energy X-ray absorptiometry: 43 hips followed for 2 years in a prospective trial. *Acta Orthop.* 2012;83(6):592 – 8.
10. Sprowson AP, Jensen C, Chambers S, Parsons NR, Aradhyula NM, Carluke I, Inman D, Reed MR. The use of high-dose dual-impregnated antibiotic-laden cement with hemiarthroplasty for the treatment of a fracture of the hip: The Fractured Hip Infection trial. *Bone Joint J.* 2016;98-B(11):1534 – 41.
11. Tyas B, Marsh M, Oswald T, Refaie R, Molyneux C, Reed M. Antibiotic resistance profiles of deep surgical site infections in hip hemiarthroplasty; comparing low dose single antibiotic versus high dose dual antibiotic impregnated cement. *J Bone Jt Infect.* 2018;3(3):123 – 129.
12. Unger AC, Dirksen B, Renken FG, Wilde E, Willkomm M, Schulz AP. Treatment of femoral neck fracture with a minimal invasive surgical approach for hemiarthroplasty - clinical and radiological results in 180 geriatric patients. *Open Orthop J.* 2014;8:225 – 31.
13. Sullivan NPT, Hughes AW, Halliday RL, Ward AL, Chesser TJS. Early complications following cemented modular hip hemiarthroplasty. *Open Orthop J.* 2015;9:15 – 9.
14. Kiran M, Johnston LR, Sripada S, Mcleod GG, Jariwala AC. Cemented total hip replacement in patients under 55 years. *Acta Orthop.* 2018;89(2):152 – 5.

Titel: SSCP PALACOS R+G

Doc.-No.: 61339

Page 13 of 17

15. Sanz-Ruiz P, Matas-Diez JA, Sanchez-Somolinos M, Villanueva-Martinez M, Vaquero-Martín J. Is the commercial Antibiotic-loaded bone cement useful in prophylaxis and cost saving after knee and hip joint arthroplasty? The transatlantic paradox. *J Arthroplasty*. 2017;32(4):1095 – 9.
16. Pallaver A, Zwicky L, Bolliger L, Bösebeck H, Manzoni I, Schädelin S, Ochsner PE, Clauss M. Long-term results of revision total hip arthroplasty with a cemented femoral component. *Arch Orthop Trauma Surg*. 2018;138(11):1609 – 16.
17. Born P, Manzoni I, Ilchmann T, Clauss M. Is cemented revision total hip arthroplasty a reasonable treatment option in an elderly population. *Orthop Rev (Pavia)*. 2019;11(3):8263.
18. Ebied AM, Ebied AA, Marei S, Smith E. Enhancing biology and providing structural support for acetabular reconstruction in single-stage revision for infection. *J Orthop Traumatol*. 2019;20(1):23.
19. Ilchmann T, Zimmerli W, Ochsner PE, Kessler B, Zwicky L, Graber P, Clauss M. One-stage revision of infected hip arthroplasty: outcome of 39 consecutive hips. *Int Orthop*. 2016;40(5):913-8.
20. Brüggemann A, Mallmin H, Hailer NP. Do dual-mobility cups cemented into porous tantalum shells reduce the risk of dislocation after revision surgery. *Acta Orthop*. 2018;89(2):156 – 62.
21. Mueller Greber P, Manzoni I, Ochsner PE, Ilchmann T, Zwicky L, Clauss M. Excellent long-term results of the Müller acetabular reinforcement ring in primary cup revision. *Acta Orthop*. 2017;88(6):619 – 26.
22. Tootsi K, Heesen V, Lohrengel M, Enz AE, Illiger S, Mittelmeier W, Lohmann CH. The use of antibiotic-loaded bone cement does not increase antibiotic resistance after primary total joint arthroplasty. *Knee Surg Sports Traumatol Arthrosc*. 2021.
23. Birkeland Ø, Espehaug B, Havelin LI, Furnes O. Bone cement product and failure in total knee arthroplasty. *Acta Orthop*. 2017;88(1):75 – 81.
24. Bendich I, Zhang N, Barry JJ, Ward DT, Whooley MA, Kuo AC. Antibiotic-laden bone cement use and revision risk after primary total knee arthroplasty in U.S. veterans. *J Bone Joint Surg Am*. 2020;102(22):1939 – 47.
25. Hauer G, Hörlesberger N, Klim S, Bernhardt GA, Leitner L, Glehr M, Leithner A, Sadoghi P. Mid-term results show no significant difference in postoperative clinical outcome, pain and range of motion between a well-established total knee arthroplasty design and its successor: a prospective, randomized, controlled trial. *Knee Surg Sports Traumatol Arthrosc*. 2021;29(3):827-831.
26. Fuchs A, Häussermann P, Hömig D, Ochs BG, Klopfer T, Müller CA, Helwig P, Konstantinidis L. 10-year follow-up of the Columbus knee prostheses system in a prospective multicenter study. *Arch Orthop Trauma Surg*. 2021.
27. Rassir R, Nolte PA, Lugt JCT, Nelissen RGHH, Sierveit IN, Verra WC. No differences in cost-effectiveness and short-term functional outcomes between cemented and uncemented total knee arthroplasty. *BMC Musculoskelet Disord*. 2020;21(1):448.
28. Sanz-Ruiz P, Villanueva-Martínez M, Matas-Diez JA, Vaquero-Martín J. Revision TKA with a condylar constrained prosthesis using metaphyseal and surface cementation: a minimum 6-year follow-up analysis. *BMC Musculoskelet Disord*. 2015;16:39.
29. Baur J, Zwicky L, Hirschmann MT, Ilchmann T, Clauss M. Metal backed fixed-bearing unicondylar knee arthroplasties using minimal invasive surgery: a promising outcome analysis of 132 cases. *BMC Musculoskelet Disord*. 2015;16:177.
30. Böhler C, Kolbitsch P, Schuh R, Lass R, Kubista B, Giurea A. Midterm results of a new rotating hinge knee implant: A 5-year follow-up. *Biomed Res Int*. 2017;2017:7532745.

Titel: SSCP PALACOS R+G

Doc.-No.: 61339

Page 14 of 17

31. Giurea A, Neuhaus H-J, Miehle R, Schuh R, Lass R, Kubista B, Windhager R. Early results of a new rotating hinge knee implant. *Biomed Res Int.* 2014;2014:948520.
32. Staats K, Wannmacher T, Weihs V, Koller U, Kubista B, Windhager R. Modern cemented total knee arthroplasty design shows a higher incidence of radiolucent lines compared to its predecessor. *Knee Surg Sports Traumatol Arthrosc.* 2019;27(4):1148 – 55.
33. Petursson G, Fenstad AM, Gøthesen, Haugan K, Dyrhovden GS, Hallan G, Röhrl SM, Aamodt A, Nilsson KG, Furnes O. Similar migration in computer-assisted and conventional total knee arthroplasty. *Acta Orthop.* 2017;88(2):166 – 72.
34. Shah N, Iqbal HJ, Brookes-Fazakerley S, Sinopidis C. Shoulder hemiarthroplasty for the treatment of three and four part fractures of the proximal humerus using Comprehensive® Fracture stem. *Int Orthop.* 2011;35(6):861 – 7.
35. Ernstbrunner L, Hingsammer A, Imam MA, Sutter R, Brand B, Meyer DC, Wieser K. Long-term results of total elbow arthroplasty in patients with hemophilia. *J Shoulder Elbow Surg.* 2018;27(1):126-32.
36. Wang X, Luo F, Huang K, Xie Z. Induced membrane technique for the treatment of bone defects due to post-traumatic osteomyelitis. *Bone Joint Res.* 2016;5(3):101-5.
37. Marais LC, Ferreira N. Bone transport through an induced membrane in the management of tibial bone defects resulting from chronic osteomyelitis. *Strategies Trauma Limb Reconstr.* 2015;10(1):27 – 33.
38. Scholz AO, Gehrman S, Glombitza M, Kaufmann RA, Bostelmann R, Flohe S, Windolf J. Reconstruction of septic diaphyseal bone defects with the induced membrane technique. *Injury.* 2015;46 Suppl 4:S121 – 4.
39. Helbig L, Bechberger M, Aldeeri R, Ivanova A, Haubruck P, Miska M, Schmidmaier G, Omlor GW. Initial peri- and postoperative antibiotic treatment of infected nonunions: results from 212 consecutive patients after mean follow-up of 34 months. *Ther Clin Risk Manag.* 2018;14:59 – 67.
40. Omlor GW, Lohnherr V, Lange J, Gantz S, Merle C, Fellenberg J, Raiss P, Lehner B. Enchondromas and atypical cartilaginous tumors at the proximal humerus treated with intralesional resection and bone cement filling with or without osteosynthesis: retrospective analysis of 42 cases with 6 years mean follow-up. *World J Surg Oncol.* 2018;16(1):139.
41. Kasha S, Rathore SS, Kumar H. Antibiotic cement spacer and induced membrane bone grafting in open fractures with bone loss: a case series. *Indian J Orthop.* 2019;53(2):237 – 45.
42. Tetsworth K, Woloszyk A, Glatt V. 3D printed titanium cages combined with the Masquelet technique for the reconstruction of segmental femoral defects: Preliminary clinical results and molecular analysis of the biological activity of human-induced membranes. *OTA Int.* 2019;2(1):016.
43. de Jong L, Klem T, Kuijper TM, Roukema GR. Factors affecting the rate of surgical site infection in patients after hemiarthroplasty of the hip following a fracture of the neck of the femur. *Bone Joint J.* 2017;99-B(8):1088-94.
44. Kienapfel H, Hildebrand R, Neumann T, Specht R, Koller M, Celik I, Mueller HH, Griss P, Klose KJ, Georg C. The effect of Palamed G bone cement on early migration of tibial components in total knee arthroplasty. *Inflamm Res.* 2004;53 Suppl 2:S159-63.
45. Yen SH, Chen JH, Lu YD, Wang JW. Perioperative complications of total knee arthroplasty in dialysis patients. *J Arthroplasty.* 2018;33(3):872-7.
46. Mjoberg B. Loosening of the cemented hip prosthesis. The importance of heat injury. *Acta Orthop Scand Suppl.* 1986;221:1-40.

Titel: SSCP PALACOS R+G

Doc.-No.: 61339

Page 15 of 17

47. Mjoberg B, Franzen H, Selvik G. Early detection of prosthetic-hip loosening. Comparison of low- and high-viscosity bone cement. *Acta Orthop Scand.* 1990;61(3):273-4.
48. Mjoberg B, Rydholm A, Selvik G, Onnerfalt R. Low- versus high-viscosity bone cement. Fixation of hip prostheses analyzed by roentgen stereophotogrammetry. *Acta Orthop Scand.* 1987;58(2):106-8.
49. Chotel F, Nguiabanda L, Braillon P, Kohler R, Berard J, Abelin-Genevois K. Induced membrane technique for reconstruction after bone tumor resection in children: a preliminary study. *Orthop Traumatol Surg Res.* 2012;98(3):301-8.
50. Berberich C, Sanz-Ruiz P. Risk assessment of antibiotic resistance development by antibiotic-loaded bone cements: is it a clinical concern? *EFORT Open Rev.* 2019;4(10):576-84.
51. Martínez-Moreno J, Merino V, Nácher A, Rodrigo JL, Climente M, Merino-Sanjuán M. Antibiotic-loaded bone cement as prophylaxis in total joint replacement. *Orthop Surg.* 2017;9(4):331-41.
52. Xu YM, Peng HM, Feng B, Weng XS. Progress of antibiotic-loaded bone cement in joint arthroplasty. *Chin Med J (Engl).* 2020;133(20):2486-94.
53. Wahlig H, Dingeldein E. Antibiotics and bone cements. Experimental and clinical long-term observations. *Acta Orthop Scand.* 1980;51(1):49-56.
54. Salvati EA, Callaghan JJ, Brause BD, Klein RF, Small RD. Reimplantation in infection. Elution of gentamicin from cement and beads. *Clin Orthop Relat Res.* 1986(207):83-93.
55. Cambon-Binder A, Revol M, Hannouche D. Salvage of an osteocutaneous thermonecrosis secondary to tibial reaming by the induced membrane procedure. *Clin Case Rep.* 2017;5(9):1471 – 1476.
56. Schmitt JW, Benden C, Dora C, Werner CML. Is total hip arthroplasty safely performed in lung transplant patients? Current experience from a retrospective study of the Zurich lung transplant cohort. *Patient Saf Surg.* 2016;10:17.
57. Grubhofer F, Imam MA, Wieser K, Achermann Y, Meyer DC, Gerber C. Staged revision with antibiotic spacers for shoulder prosthetic joint infections yields high infection control. *Clin. Orthop. Relat. Res.* 2018;476(1):146 – 152.
58. Wilairatana V, Sinlapavilawan P, Honsawek S, Limpaphayom N. Alteration of inflammatory cytokine production in primary total knee arthroplasty using antibiotic-loaded bone cement. *J Orthop Traumatol.* 2017;18(1):51 – 7.
59. Agni NR, Costa ML, Achten J, Peckham N, Dutton SJ, Png ME, Reed MR, WHITE 8 Investigators. High-dose dual-antibiotic loaded cement for hip hemiarthroplasty in the UK (WHiTE 8): a randomised controlled trial. *Lancet.* 2023;402(10397):196-202.
60. Baryeh K, Wang C, Sochart DH. Periprosthetic femoral fractures around the original cemented polished triple-tapered C-stem femoral implant: a consecutive series of 500 primary total hip arthroplasties with an average follow-up of 15 years. *Arch Orthop Trauma Surg.* 2023;143(7):4511-8.
61. Bourbotte-Salmon F, Ferry T, Cardinale M, Servien E, Rongieras F, Fessy M-H, Bertani A, Laurent F, Buffe-Lidove M, Batailler C, Lustig S, Lyon Bone and Joint Infections Study Group. Rotating hinge knee arthroplasty for revision prosthetic-knee infection: good functional outcomes but a crucial need for superinfection prevention. *Front Surg.* 2021;8:551814.
62. Crego-Vita D, Aedo-Martín D, Garcia-Cañas R, Espigares-Correa A, Sánchez-Pérez C, Berberich CE. Periprosthetic joint infections in femoral neck fracture patients treated with hemiarthroplasty - should we use antibiotic-loaded bone cement? *World J Orthop.* 2022;13(2):150-9.

Titel: SSCP PALACOS R+G

Doc.-No.: 61339

Page 16 of 17

63. Matar HE, Bloch BV, James PJ. High ten-year implant survivorship and low patellofemoral complication rate for S-ROM rotating-hinge implants in revision total knee arthroplasty: a single-centre study. *Bone Jt Open.* 2022;3(3):205-10.
64. Mukka S, Hailer NP, Möller M, Gordon M, Lazarinis S, Rogmark C, Östlund O, Sköldenberg O, Wolf O, DAICY study group. Study protocol: The DAICY trial-dual versus single-antibiotic impregnated cement in primary hemiarthroplasty for femoral neck fracture-a register-based cluster-randomized crossover-controlled trial. *Acta Orthop.* 2022;93:794-800.
65. Ott N, Müller C, Jacobs A, Paul C, Wegmann K, Müller LP, Kabir K. Outcome of geriatric proximal humeral fractures: a comparison between reverse shoulder arthroplasty versus open reduction and internal fixation. *OTA Int.* 2022;5(2 Suppl):188.
66. Reinbacher P, Wittig U, Hauer G, Draschl A, Leithner A, Sadoghi P. Impact of the COVID-19 pandemic on early clinical outcome after total knee arthroplasty: a retrospective comparative analysis. *Arch Orthop Trauma Surg.* 2023;143(6):3319-26.
67. National Joint Registry. Implant Summary Report for Ankles using PALACOS R+G Cement [in primary ankle procedures]. National Joint Registry; 2024.
68. National Joint Registry. Implant Summary Report for Ankles using PALACOS R+G Cement [in revision ankle procedures]. National Joint Registry; 2024.
69. National Joint Registry. Implant Summary Report for Elbows using PALACOS R+G Cement [in primary elbow procedures]. National Joint Registry; 2024.
70. National Joint Registry. Implant Summary Report for Elbows using PALACOS R+G Cement [in revision elbow procedures]. National Joint Registry; 2024.
71. National Joint Registry. Implant Summary Report for Hips using PALACOS R+G Cement [in primary hip procedures]. National Joint Registry; 2024.
72. National Joint Registry. Implant Summary Report for Hips using PALACOS R+G Cement [in revision hip procedures]. National Joint Registry; 2024.
73. National Joint Registry. Implant Summary Report for Knees using PALACOS R+G Cement [in primary knee procedures]. National Joint Registry; 2024.
74. National Joint Registry. Implant Summary Report for Knees using PALACOS R+G Cement [in revision knee procedures]. National Joint Registry; 2024.
75. National Joint Registry. Implant Summary Report for Shoulders using PALACOS R+G Cement [in primary shoulder procedures]. National Joint Registry; 2024.
76. National Joint Registry. Implant Summary Report for Shoulders using PALACOS R+G Cement [in revision shoulder procedures]. National Joint Registry; 2024.
77. LROI. Online LROI annual report 2023. Landelijke Registratie Orthopedische Implantaten; 2023.
78. EPRD. EPRD Data Summary for PALACOS® R+G (B1.1 Bone Cement). EPRD Deutsche Endoprothesenregister gGmbH; 2023.
79. Angelomenos V, Mohaddes M, Kärrholm J, Malchau H, Shareghi B, Itayem R. A prospective randomized study of Refobacin Bone Cement R versus Palacos R + G. *Bone Joint J.* 2024;106-B(5):435-41.
80. Duncan ST, Sabatini F. The use of calcium sulfate/hydroxyapatite bone graft substitute to restore acetabular bone loss in revision total hip arthroplasty. *Arthroplast Today.* 2023;23:101217

Titel: SSCP PALACOS R+G

Doc.-No.: 61339

Page 17 of 17

81. Giovanoulis V, Koutserimpas C, Lepidas N, Vasiliadis AV, Batailler C, Ferry T, Lustig S. Restoring the anatomy of long bones with large septic non-union defects with the masquelet technique. *Maedica (Bucur)*. 2023;18(3):413-9.
82. Holzer LA, Finsterwald MA, Sobhi S, Jones CW, Yates PJ. Application of bone cement directly to the implant in primary total knee arthroplasty. Short-term radiological and clinical follow-up of two different cementing techniques. *Arch Orthop Trauma Surg*. 2024;144(1):333-340.
83. Howie C, McCrosson M, Padgett AM, Sanchez T, McGwin G, Naranje S. The economic and clinical impact of fast- versus slow-setting cement in primary total knee arthroplasty. *Arch Orthop Trauma Surg*. 2024;144(1):15-21.
84. Matar HE, Bloch BV, James PJ. Satisfactory medium- to long-term outcomes of cemented rotating hinge prosthesis in revision total knee arthroplasty. *Bone Jt Open*. 2023;4(10):776-81.
85. O'Donovan P, McAleese T, Harty J. Does lucency equate to revision? A five-year retrospective review of Attune and Triathlon total knee arthroplasty. *Knee Surg Sports Traumatol Arthrosc*. 2023;31(11):4773-4781.
86. Png ME, Costa M, Nickil A, Achten J, Peckham N, Reed MR, WHITE-8 Investigators. Cost-utility analysis of dual-antibiotic cement versus single-antibiotic cement for the treatment of displaced intracapsular hip fractures in older adults. *Bone Joint J*. 2023;105-B(10):1070-1077.
87. Saccomanno MF, Lädemann A, Collin P. Two-stage exchange arthroplasty for periprosthetic reverse shoulder arthroplasty infection provides comparable functional outcomes to primary reverse shoulder arthroplasty. *J Clin Med*. 2024;13(3).
88. Shi J, Yang X, Song M, Zhang X, Xu Y. Clinical effects of early debridement, internal fixation, and Masquelet technique for childhood chronic haematogenous osteomyelitis of long bones. *J Orthop Surg Res*. 2023;18(1):11
89. Swedish Arthroplasty Register. Annual Report 2023. Swedish Arthroplasty Register; 2023.
90. Turnbull G, Blacklock C, Akhtar A, Dunstan E, Ballantyne JA. Experience of an anatomic femoral stem in a UK orthopaedic centre beyond 20 years of follow-up. *Eur J Orthop Surg Traumatol*. 2024;34(4):2155-62.
91. Turnbull GS, Akhtar MA, Dunstan ERR, Ballantyne JA. Experience of an anatomic femoral stem in a United Kingdom center - excellent survivorship and negligible periprosthetic fracture rates at mean 12 years following primary total hip arthroplasty. *J Arthroplasty*. 2024;39(1):187-192.
92. Wagh A, Tandel J, Ballyapally D, Jagtap K, Bharadwaj B. Can "aseptic" looking tibia non-union be result of an unrecognized subclinical infection. *J Orthop Case Rep*. 2023;13(12):75-9.
93. Yang X, Xu X, Li J, Song M, Sun H, Zhang H, Zhang X, Xu Y, Shi J. Management of infected bone defects of the femoral shaft by Masquelet technique: sequential internal fixation and nail with plate augmentation. *BMC Musculoskelet Disord*. 2024;25(1):552.