
Heraeus

PALACOS[®]R+G

High viscosity, radiopaque bone cement
containing Gentamicin

PALACOS® R+G

EN	Instructions for use Working times (Revision status: 2019-07)	5 137
EN (CA)	Note: Applicable in Canada! Instructions for use Working times (Revision status: 2019-07)	9 137
AR	تعليمات الاستخدام أوقات العمل (تاريخ المراجعة: 2019-07)	13 137
BG	Инструкции за употреба Време за обработка при (Дата на актуализиране на текста: 2019-07)	17 137
CS	Návod k použití Doby zpracování (Poslední revise: 2019-07)	21 137
DA	Brugsanvisning Arbejdstider (Aktualiseret: 2019-07)	25 137
DE	Gebrauchsanweisung Verarbeitungszeiten (Revisionsstand: 2019-07)	29 137
EL	Οδηγίες χρήσης Χρόνοι εργασίας (Ημερομηνία αναθεώρησης: 2019-07)	33 137
ES	Instrucciones de uso Tiempos de manipulación (Última revisión: 2019-07)	37 137
ET	Kasutusjuhend Tööjad (Viimati muudetud: 2019-07)	41 137
FI	Käyttöohje Työstöajat (Revisiostatus: 2019-07)	45 137
FR	Notice d'utilisation Temps de travail (Date de révision: 2019-07)	49 137
FR (CA)	Nota: Applicable au Canada! Notice d'utilisation Temps de travail (Date de révision: 2019-07)	53 137

HR	Upute za uporabu Vremena obrade (Stanje revizije: 2019-07)	57 137
HU	Használati utasítás Mégmunkálási (Aktualizálás időpontja: 2019-07)	61 137
IT	Istruzioni per l'uso Tempi di lavorazione (Ultima revisione: 2019-07)	65 137
KO	사용 설명서 작업 시간 (2019년 07월 수정)	69 137
LT	Naudojimo instrukcija Darbo laikas (Tikrinta: 2019-07)	73 137
LV	Lietošanas instrukcija Apstrādes laiks (Pārbaudes statuss: 2019-07)	77 137
NL	Gebruiksaanwijzing Verwerkingstijden (Revisiestatus: 2019-07)	81 137
NO	Brugsanvisning Virketider (Revisjonsstatus: 2019-07)	85 137
PL	Instrukcja użytkowania Czasy obróbki (Ostatnia aktualizacja: 2019-07)	89 137
PT	Instruções de utilização Fases de aplicação (Data da última actualização: 2019-07)	93 137
PT (BR)	Nota: Apenas para fins brasileiros! Instruções de utilização Fases de aplicação (Data da última actualização: 2019-07)	97 137
RO	Instrucțiuni de utilizare Timpi de prelucrare (Stadiu revizie: 2019-07)	101 137
RU	Инструкция по применению Временные фазы для смешивания (Ревизионный статус: 2019-07)	105 137

SK	Návod na použitie	109
	Časy spracovateľnosti (Posledná revízia textu: 2019-07)	137
SL	Navodilo za uporabo	113
	Časovne opredelitve (Datum revizije: 2019-07)	137
SV	Bruksanvisning	117
	Hanteringstider (Revision status: 2019-07)	137
TH	คำแนะนำวิธีใช้	121
	เวลาในการใช้งาน (แก้ไขปรับปรุงครั้งสุดท้ายเมื่อ: 2019-07)	137
TR	Kullanma Talimatı	125
	İşleme süreleri (Bilgilerin son güncellenme tarihi: 2019-07)	137
UK	Інструкція для застосування	129
	Робочий час (Редакція: 2019-07)	137
ZH	使用说明书	第133页
	加工时间 (修订状态: 2019 年 07 月)	第137页

Properties

PALACOS® R+G is a fast-curing, radiopaque, poly(methyl methacrylate)-based bone cement. It contains the aminoglycoside antibiotic gentamicin to protect the cured cement and contiguous tissue against contamination by microbes sensitive to gentamicin.

PALACOS® R+G contains the X-ray contrast medium zirconium dioxide. To improve visibility in the surgical field **PALACOS® R+G** has been coloured with chlorophyll (E141). The bone cement is prepared directly before use by mixing a polymer powder component with a liquid monomer component. A ductile dough forms which cures within a few minutes.

Composition

A pack of **PALACOS® R+G** contains one/two bag/s of gentamicin-containing cement powder (polymer powder) and one/two amber glass ampoule/s (monomer liquid).

Cement powder ingredients:

Gentamicin (in the form of gentamicin sulphate), poly(methylacrylate, methyl methacrylate), zirconium dioxide, benzoyl peroxide and colorant E141.

Ingredients

Pack size	10	20	40	60
Gentamicin base	0.13 g	0.25 g	0.5 g	0.75 g
Cement powder	10.2 g	20.4 g	40.8 g	61.2 g
Monomer liquid	5 ml	10 ml	20 ml	30 ml

The mass ratio is about 32 percent monomer liquid to 68 percent cement powder.

Monomer liquid ingredients:

Methyl methacrylate, N,N-dimethyl-p-toluidine, hydroquinone and colorant E141.

The cement powder is triple packaged. The outer, unsterile protective aluminium bag contains a polyethylene paper bag (peel-off) which is unsterile on the outside and sterile on the inside. Inside this there is another sterile polyethylene paper bag containing the cement powder. The brown glass ampoule with the sterile-filtered monomer liquid is also sterile packed in an ethylene oxide sterilised individual blister pack.

Composition**The cement powder contains:**

poly(methyl acrylate, methyl methacrylate)	82 %
zirconium dioxide	15 %
hydrous benzoyl peroxide	1 %
gentamicin sulphate	2 %

The monomer liquid contains:

methyl methacrylate	98 %
N,N-dimethyl-p-toluidine	2 %

Other constituents:

In the cement powder: chlorophyll VIII (colorant E141)
In the monomer liquid: chlorophyll VIII (colorant E141)
in an oily solution, hydroquinone

Indication

The use of **PALACOS® R+G** with gentamicin is indicated for

- stable anchoring of total or partial joint endoprostheses in bone or filling and stabilising bone defects within the scope of internal fixation treatment or for endoprostheses revision surgery
- primary and secondary plastic reconstruction of bone defects.

Contraindication

PALACOS® R+G with gentamicin must not be used in the following cases

- In case of suspected or proven hypersensitivity to components of the bone cement
- during pregnancy or breast-feeding
- in cases of severe renal insufficiency
- in the presence of an active or incompletely treated infection at the bone site caused by gentamicin non-sensitive strains.

Target population

There is very little evidence with children and adolescents, therefore it is not recommended to use **PALACOS® R+G**. If no other option is available, for example with present surgical trauma, the decision whether to use **PALACOS® R+G** lies with the attending surgeon.

Target user group

Healthcare professionals in a clinical context.

Side effects

General

Rarely there is a temporary drop in blood pressure after preparation of the prosthesis bed or directly after the implantation of bone cements and the endoprosthesis.

In individual cases severe complications can occur such as severe allergic reactions which are associated with cardiac arrest, anaphylactic shock or even sudden death.

In order to avoid pulmocardiovascular complications such as pulmonary embolism and cardiac arrest it is recommended that the implantation site be irrigated thoroughly with an isotonic solution (pulse lavage application) before the bone cement is introduced. In the event of pulmonary or cardiovascular events it is necessary to monitor blood volume and possibly increase it.

In the case of acute respiratory insufficiency anaesthesiological measures should be taken.

The following adverse effects have been observed when using poly(methyl methacrylate) bone cements: thrombophlebitis, haemorrhage, trochanteric bursitis.

Other side effects that were observed are heterotopic bone formation, osteolysis due to cement fragments, myocardial infarction, brief cardiac arrhythmia, cerebrovascular incident.

Gentamicin

Due to the gentamicin contained in **PALACOS® R+G**, side effects that are typical for gentamicin can occur:

- impairment to auditory and vestibular nerves
- renal toxicity: Adverse renal effects have been reported, as demonstrated by the presence of cells or protein in the urine or by rising of serum creatinine or oliguria. They occur more frequently in patients with a history of renal impairment. Application of clinical doses of gentamicin has led to the occasional appearance of Fanconi's syndrome or a Bartter's-like syndrome
- neuromuscular blockade/ neurotoxicity: Serious adverse effects on both vestibular and auditory branches of the 8th cranial nerves have been reported, primarily in patients with renal impairment. Symptoms include dizziness, vertigo, tinnitus, roaring in the ears and hearing loss, which, as with the other aminoglycosides, may be irreversible. Other factors which may increase the risk of toxicity include dehydration and previous exposure to other ototoxic drugs
- rare cases of paresthesia, tetany and amyasthenia
- rare cases of allergic reactions (exanthema, urticaria, anaphylactic reactions).

Like all aminoglycosides gentamicin is also potentially nephrotoxic and/or ototoxic. In most cases, damage to hearing caused by gentamicin is irreversible, whereas renal damage is reversible.

Other reported adverse reactions possibly related to gentamicin include: respiratory depression, lethargy, confusion, depression, visual disturbances, decreased appetite, weight loss, and hypotension or hypertension; rash, itching, urticaria, generalised burning, laryngeal oedema, anaphylactoid reactions, fever, and headache, nausea, vomiting, increased salivation, and stomatitis; purpura, pseudotumor cerebri, acute organic brain syndrome, pulmonary fibrosis, alopecia, joint pain, transient hepatomegaly, and splenomegaly.

The use of **PALACOS® R+G**, results in a desired high bactericidal but locally limited concentration of gentamicin.

Therefore, the occurrence of the afore-mentioned side effects is highly unlikely. Care should be taken and the level of gentamicin should be monitored in patients with impaired renal function, especially if treated with parenteral aminoglycosides or muscle relaxants simultaneously. The same applies to patients with previous neuromuscular disorders (e.g. myasthenia gravis, Parkinson's disease). Allergic reactions can occur dosage-independent.

Interactions

Gentamicin

Owing to the administration of muscle relaxants and ether, the neuromuscular blocking properties of gentamicin may be intensified. However, in view of the very low serum levels this is relatively unlikely in patients with healthy kidneys. The probability of interactions occurring increases in proportion to the serum levels of gentamicin, e.g. in patients with impaired renal function.

Interactions of gentamicin have been reported to occur with the following groups of drugs:

Concurrent or sequential use of neurotoxic and / or nephrotoxic pharmaceuticals, including other aminoglycosides, polymyxin B, colistin, cisplatin, ciclosporines, foscarnet, vancomycin, amphotericin B, clindamycin and cephalosporins.

Potent diuretics: Ethacrynic acid, furosemide, or other potent diuretics which may themselves cause ototoxicity or enhance aminoglycoside toxicity by altering antibiotic concentrations in serum and tissue.

Muscle relaxants: in particular succinyl choline and tubocurarine, decamethonium, halogenated hydro-carbon inhalation anaesthetics or opioid analgesics. Gentamicin may interact with these medications and result in skeletal muscle weakness and respiratory depression (apnea).

Concurrent application of these medications and gentamicin during surgery or in the post-operative period should be monitored with caution, especially if there is a possibility of incomplete reversal of neuromuscular blockade post-operatively.

Moreover, complications are possible that can occur during any surgical procedure.

Precautions

Use by operating staff

Before using **PALACOS® R+G** the user should be well-acquainted with its properties, handling and application. The user is advised to practise the entire procedure of mixing, handling and introducing **PALACOS® R+G** before using it for the first time. Detailed knowledge is necessary even if mixing systems and syringes are being used for application of the cement.

The liquid monomer is highly volatile and flammable. Ignition of monomer fumes caused by use of electrocautery devices in surgical sites near freshly implanted bone cements has been reported. The monomer is also a strong lipid solvent and should not come into direct contact with the body.

When handling the monomer or the prepared cement **PALACOS® R+G** gloves must be worn which provide the necessary protection against penetration of the skin by the monomer methyl methacrylate.

Gloves made of PVP (three-ply polyethylene, ethylene vinyl alcohol copolymer, polyethylene) and Viton®/butyl gloves have proved to provide good protection over a lengthy period. For safety's sake it is recommended that two pairs of gloves be worn over one another, e.g. one polyethylene surgical glove over an inner pair of standard latex surgical gloves.

The use of latex or polystyrene-butadiene gloves on their own is inadequate. Please make enquiries with your supplier to establish which gloves are suitable for such an application. The monomer vapours can irritate respiratory passages and eyes and possibly cause damage to the liver. Skin irritations have been described which are due to contact with the monomer.

Manufacturers of soft contact lenses recommend removing the lenses in the presence of damaging or irritant vapours. Since soft contact lenses are permeable to liquids and gases, they should not be worn in the operating theatre if methyl methacrylate is being used.

Use on the patient

Blood pressure, pulse and breathing must be monitored carefully during and directly after introduction of the bone cement. Any significant change in these vital signs must be eliminated without delay by taking appropriate action.

When using **PALACOS® R+G** the prepared bone should be carefully cleaned, aspirated and dried just before the bone cement is placed.

Incompatibilities

Aqueous solutions (e.g. ones containing antibiotics) must not be added to the bone cement because they have a considerable detrimental effect on the physical and mechanical properties of the cement.

Amount required

After mixing the cement powder with the monomer liquid, a fast-curing ductile dough develops which is introduced to bone cavities for anchorage and/or filling purposes.

PALACOS® R+G is coloured green in order to make the cement clearly visible in the surgical field.

A dose is prepared by mixing the entire contents of the bag of cement powder with all the monomer liquid of an ampoule. The quantity of cement dough required depends on the specific surgical intervention and on the technique being used.

At least one additional dose of **PALACOS® R+G** should be available before commencing the operation.


Preparation

Preparation:

Before opening the non-sterile aluminium protective pouch (see above) move the contents down (1 polyethylene-paper bag) by shaking or tapping in order to ensure that when the bag is cut open at the top the content is not damaged.

The polyethylene-paper bag and the ampoule may only be opened under sterile conditions. For this purpose the sterile components (inner PE-paper bag and glass ampoule) are sterile when they are delivered.

Opening under sterile conditions:

The opening flaps from the bag top assist to detach the foil from the paper.

In order to grip as much of the opening flaps as possible, the side of the paper/foil should be kept between thumb, index and middle finger.

Please use the whole thumb surface to grip foil and paper side and take off each side evenly.

Open the outer PE-paper bag at the special point under sterile conditions so that the inner PE-paper bag remains sterile when it is removed. Also open the blister pack at the special point under sterile conditions so that the glass ampoule remains sterile when it is removed.

Before opening the inner PE-paper bag, move the contents down by shaking or tapping in order to ensure that when the bag is cut open at the top no powder is lost. To make it easier to open the glass ampoule the latter is provided with a predetermined breaking point at the transition to the head of the ampoule.

Do not open the ampoule over the mixing device to prevent contamination of the cement with glass fragments.

The 10ml- and 20ml- ampoule is provided with a snapping-off device (tube) in order to facilitate the opening procedure. If this is the case, take hold of the fitted snapping-off device instead of the ampoule head and break off the ampoule head over it. When the ampoule head has snapped off, it remains inside the tube.

Mixing the components:

It is advisable to first measure out the liquid and then add the powder. If this order is reversed, powder nests are more likely to form as a result of polymerisation commencing immediately at the surface. Both components, i.e. the relative proportions of powder and monomer, are precisely matched. The pouch and ampoule must therefore be emptied completely if an optimal mix is to be achieved. The components can be mixed with a vacuum mixing system or by hand.

The mixing, waiting, working and curing times of **PALACOS® R+G** are shown in the diagrams at the end of the instructions for use. Please note that these are stated for guidance only because working time and curing time depend on temperature, mixing and humidity, whereby the direct ambient temperatures are important, e.g. of cement powder, mixing system, bench and hands. A high temperature accelerates the waiting, working and curing times.

Preparation with a vacuum mixing system

In order to obtain a bone cement with a reduced number of air inclusions, the liquid and the powder are mixed in a vacuum. For this purpose an airtight mixing system has to be used, ensuring that a sufficient vacuum builds up in the mixing vessel quickly (pressure approx. 200 mbar absolute). It is recommended that the cement components be precooled at 4–7°C for at least 24 hours. The cement components should only be removed from the cooling vessel and filled into the mixing vessel just before mixing.

Perform filling and mixing under sterile conditions. The mixing time is 30 seconds unless otherwise recommended. As a result of prior cooling the working time and curing time are longer. Initial viscosity is slightly reduced for **PALACOS® R+G**, compared to not pre-chilled cement. For details of mixing technique please refer to the instructions for the mixing system being used. Always mix the entire contents of a bag with the entire contents of an ampoule of monomer liquid.

Preparation mixing by hand

The cement components should only be filled into the mixing vessel just before mixing. Perform filling and mixing under sterile conditions – mixing time is 30 seconds. During this time the two components are mixed with one another while stirring evenly. Always mix the entire contents of a bag with the entire contents of an ampoule of monomer liquid.

Use of the bone cement

The bone cement can be applied as soon as the doughy bone cement no longer adheres to the gloves. The application period depends on the temperature of the material and the room temperature. To ensure adequate fixation the prosthesis must be introduced and held within the time window allowed for working, until the bone cement has set completely. Remove any surplus cement as long as it is still soft.

Cranial vault reconstruction

In the management of major skull defects, after careful preparation of the bone gap the dura mater is first of all covered with moist cotton wool or cellulose. A thin plastic or aluminium foil is laid over it for further protection. The mixed pasty dough is introduced to the prepared bone gap and moulded at the bone margins to the required thickness

of 4–5 mm. While it is setting, irrigate with normal saline solution in order to dissipate the heat caused by polymerisation. When the reconstruction has almost hardened it is then removed, corrected at the margins and provided with perforations, through which epidural fluid can drain away and connective tissue can grow. After removal of the cotton wool/cellulose and plastic foil the prosthesis is fixed in place at three or four points with non-absorbable suture.

Storage

Do not store above 25°C (77°F).

Shelf life/sterility

The shelf life is printed on the folding box, the protective aluminium bag and the inner bag. Do not use **PALACOS® R+G** if this date indicated has expired. The information on the blister pack may be different for production reasons.

The contents of opened or damaged aluminium sachets or ampoule blisters must not be re-sterilized and thus must be discarded. If the cement powder has turned yellow, do not use **PALACOS® R+G**. The cement powder and packagings are sterilised by gassing with ethylene oxide. The monomer liquid has been sterilised by filtration.

Disposal

Single components of the bone cement, cured solid material as well as (uncleaned) packaging material must be disposed by following the regulations of the local authorities. Dispose of the polymer component in an authorized waste facility. The liquid component should be evaporated under a well-ventilated hood or absorbed by an inert material and transferred in a suitable container for disposal.

EN (CA)

PALACOS® R+G

Note: Applicable in Canada!

Properties

PALACOS® R+G is a fast-curing, radiopaque, poly(methyl methacrylate)-based bone cement. It contains the aminoglycoside antibiotic gentamicin to protect the cured cement and contiguous tissue against contamination by microbes sensitive to gentamicin.

PALACOS® R+G contains the X-ray contrast medium zirconium dioxide. To improve visibility in the surgical field **PALACOS® R+G** has been coloured with chlorophyll (E141). The bone cement is prepared directly before use by mixing a polymer powder component with a liquid monomer component. A ductile dough forms which cures within a few minutes.

Composition

A pack of **PALACOS® R+G** contains one/two bag/s of gentamicin-containing cement powder (polymer powder) and one/two amber glass ampoule/s (monomer liquid).

Cement powder ingredients:

Gentamicin (in the form of gentamicin sulphate), poly(methylacrylate, methyl methacrylate), zirconium dioxide, benzoyl peroxide and colorant E141.

Ingredients

Pack size	10	20	40	60
Gentamicin base	0.13g	0.25g	0.5g	0.75g
Cement powder	10.2g	20.4g	40.8g	61.2g
Monomer liquid	5ml	10ml	20ml	30ml

The mass ratio is about 32 percent monomer liquid to 68 percent cement powder.

Monomer liquid ingredients:

Methyl methacrylate, N,N-dimethyl-p-toluidine, hydroquinone and colorant E141.

The cement powder is triple packaged. The outer, unsterile protective aluminium bag contains a polyethylene paper bag (peel-off) which is unsterile on the outside and sterile on the inside. Inside this there is another sterile polyethylene paper bag containing the cement powder. The brown glass ampoule with the sterile-filtered monomer liquid is also sterile packed in an ethylene oxide sterilised individual blister pack.

Composition

The cement powder contains:

poly(methyl acrylate, methyl methacrylate)	82%
zirconium dioxide	15%
hydrous benzoyl peroxide	1%
gentamicin sulphate	2%

The monomer liquid contains:

methyl methacrylate	98%
N,N-dimethyl-p-toluidine	2%

Other constituents:

In the cement powder: chlorophyll VIII (colorant E141)
In the monomer liquid: chlorophyll VIII (colorant E141)
in an oily solution, hydroquinone

Indication

- The use of **PALACOS® R+G** with gentamicin is indicated for
- stable anchoring of total or partial joint endoprostheses in bone or filling and stabilising bone defects within the scope of internal fixation treatment or for endoprostheses revision surgery
 - primary and secondary coverage of skull bone defects.

Contraindication

PALACOS® R+G with gentamicin must not be used in the following cases

- In case of suspected or proven hypersensitivity to components of the bone cement
- during pregnancy or breast-feeding
- in cases of severe renal insufficiency
- in the presence of an active or incompletely treated infection at the bone site caused by gentamicin non-sensitive strains.

Target population

There is very little evidence with children and adolescents, therefore it is not recommended to use **PALACOS® R+G**. If no other option is available, for example with present surgical trauma, the decision whether to use **PALACOS® R+G** lies with the attending surgeon.

Target user group

Healthcare professionals in a clinical context.

Side effects

General

Rarely there is a temporary drop in blood pressure after preparation of the prosthesis bed or directly after the implantation of bone cements and the endoprosthesis.

In individual cases severe complications can occur such as severe allergic reactions which are associated with cardiac arrest, anaphylactic shock or even sudden death.

In order to avoid pulmocardiovascular complications such as pulmonary embolism and cardiac arrest it is recommended that the implantation site be irrigated thoroughly with an isotonic solution (pulse lavage application) before the bone cement is introduced. In the event of pulmonary or cardiovascular events it is necessary to monitor blood volume and possibly increase it.

In the case of acute respiratory insufficiency anaesthesiological measures should be taken.

The following adverse effects have been observed when using poly(methyl methacrylate) bone cements: thrombophlebitis, haemorrhage, trochanteric bursitis.

Other side effects that were observed are heterotopic bone formation, osteolysis due to cement fragments, myocardial infarction, brief cardiac arrhythmia, cerebrovascular incident.

Gentamicin

Due to the gentamicin contained in **PALACOS® R+G**, side effects that are typical for gentamicin can occur:

- impairment to auditory and vestibular nerves
- renal toxicity: Adverse renal effects have been reported, as demonstrated by the presence of cells or protein in the urine or by rising of serum creatinine or oliguria. They occur more frequently in patients with a history of renal impairment. Application of clinical doses of gentamicin has led to the occasional appearance of Fanconi's syndrome or a Bartter's-like syndrome
- neuromuscular blockade/ neurotoxicity: Serious adverse effects on both vestibular and auditory branches of the 8th cranial nerves have been reported, primarily in patients with renal impairment. Symptoms include dizziness, vertigo, tinnitus, roaring in the ears and hearing loss, which, as with the other aminoglycosides, may be irreversible. Other factors which may increase the risk of toxicity include dehydration and previous exposure to other ototoxic drugs
- rare cases of paresthesia, tetany and amyosthenia
- rare cases of allergic reactions (exanthema, urticaria, anaphylactic reactions).

Like all aminoglycosides gentamicin is also potentially nephrotoxic and/or ototoxic. In most cases, damage to hearing caused by gentamicin is irreversible, whereas renal damage is reversible.

Other reported adverse reactions possibly related to gentamicin include: respiratory depression, lethargy, confusion, depression, visual disturbances, decreased appetite, weight loss, and hypotension or hypertension; rash, itching, urticaria, generalised burning, laryngeal oedema, anaphylactoid reactions, fever, and headache, nausea, vomiting, increased salivation, and stomatitis; purpura, pseudotumor cerebri, acute organic brain syndrome, pulmonary fibrosis, alopecia, joint pain, transient hepatomegaly, and splenomegaly.

The use of **PALACOS® R+G**, results in a desired high bactericidal but locally limited concentration of gentamicin.

Therefore, the occurrence of the afore-mentioned side effects is highly unlikely. Care should be taken and the level of gentamicin should be monitored in patients with impaired renal function, especially if treated with parenteral aminoglycosides or muscle relaxants simultaneously. The same applies to patients with previous neuromuscular disorders (e.g. myasthenia gravis, Parkinson's disease). Allergic reactions can occur dosage-independent.

Interactions

Gentamicin

Owing to the administration of muscle relaxants and ether, the neuromuscular blocking properties of gentamicin may be intensified. However, in view of the very low serum levels this is relatively unlikely in patients with healthy kidneys. The probability of interactions occurring increases in proportion to the serum levels of gentamicin, e.g. in patients with impaired renal function.

Interactions of gentamicin have been reported to occur with the following groups of drugs:

Concurrent or sequential use of neurotoxic and / or nephrotoxic pharmaceuticals, including other aminoglycosides, polymyxin B, colistin, cisplatin, ciclosporines, foscarnet, vancomycin, amphotericin B, clindamycin and cephalosporins.

Potent diuretics: Ethacrynic acid, furosemide, or other potent diuretics which may themselves cause ototoxicity or enhance aminoglycoside toxicity by altering antibiotic concentrations in serum and tissue.

Muscle relaxants: in particular succinyl choline and tubocurarine, decamethonium, halogenated hydro-carbon inhalation anaesthetics or opioid analgesics. Gentamicin may interact with these medications and result in skeletal muscle weakness and respiratory depression (apnea).

Concurrent application of these medications and gentamicin during surgery or in the post-operative period should be monitored with caution, especially if there is a possibility of incomplete reversal of neuromuscular blockade post-operatively.

Moreover, complications are possible that can occur during any surgical procedure.

Precautions

Use by operating staff

Before using **PALACOS® R+G** the user should be well-acquainted with its properties, handling and application. The user is advised to practise the entire procedure of mixing, handling and introducing **PALACOS® R+G** before using it for the first time. Detailed knowledge is necessary even if mixing systems and syringes are being used for application of the cement.

The liquid monomer is highly volatile and flammable. Ignition of monomer fumes caused by use of electrocautery devices in surgical sites near freshly implanted bone cements has been reported. The monomer is also a strong lipid solvent and should not come into direct contact with the body.

When handling the monomer or the prepared cement **PALACOS® R+G** gloves must be worn which provide the necessary protection against penetration of the skin by the monomer methyl methacrylate.

Gloves made of PVP (three-ply polyethylene, ethylene vinyl alcohol copolymer, polyethylene) and Viton®/butyl gloves have proved to provide good protection over a lengthy period. For safety's sake it is recommended that two pairs of gloves be worn over one another, e.g. one polyethylene surgical glove over an inner pair of standard latex surgical gloves.

The use of latex or polystyrene-butadiene gloves on their own is inadequate. Please make enquiries with your supplier to establish which gloves are suitable for such an application. The monomer vapours can irritate respiratory passages and eyes and possibly cause damage to the liver. Skin irritations have been described which are due to contact with the monomer.

Manufacturers of soft contact lenses recommend removing the lenses in the presence of damaging or irritant vapours. Since soft contact lenses are permeable to liquids and gases, they should not be worn in the operating theatre if methyl methacrylate is being used.

Use on the patient

Blood pressure, pulse and breathing must be monitored carefully during and directly after introduction of the bone cement. Any significant change in these vital signs must be eliminated without delay by taking appropriate action.

When using **PALACOS® R+G** the prepared bone should be carefully cleaned, aspirated and dried just before the bone cement is placed.

Incompatibilities

Aqueous solutions (e.g. ones containing antibiotics) must not be added to the bone cement because they have a considerable detrimental effect on the physical and mechanical properties of the cement.

Amount required

After mixing the cement powder with the monomer liquid, a fast-curing ductile dough develops which is introduced to bone cavities for anchorage and/or filling purposes.

PALACOS® R+G is coloured green in order to make the cement clearly visible in the surgical field.

A dose is prepared by mixing the entire contents of the bag of cement powder with all the monomer liquid of an ampoule. The quantity of cement dough required depends on the specific surgical intervention and on the technique being used.

At least one additional dose of **PALACOS® R+G** should be available before commencing the operation.


Preparation

Preparation:

Before opening the non-sterile aluminium protective pouch (see above) move the contents down (1 polyethylene-paper bag) by shaking or tapping in order to ensure that when the bag is cut open at the top the content is not damaged.

The polyethylene-paper bag and the ampoule may only be opened under sterile conditions. For this purpose the sterile components (inner PE-paper bag and glass ampoule) are sterile when they are delivered.

Opening under sterile conditions:

The opening flaps from the bag top assist to detach the foil from the paper.

In order to grip as much of the opening flaps as possible, the side of the paper/foil should be kept between thumb, index and middle finger.

Please use the whole thumb surface to grip foil and paper side and take off each side evenly.

Open the outer PE-paper bag at the special point under sterile conditions so that the inner PE-paper bag remains sterile when it is removed. Also open the blister pack at the special point under sterile conditions so that the glass ampoule remains sterile when it is removed.

Before opening the inner PE-paper bag, move the contents down by shaking or tapping in order to ensure that when the bag is cut open at the top no powder is lost. To make it easier to open the glass ampoule the latter is provided with a predetermined breaking point at the transition to the head of the ampoule.

Do not open the ampoule over the mixing device to prevent contamination of the cement with glass fragments.

The 10ml- and 20ml- ampoule is provided with a snapping-off device (tube) in order to facilitate the opening procedure. If this is the case, take hold of the fitted snapping-off device instead of the ampoule head and break off the ampoule head over it. When the ampoule head has snapped off, it remains inside the tube.

Mixing the components:

It is advisable to first measure out the liquid and then add the powder. If this order is reversed, powder nests are more likely to form as a result of polymerisation commencing immediately at the surface. Both components, i.e. the relative proportions of powder and monomer, are precisely matched. The pouch and ampoule must therefore be emptied completely if an optimal mix is to be achieved. The components can be mixed with a vacuum mixing system or by hand.

The mixing, waiting, working and curing times of **PALACOS® R+G** are shown in the diagrams at the end of the instructions for use. Please note that these are stated for guidance only because working time and curing time depend on temperature, mixing and humidity, whereby the direct ambient temperatures are important, e.g. of cement powder, mixing system, bench and hands. A high temperature accelerates the waiting, working and curing times.

Preparation with a vacuum mixing system

In order to obtain a bone cement with a reduced number of air inclusions, the liquid and the powder are mixed in a vacuum. For this purpose an airtight mixing system has to be used, ensuring that a sufficient vacuum builds up in the mixing vessel quickly (pressure approx. 200 mbar absolute). It is recommended that the cement components be precooled at 4–7°C for at least 24 hours. The cement components should only be removed from the cooling vessel and filled into the mixing vessel just before mixing.

Perform filling and mixing under sterile conditions. The mixing time is 30 seconds unless otherwise recommended. As a result of prior cooling the working time and curing time are longer. Initial viscosity is slightly reduced for **PALACOS® R+G**, compared to not pre-chilled cement. For details of mixing technique please refer to the instructions for the mixing system being used. Always mix the entire contents of a bag with the entire contents of an ampoule of monomer liquid.

Preparation mixing by hand

The cement components should only be filled into the mixing vessel just before mixing. Perform filling and mixing under sterile conditions – mixing time is 30 seconds. During this time the two components are mixed with one another while stirring evenly. Always mix the entire contents of a bag with the entire contents of an ampoule of monomer liquid.

Use of the bone cement

The bone cement can be applied as soon as the doughy bone cement no longer adheres to the gloves. The application period depends on the temperature of the material and the room temperature. To ensure adequate fixation the prosthesis must be introduced and held within the time window allowed for working, until the bone cement has set completely. Remove any surplus cement as long as it is still soft.

Cranial vault reconstruction

In the management of major skull defects, after careful preparation of the bone gap the dura mater is first of all covered with moist cotton wool or cellulose. A thin plastic or aluminium foil is laid over it for further protection. The mixed pasty dough is introduced to the prepared bone gap and moulded at the bone margins to the required thickness

of 4–5 mm. While it is setting, irrigate with normal saline solution in order to dissipate the heat caused by polymerisation. When the reconstruction has almost hardened it is then removed, corrected at the margins and provided with perforations, through which epidural fluid can drain away and connective tissue can grow. After removal of the cotton wool/cellulose and plastic foil the prosthesis is fixed in place at three or four points with non-absorbable suture.

Storage

Do not store above 25°C (77°F).

Shelf life/sterility

The shelf life is printed on the folding box, the protective aluminium bag and the inner bag. Do not use **PALACOS® R+G** if this date indicated has expired. The information on the blister pack may be different for production reasons.

The contents of opened or damaged aluminium sachets or ampoule blisters must not be re-sterilized and thus must be discarded. If the cement powder has turned yellow, do not use **PALACOS® R+G**. The cement powder and packagings are sterilised by gassing with ethylene oxide. The monomer liquid has been sterilised by filtration.

Disposal

Single components of the bone cement, cured solid material as well as (uncleaned) packaging material must be disposed by following the regulations of the local authorities. Dispose of the polymer component in an authorized waste facility. The liquid component should be evaporated under a well-ventilated hood or absorbed by an inert material and transferred in a suitable container for disposal.

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


Working times for manual mixing (not pre-chilled bone cement) [A]

Working times for vacuum mixing (pre-chilled bone cement) [B]

Working times for vacuum mixing (not pre-chilled bone cement) [C]

■ Mixing ① ■ Waiting ② ■ Application ③ ■ Setting ④

Test conditions: Not pre-chilled vacuum mixing system PALAMIX®, 55% humidity. [D]

	①	②	③	④	⑤	⑥
EN	Mixing	Waiting	Application	Setting	Time (min)	Room temperature (°C)
AR	أخلط	انتظر	الاستخدام	تصلب	(الوقت (min)	درجة الحرارة (°C)
BG	Смесване	Изчакване	Приложение	Втвърдяване	Време (мину)	Температура на помещението (°C)
CS	Míchání	Čekání	Použití	Tuhnutí	Čas (min)	Pokojoiná teplota (°C)
DA	Blanding	Ventetid	Anvendelse	ærdning	Tid (min)	Rumtemperatur (°C)
DE	Mischen	Warten	Applikation	Aushärtung	Zeit (min)	Raumtemperatur (°C)
EL	Ανάμειξη	Αναμονή	Εφαρμογή	Πήξη	Χρόνος (λεπτά)	Θερμοκρασία δωματίου (°C)
ES	Mezcla	Espere	Aplicación	Endurecimiento	Tiempo (min.)	Temperatura ambiente (°C)
ET	Segamine	Ootamine	Kasutamine	Kõvastumine	Aeg (minutid)	Operatsioonisaali temperatuur (°C)
FI	Sekoittaminen	Odotusaika	Applikointi	Kovettuminen	Aika (min)	Huoneenlämpötila (°C)
FR	Mélange	Attente	Application	Prise	Temps (min)	Température ambiante (°C)
HR	Miješanje	Čekanje	Nanošenje	Stvrđnjavanje	Vrijeme (min)	Temperatura prostorije (°C)
HU	Keverés	Várakozás	Alkalmazás	Kötés	Idő (min)	Szobahőmérséklet (°C)
IT	Miscelazione	Attesa	Applicazione	Indurimento	Tempo (min)	Temperatura ambiente (°C)
KO	혼합	대기	적용	세팅	시간 (min)	온도 (°C)
LT	Maišymas	Laukimas	Taikymas	Kietėjimas	Laikas (min.)	Patalpos temperatūra (°C)
LV	Sajaukšana	Nogaidīšana	Lietošana	Fiksēšanās	Laiks (min)	Istabas temperatūra (°C)
NL	Mengfase	Wachtfase	Applicatiefase	Uithardingsfase	Tijd (min)	Kamertemperatuur (°C)
NO	Blanding	Ventetid	Applisering	Herding	Tid (min)	Romtemperatur (°C)
PL	Mieszanie	Czekanie	Aplikacja	Wiązanie	Czas (min)	Temperatura pokojowa (°C)
PT	Mistura	Espere	Aplicação	Polimerização	Tempo (min)	Temperatura ambiente (°C)
RO	Amestecare	Așteptare	Aplicare	Întărire	Timp (min)	Temperatura încăperii (°C)
RU	Смешивание	Ожидание	Аппликация	Отверждение	Время (мин.)	Комнатная температура (°C)
SK	Miešanie	Čakanie	Aplikácia	Tvrdenutie	Čas (min)	Izbová teplota (°C)
SL	Mešanje	Čakanje	Vnos	Strjevanje	Čas (min)	Sobna temperatura (°C)
SV	Blandning	Väntetid	Applisering	Hårdning	Tid (min)	Rumtemperatur (°C)
TH	การผสม	รอ	การใช้	การแข็งตัว	เวลา (min)	อุณหภูมิ (°C)
TR	Karıştırma	Bekleme	Uygulama	Sertleşme	Süre (dk)	Oda sıcaklığı (°C)
UK	Смешивание	Ожидание	Аппликация	Отверждение	Время (мин.)	Комнатная температура (°C)
ZH	混合	等待	应用	固化	时间 (分钟)	室温 (°C)

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Do not re-use

Do not re-sterilize

Do not use if the product sterile barrier system or its packaging is compromised

Catalogue number

Use by date

Batch code

Flammable liquid – Flashpoint 10 °C

Causes skin irritation

- EN Do not use if the product sterile barrier system or its packaging is compromised.
- AR لا تستخدم المنتج إذا كان هناك ضرراً في نظام حاجز التعقيم أو التغليف الخاص بالمنتج به تالفاً.
- BG Не използвайте при нарушена стерилна преградна система на продукта или неговата опаковка.
- CS Nepoužívejte, pokud je narušena sterilní bariéra systému nebo obalu.
- DA Må ikke anvendes, hvis produktets sterile barriersystem eller dets emballage er beskadiget.
- DE Bei beschädigtem Sterilbarriersystem des Produktes oder seiner Verpackung nicht verwenden.
- EL Μην το χρησιμοποιείτε εάν η αποστειρωμένη συσκευασία έχει ανοιχτεί ή έχει υποστεί ζημιά.
- ES No utilizar el producto en caso de daños en el sistema de barrera estéril o en el envase.
- ET Mitte kasutada, kui toote steriiliden suojausjärjestelmä või pakend on kahjustunud.
- FI Älä käytä tuotetta, jos sen steriilijärjestelmä tai pakkaus on vahingoittunut.
- FR Ne pas utiliser en cas de détérioration du système de barrière stérile du produit ou de son emballage.
- HR Ne koristiti ako su sustav sterilnih barijera proizvoda ili njegova ambalaža oštećeni.
- HU Ne használja, ha a termék sterilgát-rendszer vagy csomagolóasa sérült.
- IT Non utilizzare in caso di danneggiamento alla barriera sterile del prodotto o alla confezione.
- KO 제품 멸균 장벽 또는 그 포장에 손상된 경우 사용하지 마십시오.
- LT Nenaudokite, jei produkto sterilumo barjero sistema ar pakuotė yra pažeista.
- LV Neizmantot, ja produkta sterilā barjera vai iepakojums ir bojāts.
- NL Niet meer gebruiken als het steriele barriersysteem van het product of de verpakking ervan beschadigd is.
- NO Må ikke brukes om produktets sterile barriere er brutt, eller pakningen er ødelagt.
- PL Nie stosować w przypadku uszkodzenia systemu bariery sterylnej produktu lub jego opakowania.
- PT Não utilizar em caso de danos no sistema de barreira estéril do produto ou da respetiva embalagem.
- RO A nu se utilizeza dacă sistemul de barieră sterilă sau ambalajul produsului este deteriorat.
- RU Не использовать при повреждении стерильной барьерной системы изделия или его упаковки.
- SK Nepoužívať, ak došlo k strate sterility produktu alebo k poškodeniu jeho obalu.
- SL Ne uporabljajte, če je ovojnina poškodovana oziroma če izdelek ni sterilen.
- SV Får ej användas om produktens sterilbarriärsystem eller produktförpackningen är skadad.
- TH ห้ามใช้ หากระบบป้องกันการปลอดเชื้อของผลิตภัณฑ์หรือบรรจุภัณฑ์ไม่อยู่ในสภาพที่สมบูรณ์.
- TR Ürün steril bariyer sisteminin veya ambalajının hasarlı olması durumunda kullanmayın.
- UK Не использовать при повреждении стерильной барьерной системы изделия или его упаковки.
- ZH 如果产品无菌屏障系统或其包装已受损，请不要使用

Heraeus Medical GmbH
 Philipp-Reis-Straße 8/13
 61273 Wehrheim, Germany
 Phone: +49 (0)6181 353399
www.heraeus-medical.com