
Heraeus

# OSTEOPAL<sup>®</sup> V

Radiopaque bone cement  
for filling and stabilising vertebral bodies

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### Properties

**OSTEOPAL®V** is a radiopaque, low-viscosity bone cement, based on polymethylmethacrylate, that is used to fill and stabilise vertebral bodies.

**OSTEOPAL®V** contains zirconium dioxide as an X-ray contrast medium. The colourant chlorophyll (E 141) is mixed with **OSTEOPAL®V** to improve visibility in the surgical field. The bone cement is prepared immediately before using by mixing the polymer powder component and the liquid monomer component. This produces a low-viscosity paste that is inserted into the vertebral body using an application system and allowed to set.

### Composition

One package of **OSTEOPAL®V** contains one pouch with bone cement powder (polymer powder) and one glass ampoule (monomer liquid).

#### Cement powder ingredients:

Poly(methyl acrylate, methyl methacrylate), zirconium dioxide, benzoyl peroxide, colourant E 141.

#### Monomer liquid ingredients:

Methyl methacrylate, N,N-dimethyl-p-toluidine, hydroquinone, colourant E 141.

The cement powder is triple packaged. The outer sterile protective aluminium pouch contains a polyethylene paper bag (peelable) that is sterile on the outside and sterile on the inside. Inside this polyethylene paper bag is another sterile polyethylene paper bag that contains the cement powder.

The glass ampoule containing the sterile-filtered monomer liquid is also sterile packed in an ethylene oxide sterilised ampoule blister pack.

### Indications

**OSTEOPAL®V** is indicated for augmentation and stabilisation of vertebral bodies in percutaneous vertebroplasty and balloon kyphoplasty

- for compression fractures of the vertebral body
- for vertebral body tumours (metastases and myelomas)
- for symptomatic vertebral hemangioma.

The vertebroplasty or kyphoplasty procedure must only be considered a palliative treatment to stabilise the vertebral body. The procedure does not treat the underlying disease (osteoporosis, cancer).

### Contraindications

Haemorrhagic diathesis and infections are absolute contraindications. Lesions of the vertebral body with epidural extension are relative contraindications due to the risk of spinal cord compression.

**OSTEOPAL®V** must not be used

- in patients with suspected or known hypersensitivity to constituents of the bone cement
- during pregnancy or breast-feeding.

### Target patient population

Adult population.

There is very little evidence with children and adolescents. Therefore, it is not recommended to use **OSTEOPAL®V**. If no other option is available, the decision whether to use **OSTEOPAL®V** lies with the attending surgeon.

### Side effects

A temporary drop in blood pressure is commonly observed immediately after implantation of the bone cement. In isolated cases, serious complications can develop which may be associated with cardiac arrest, anaphylactic shock or even sudden death.

These cardiovascular and respiratory side effects, also known as implantation syndrome or bone cement syndrome, originate from an influx of bone marrow components into the venous vascular system. In case of pulmonary or cardiovascular complications, it may be necessary to monitor and possibly increase the blood volume. In the presence of acute respiratory failure, anaesthetic management measures must be taken. The following additional adverse effects have been observed as a result of the use of poly(methyl acrylate) bone cements: thrombophlebitis, haemorrhage, trochanteric bursitis. Other side effects observed: myocardial infarction, temporary cardiac arrhythmia, cerebrovascular incident.

Those complications that may occur with any surgical procedure are also possible.

The requirements for a contemporary cementing technique must be considered when using **OSTEOPAL®V** in order to limit undesirable side effects and to safeguard stable and lasting anchoring of the bone cement in the vertebral body.

### Interactions

None reported to date.

### Precautions

#### Use by operating staff

Before using **OSTEOPAL®V**, the user should be well acquainted with its properties, handling and application.

The user is advised to practise the entire procedure of mixing, handling and inserting **OSTEOPAL®V** before using it for the first time. In-depth knowledge is also required if mixing systems and syringes are used to apply the bone cement.

The liquid monomer is highly volatile and flammable. Ignition of monomer fumes caused by use of electrocautery devices in surgical sites near freshly implanted bone cements has been reported. The monomer is also a strong lipid solvent and should not come into direct contact with the body.

When handling the monomer or the prepared bone cement **OSTEOPAL®V**, gloves must be worn that provide the requisite protection against penetration of the skin by the monomer methyl methacrylate. PVP gloves (three-layered polyethylene, ethylene vinyl alcohol copolymer and polyethylene) and Viton®/Butyl gloves are highly recommended over a long period of use. As an additional precaution, it is recommended to wear two pairs of gloves over one

another, e.g. one polyethylene surgical glove over an inner pair of standard latex surgical gloves.  
The use of latex or polystyrene-butadiene gloves alone is insufficient. Please ask your supplier which gloves are suitable for this application.

The monomer vapours can irritate airways and eyes and may damage the liver. Skin irritations have been described that can be attributed to contact with the monomer.  
Manufacturers of soft contact lenses recommend removing the lenses in the presence of damaging or irritant vapours. As soft contact lenses are permeable to liquids and gases, they should not be worn in the operating theatre if methyl methacrylate is used.

#### Use in patients

Blood pressure, pulse and breathing must be monitored carefully during and immediately after introduction of the bone cement. Any significant change in these vital signs must be immediately remedied using appropriate measures. Leakage of bone cement can cause damage to the paravertebral structures. Complications such as spinal cord compression, intercostal neuralgia, cement escaping into the intervertebral space, perivertebral filling of veins and arteries (danger of embolism), infections and post-procedural pain are possible. To prevent bone cement escaping and to detect unwanted occurrences in good time, application must be performed using imaging techniques (real-time imaging). Immediate surgical intervention may also be necessary to counteract the complications described. Prior to surgery, careful radiological examination must be performed to assess possible risks (e.g. vertebral body lesions, vascular supply of the vertebral body or oedema). Incomplete filling of the vertebral body with bone cement may result in an insufficient reduction in acute pain and reduced long-term stability of the treated vertebral body.

#### Incompatibilities

Aqueous solutions (e.g. containing antibiotics) must not be mixed with the bone cement as they have a significant negative effect on the physical and mechanical properties of the bone cement.

#### Dosage

One dose is prepared by mixing the entire contents of a cement powder pouch with all the monomer liquid from a single glass ampoule. The quantity to be applied depends on the particular anatomical conditions. There is no information available on the maximum amount of bone cement that can be introduced and the maximum number of vertebral bodies to be treated for a vertebroplasty or kyphoplasty. Generally, one or two doses are used. However, this depends on the surgical technique used as well as the size of the particular vertebral body and the defect.

At least one additional dose of **OSTEOPAL®V** should be available before starting the operation.


#### Procedure for use

##### Preparation:

Before opening the unsterile protective aluminium pouch, the contents should be moved down to the bottom of the sachet by shaking or tapping so that the contents are not damaged when opening. The polyethylene paper bag (PE paper bag) and the glass ampoule may only be opened

under sterile conditions. The sterile components (inner polyethylene paper bag and glass ampoule) are therefore supplied sterile.

##### Opening under sterile conditions:

The opening flaps from the bag top assist to detach the foil from the paper. In order to grip as much of the opening flaps as possible, the side of the paper/foil should be kept between thumb, index and middle finger. Please use the whole thumb surface to grip foil and paper side and take off each side evenly.

Open the outer PE paper bag at the point provided under sterile conditions so that the inner PE paper bag remains sterile when it is removed. The blister packaging is also opened under sterile conditions and at the point provided so that the glass ampoule can be removed and kept sterile. Before opening the inner polyethylene paper bag, the contents are moved down by shaking or tapping to ensure that there is no loss of cement powder when the top edge is cut open.

To facilitate opening of the glass ampoule, it has a preset breaking point at the transition to the tip of the glass ampoule. Do not open the glass ampoule over the mixing vessel to prevent contamination of the bone cement with glass fragments. A snapping-off device (tube) is attached to the glass ampoule to facilitate opening of the glass ampoule. Rather than the tip of the glass ampoule, the attached snapping-off device is held and the tip of the glass ampoule is broken off using the device. The head of the glass ampoule remains in the snapping-off device after being broken off.

##### Mixing the components:

The liquid from the glass ampoule is placed in a sterile mixing vessel. Then the powder is added from the open inner sachet. The mixture is stirred with a sterile spatula or spoon until a homogeneous paste is produced. The mixture should be stirred for 30 seconds regardless of the ambient temperature. Strict adherence to the instructions for mixing the cement powder and monomer liquid components can minimise complications.

The cement components can also be mixed in a mixing system with or without vacuum. The mixing times for mixing with or without vacuum are also 30 seconds, independent of the ambient temperature. Details on the mixing systems are provided in their instructions for use.

The processing time and polymerisation is greatly dependent on the temperature of the components and the environment. Higher temperatures accelerate and lower temperatures slow down the hardening time.

##### Processing

The viscosity increases with progressing polymerisation, i.e. as the processing phase progresses. The bone cement paste should be placed in an application system immediately after mixing because its viscosity is still low and it is easy to aspirate at this point. To prevent vascular leakage of bone cement, the bone cement should be applied when it is a paste-like consistency. **OSTEOPAL®V** can be introduced

into the vertebral body using an application system that has been approved for percutaneous vertebroplasty or kyphoplasty and enables steady and controlled injection. The handling of the system is described in the manufacturer's instructions for use. During intravertebral application, systematic real-time fluoroscopy (latero-lateral) is necessary. In the event of paravertebral cement leakage, cement injection must be immediately stopped and can be continued after the cement's viscosity has been increased. If the vertebral filling is not sufficient, additional contralateral access is viable. After augmentation, a mandrel should be inserted into the hypodermic needle so that no traces of bone cement remain in the soft tissue after removal of the hypodermic needle.

The mixing, waiting, application and setting times for **OSTEOPAL®V** are shown in the diagrams at the end of the instructions for use. The values are valid for application using a 3-mL single-use syringe and a cannula with a diameter of 3.0 mm and a length of 10 cm. The processing may be different for other application systems. Cannulae with a diameter of less than 1.8 mm (13G) should not be used. The patient must remain immobilised until the bone cement has fully set.

#### **Storage**

**OSTEOPAL®V** must be stored unopened and protected from light at a maximum temperature of 25°C (77°F) in a dry, clean place in the original packaging.

#### **Shelf life / sterility**

The shelf life is printed on the folding box, the protective aluminium pouch and the polyethylene paper bag. Do not use **OSTEOPAL®V** after the expiry date. For production reasons, the shelf lives of the individual components may differ from that specified on the folding box. The contents of opened or damaged protective aluminium pouch or ampoule blisters may not be re-sterilised and must therefore be discarded. If the cement powder has turned yellow, do not use **OSTEOPAL®V**. The cement powder, PE pouch and ampoule blister are sterilised by gassing with ethylene oxide. The monomer liquid has been sterilised by filtration.

#### **Disposal**

Single components of the bone cement, cured solid material as well as (uncleaned) packaging material must be disposed by following the regulations of the local authorities. Dispose of the polymer component in an authorized waste facility. The liquid component should be evaporated under a well-ventilated hood or absorbed by an inert material and transferred in a suitable container for disposal.

**Uses and properties**

**OSTEOPAL®V** is a setting polymer for use in vertebroplasty. Mixing of the two component system, consisting of a powder and a liquid, initially produces a paste, which is injected in the vertebral body for stabilization. Insoluble zirconium(IV) oxide is included in the cement powder as an X-ray contrast medium. The chlorophyll additive serves as optical marking of the bone cement at the site of the operation.

**Indications**

**OSTEOPAL®V** is used for augmentation of vertebral bodies:

- following painful vertebral body compression fractures in osteoporosis
- at painful vertebral body tumors (metastasis or myeloma)
- symptomatic vertebral hemangioma

In all named indications the percutaneous vertebroplasty is only a palliative treatment that stabilizes the vertebral body; a therapy of the systemic illness (osteoporosis, tumor) is not achieved by the percutaneous vertebroplasty.

**Contraindications**

- The use in cases of haemorrhaging diathesis and infections are absolutely contraindicated. Lesions of the vertebral body with epidural extension are relative contraindications (danger of spinal cord compression).
- **OSTEOPAL®V** must not be used in cases of known hypersensitivity to constituents of the bone cement.
- **OSTEOPAL®V** must not be used during pregnancy or the nursing period.
- Allergies or sensitivity to any of its chemical components
- Cases of active or incompletely treated infections
- Coagulation disorders, or with severe cardiopulmonary disease
- Spinal stenosis (> 20 % by retropulsed fragments)

- Compromised of the vertebral fractures due to posterior involvement
- Patient clearly improving on medical therapy
- Prophylaxis in metastatic or osteoporotic patients with no evidence of acute fracture
- Non-pathological, acute traumatic fractures of the vertebra
- Vertebral plana (collapse >90 %)
- Compromise of the vertebral body or the walls of the pedicles
- Unstable vertebral fractures due to posterior involvement

**Target patient population**

Adult population.

There is very little evidence with children and adolescents.

Therefore, it is not recommended to use **OSTEOPAL®V**.

If no other option is available, the decision whether to use **OSTEOPAL®V** lies with the attending surgeon.

**Warnings**

- Read and understand instructions
- Have specific training and experience to be thoroughly familiar with the properties, handling characteristics, and application of the product and percutaneous cement delivery.
- Heraeus Medical as a manufacture does not recommend a surgical technique, it is the responsibility of the physician to determine the appropriateness of this product and specific technique for each patient.
- Hypotensive reactions have occurred between 10 and 165 seconds. They have lasted from 30 seconds to 5 minutes. Some patients have progressed to cardiac arrest. Monitor patients carefully for any change in blood pressure during and immediately following the application of

Composition	OSTEOPAL®V 1x20
	<b>1 pouch of 26.0g powder contains:</b>
poly (methyl acrylate, methyl methacrylate)	14.2g
zirconium dioxide	11.7g
benzoyl peroxide	0.1g
	<b>1 ampoule of 10 ml liquid contains:</b>
methyl methacrylate	9.2g
N,N-dimethyl-p-toluidine	0.2g
<b>Other constituents:</b>	In the powder: chlorophyll VIII In the liquid: chlorophyll VIII in an oily solution, hydroquinone

bone cement. Adverse patient reactions affecting the cardiovascular system have been associated with the use of bone cements.

- MMA can be hypersensitivity in susceptible persons, which may result in anaphylactic response.
- Use in pregnancy and children has not been established.
- Not recommended in patients that do not exhibit a pathological condition, such as primary or secondary osteoporosis or a tumor, which would impair the ability of the patient to heal using conservative treatment methods.
- Strict adherence to good surgical principles and techniques. Deep wound infection is a serious postoperative complication and may require total removal of the embedded cement. Deep wound infection may be latent and not manifest itself even for several years postoperatively.
- Follow handling, mixing and preparation instructions carefully.
- Do not resterilize. Single patient use only. Sterile only if package is unopened and undamaged.
- Caution should be used to prevent excessive exposure to the concentrated vapors of the monomer, which may produce irritation of the respiratory tract, eyes and possible liver.
- Do not use monomer if the liquid component shows any sign of thickening or premature polymerization.
- Do not allow the liquid component to contact rubber or latex gloves. The liquid component is a powerful lipid solvent. Should contact occur, the gloves may dissolve and tissue damage may occur. Wearing a second pair of gloves may diminish the possibility of hypersensitivity reactions.
- Do not allow personnel wearing contact lenses to be near or involved in mixing of the bone cement.
- Use appropriate imaging techniques to verify correct needle placement, absence of damage to surrounding structures, and appropriate location of the injected bone cement. Use imaging, such as fluoroscopy, to assess the ability of the vertebra to contain the injected bone cement.
- Avoid over pressurization of the bone cement because this may lead to leakage of the bone cement beyond the site of its intended application.
- Cement leakage may cause tissue damage, nerve or circulatory problems, and other serious adverse events.
- Leaks can also occur when injecting if the needle is in a vein or if unseen microfracture are prevalent.
- If bone cement is seen outside of the vertebral body or in the circulatory system during the procedure, immediately stop the injection.
- Inadequate fixation or unanticipated postoperative events may affect the cement-bone interface and lead to micro-motion of the cement against bone surface. A fibrous tissue layer may develop between the cement and the bone. Long-term follow-up is advised for all patients on a regular scheduled basis.
- The completion of polymerization occurs in the patient and is an exothermic reaction with considerable liberation of heat. The released heat may damage bone or other tissues surrounding the implant. According to ISO 5833 standard, the temperature can be as high as 95 °C. The long term effect of the heat produced in situ, positioning of the patient should maintained securely to obtain proper

fixation. For proper fixation, 1–2 hours or longer may be required, as determined by the patient's medical condition and the attending physician.

- The long term effects of the bone cement in the spine have not yet been established.

#### Precautions

- Do not use this product after expiration date printed on the package. This device may not be safe or effective beyond its expiration date.
- Follow the handling and mixing instructions carefully to avoid contact dermatitis. Strict adherence to the instructions for mixing the powder and liquid components may reduce the incidence of this complication.
- Adequately ventilate the operating room to eliminate as much monomer vapor as possible. The liquid monomer is highly volatile and flammable. Ignition of monomer fumes caused by the use of electrocautery devices in surgical sites near freshly implanted bone cements has been reported.
- When working with the monomer or the cement, gloves must be worn to assure adequate protection against the penetration of the monomer (methyl methacrylate) into the skin. PVP (three-layer polyethylene, ethylene-vinyl alcohol-copolymer, polyethylene) and Viton/Butyl gloves have proved to provide good protection over an extended period. Putting on two pairs of gloves – one pair of polyethylene surgical gloves over a pair of standard latex surgical gloves has also proved to offer adequate protection. The use of latex or polystyrene-butadiene gloves, however, must be avoided. Please request the confirmation of your glove supplier whether the respective gloves are suitable for the use of this cement.

#### Potential adverse events associated with vertebroplasty procedure include

- Pneumonia
- Intercostal neuralgia
- Collapse of a vertebra adjacent to the one injected, due to an osteoporotic disease
- Pneumothorax
- Cement extrusion into soft tissue
- Fracture of a pedicle
- Rib fracture in patients with diffuse osteopenia, especially during thoracic vertebroplasty procedures, due to the significant downward force exerted during needle insertion
- Compression of the spinal cord with paralysis or loss of feeling
- Cement leakage between intervertebral discs.
- Serious adverse events, some with fatal outcome, associated with the use of acrylic bone cements for vertebroplasty or kyphoplasty include myocardial infarction, cardiac arrest, cerebrovascular accident, pulmonary embolism, and cardiac embolism. Although the majority of these adverse events present early within the post-operative period, there have been some reports of diagnoses beyond a year or more after the procedure. Other reported adverse events for acrylic bone cements intended for vertebroplasty or kyphoplasty include: Leakage of the bone cement beyond the site of its intended application with introduction into the vascular

system resulting in embolism of the lung and/or heart or other clinical sequelae.

- The most frequent adverse reactions reported with acrylic bone cements are a transitory fall in blood pressure, thrombophlebitis, hemorrhage and hematoma, loosening or displacement of the prosthesis, superficial or deep wound infection, trochanteric bursitis, and short-term cardiac conduction irregularities. Other reported adverse events for acrylic bone cements include pyrexia due to an allergy to the bone cement, hematuria, dysuria, bladder fistula, and delayed sciatic nerve entrapment due to extrusion of the bone cement beyond the region of its intended application.

#### Important Physician Information

- Percutaneous vertebroplasty procedures should only be performed in medical settings in which emergency decompressive surgery is available.
- Adverse reactions affecting the cardiovascular system have been attributed to recent data that indicate the monomer undergoes rapid hydrolysis to methacrylic acid, and that a significant fraction of the circulating methacrylate is in the form of circulating concentrations of methyl methacrylate/methacrylic acid and changes in blood pressure has not been established.
- The physician is responsible for any complication or harmful consequences, which may result from an erroneous indication or operating technique, from inappropriate use of the material, or non-observation of the safety instructions that appear in the directions for use.
- Additives (such as antibiotics) are not to be mixed with the bone cement, as this will alter cement properties.

#### Interactions

Not known.

#### Incompatibilities

Aqueous (e.g. antibiotic-containing) solutions must not be mixed with the bone cement, as this reduces the strength considerably.

#### Dosage

The amount of bone cement required depends on the patient's anatomy. There is no study about the maximal amount of bone cement or the maximal amount of treatable vertebral bodies by vertebroplasty. This is subject to the surgeon's responsibility. One or more complete units (the contents of one sachet and one ampule) must always be mixed together.


#### Notes on use

Before using **OSTEOPAL® V** the surgeon must be thoroughly familiar with its properties, handling and use and must have read the relevant literature. In particular, the surgeon must be proficient in the special application techniques required for low viscosity cements. An appropriate up to date cementing technique must be used with **OSTEOPAL® V** to limit the side effects and to achieve a stable and long-lasting anchorage of the cement in the vertebral body. The application must be performed with an application instrument certified for vertebroplasty. For special application techniques, the relevant instructions of the application system must be consulted by the surgeon. In order to recog-

nize immediately cement leakage, an imaging procedure (real time) during the application of the cement is recommended.

For recognition of a vertebral extension a phlebography before injection of the cement is recommended. The protective aluminium packaging and outer non-sterile polyethylene sachet („peel-off“ package) and the ampule packaging should be opened maintaining sterility by the circulating nurse. The ampule and powder sachet are then taken out under sterile conditions and placed on a sterile working surface. The ampule is opened by breaking the neck and the inner sachet is cut open with sterile scissors.

#### Opening under sterile conditions:

The opening flaps from the bag top assist to detach the foil from the paper.

In order to grip as much of the opening flaps as possible, the side of the paper/foil should be kept between thumb, index and middle finger.

Please use the whole thumb surface to grip foil and paper side and take off each side evenly.

#### Mixing

1. Pour the liquid into a sterile mixing bowl and add all the powder. Stir carefully with a sterile stirring rod until a homogeneous mass is obtained. The mixture should not be stirred for longer than 30 sec, irrespective of the room temperature.
2. Mix the liquid and powder components for a minimum of 1 minute and until the power is completely saturated with the liquid. The liquid monomer and the powder component should be vigorously mixed to create a homogenous mixture, which is essential to obtain consistent cement. Follow the instructions provided by the mixing device manufacturer.
3. Care must be taken to minimize air entrapment within the cement when using a cement delivery system.
4. The physician should, by specific training and experience, be thoroughly familiar with the properties, handling characteristics, application of the bone cements, and amount of cement required per procedure. Because the handling and curing characteristics of this bone cement may vary with temperature, humidity, and mixing technique, they are best determined by the physician's actual experience.
5. Do not apply the cement until two minutes after completely mixing the liquid component to the powder component.
6. Select an appropriate method to apply the cement to bone. If applicable, follow the device manufacturer's instructions.
7. The injection time and set time of radiopaque bone cement will vary with temperature.

#### Working procedure

The working time and the rate of polymerization are strongly dependent on ambient temperature and on temperature of the components. The hardening time is shorter at high temperatures and longer at low temperatures. Viscosity increases as polymerization progresses, i.e. as the cement hardens. After mixing, the cement is poured directly into a

cement syringe, since at this time **OSTEOPAL®V** still has a low viscosity and is easily pourable.

The application of the cement should be performed during a pasty condition phase of the cement to avoid vascular migration. **OSTEOPAL®V** can be applied with a percutaneous vertebroplasty certified cement syringe into the vertebral body which allows a controlled constant injection. The method of application is given in the instructions for the system used by the manufacturer and must be consulted. During intravertebral application a strict latero-lateral real time radioscopy is necessary. In case of paravertebral leakage the injection must be stopped immediately and can be proceeded when the viscosity of the cement is increased. In cases of insufficient vertebral filling a contralateral approach can be performed. After vertebral augmentation a mandrin should be placed into the vertebroplasty injection needle to avoid cement residues being deposited in the soft tissues on withdrawal of the application needle.

At a room and a material temperature of 20 °C (68 °F) the schedule of the temperature/time diagram applies:

The reference data belong to the CementoSet (Optimed GmbH, 76275 Ettlingen, Germany) with a canula diameter of 3.2mm. This time may vary depending on the syringe system used. No canula with a diameter less than 1.8mm should be used. The patients should be kept immobilised until the end of the setting time as illustrated in the graph.

#### **Other uses**

**OSTEOPAL®V** is not certified for other indications than percutaneous vertebroplasty. The application in another indication (i.e. endoprosthesis) is not investigated.

#### **Storage**

Do not store above 25 °C (77 °F).

#### **Shelf life/sterility**

The expiry date is printed on the outer carton, on the protective aluminium packaging and on the inner sachet.

**OSTEOPAL®V** must not be used after the expiry date. Once removed from the protective packaging, the powder sachet cannot be resterilized and must therefore be discarded.

**OSTEOPAL®V** is sterilized with ethylene oxide gas and must not be resterilized. If the polymer powder shows a yellow discoloration, it must not be used.

#### **Disposal**

Single components of the bone cement, cured solid material as well as (uncleaned) packaging material must be disposed by following the regulations of the local authorities.

Dispose of the polymer component in an authorized waste facility. The liquid component should be evaporated under a well-ventilated hood or absorbed by an inert material and transferred in a suitable container for disposal.

## Properties – Composition

**OSTEOPAL®V** is a radiopaque, low-viscosity, quicksetting bone cement for filling and stabilizing vertebral bodies. It is obtained by mixing a polymer powder component with a liquid monomer component. Zirconium oxide has been added to the cement powder as an X-ray contrast medium. The sterile-filtrated monomer component is supplied in an amber glass ampoule and comes in a sterile blister pack. The polymer powder component is supplied in a double sterile packaging. The inner polyethylene sachet which contains the powder component is wrapped in an additional polyethylene sachet; both sachets were sterilized with ethylene oxide. The polyethylene sachets are contained in a non-sterile protective aluminium packaging. Chlorophyll has been used to obtain the green color of **OSTEOPAL®V** in order to ensure clear visibility of the cement at the operating site.

After mixing, a low-viscosity dough is obtained which is filled into the vertebral bodies as a stabilization medium and allowed to harden.

## Indications

**OSTEOPAL®V** is used for augmentation of vertebral bodies

- for painful compression fractures of the vertebral body
  - for painful tumors in the vertebral body (metastatic carcinomas or myelomas)
  - for symptomatic vertebral hemangiomas
- Percutaneous vertebroplasty is only a palliative treatment for stabilizing the vertebral bodies. It does not treat the underlying illness (osteoporosis, tumor-related illness).

## Contraindications

Hemorrhagic diatheses and infections are absolute contraindications. Lesions of the vertebral body with epidural extension are relative contraindications because of the danger of spinal cord compression.

**OSTEOPAL®V** must not be used in cases of known hypersensitivity to the constituents of the bone cement, nor during pregnancy or nursing.

## Target patient population

Adult population.

There is very little evidence with children and adolescents. Therefore, it is not recommended to use **OSTEOPAL®V**. If no other option is available, the decision whether to use **OSTEOPAL®V** lies with the attending surgeon.

## Warning information – Side effects

Prior to using **OSTEOPAL®V** the surgeon should be familiar with its properties, handling and application. It is also recommended for surgeons to practice mixing, handling and application prior to the use. Precise knowledge is also required, if mixing systems and syringes are used for the application of the cement. The liquid monomer is highly volatile and flammable. Ignition of monomer fumes caused by use of electrocautery devices in surgical sites near freshly implanted bone cements has been reported. The monomer is also a powerful lipid solvent and should not come into direct contact with the body. When working with the monomer or the cement, gloves must be worn to ensure adequate protection against the penetration of the monomer

Composition	OSTEOPAL®V 1x20
	<b>1 sachet of 26 g powder contains:</b>
poly (methyl acrylate, methyl methacrylate)	14.2 g
zirconium dioxide	11.7 g
benzoyl peroxide	0.1 g
	<b>1 ampoule of 10 ml liquid contains:</b>
methyl methacrylate	9.2 g
N,N-dimethyl-p-toluidine	0.2 g
<b>Other constituents:</b>	In the powder: colourant E141 In the liquid: colourant E141, hydroquinone

(methyl methacrylate) into the skin. PVP (three-layer polyethylene, ethylene-vinyl alcohol-copolymer, polyethylene) and Viton/Butyl gloves have proved to provide good protection over an extended period. Putting on two pairs of gloves has also proved to offer adequate protection.

The use of latex or polystyrene-butadiene gloves alone, however, must be avoided. Please request the confirmation of your glove supplier whether the respective gloves are suitable for the use with **OSTEOPAL®V**.

The monomer vapors may irritate the respiratory tract and the eyes and possibly damage the liver. Skin irritations have been reported, which must be attributed to the contact with the monomer.

Manufacturers of soft contact lenses recommend the removal of the lenses in the presence of harmful or irritating vapors. Since contact lenses are permeable to liquids and gases, they should not be worn in the operating room if methyl methacrylate is used.

### Precautionary measures for application

Directly after implanting the bone cement into the vertebral body, the increase in pressure in the vertebral bodies can lead to a temporary drop in blood pressure. In rare cases, the hypotension may be accompanied by pulmonary embolism, cardiac arrest and death. These cardiovascular and respiratory side effects, also known as implantation syndrome or bone cement syndrome, result from an influx of bone marrow constituents into the venous vascular system.

If pulmonary, cardiovascular complications arise, monitoring and – in some cases – even increasing the blood volume may be required. In cases of acute respiratory insufficiency, anaesthesiological measures should be taken.

In addition, paravertebral structures may be damaged by cement escaping. Complications such as spinal cord compression, intercostal neuralgia, cement escaping into the intervertebral space, perivertebral filling of veins and arteries (danger of cement embolism to heart and lungs), infections and post-procedural pain are possible. Some evidence suggests that release of monomer into the circulation could also be a factor. These serious adverse events appear to be more frequently associated with multiple level vertebroplasty/ kyphoplasty procedures. High quality real time imaging (Biplanar fluoroscopy or CT Guidance) is needed to monitor the injection of the cement, and to minimize extravertebral escape of the cement.

In order to detect a vertebral extension, a phlebography must be performed before applying the cement. Immediate surgical intervention may also be necessary to counteract the described complications. Prior to surgery, a careful radiological investigation must be performed to assess possible risks (e.g. vertebral-body lesions, vascular supply of the vertebral body or edema). When treating hemangiomas, preoperative angiosterosclerosis through percutaneous application of alcohol may help prevent cement entering the vascular structures. Incomplete filling of the vertebral body with bone cement may result in an insufficient reduction in acute pain and reduced long-term stability of the treated vertebral body.

### Interactions

No known interactions.

### Incompatibilities

Aqueous (e.g containing antibiotics) solutions must not be added to the bone cement since they considerably impair the mechanical properties of the cement.

### Dose

The amount of bone cement required depends on the particular anatomical conditions. No information is available on the maximum amount of bone cement that can be used and the maximum number of vertebral bodies to be treated by a vertebroplasty procedure. One or more units (sachet and ampoule) must always be mixed together.

### Notes on application

All patients with osteoporotic compression fractures should have had an adequate trial of conservative management before considering vertebroplasty or kyphoplasty. Only physicians trained to do vertebroplasty and kyphoplasty should be performing these procedures.

When applying **OSTEOPAL®V**, the requirements of current cement techniques must be observed in order to prevent unwanted side effects and to ensure stable, long-term fixation of the bone cement in the vertebral body. Application must be performed with an instrument authorized for filling vertebral bodies, and the operator must be familiar with the use of this instrument.

### Preparation

A single dose is prepared by mixing the entire content of a powder sachet with one ampoule. In most cases, one or two portions suffice, although this depends on the particular surgical operation and the technique used. Before the beginning of the operation, at least one additional dose of **OSTEOPAL®V** should be readily available. Each dose is prepared separately.

The following is required to prepare the bone cement: Sterile working surface, porcelain or stainless-steel bowls, sterile mixing spoons or porcelain or stainless-steel spatulas. The protective aluminium packaging, the outer nonsterile polyethylene sachet and the blister pack of the ampoule should be opened by an assistant in a way that maintains sterility. The sterile polyethylene sachet and the ampoule are placed under aseptic conditions on a sterile desk. Sterile conditions must be ensured when opening the polyethylene sachet and the ampoule.

### Opening under sterile conditions:

The opening flaps from the bag top assist to detach the foil from the paper.

In order to grip as much of the opening flaps as possible, the side of the opening paper/foil should be kept between thumb, index and middle finger.

Please use the whole thumb surface to grip foil and paper side and take off each side evenly.

### Processing

Processing time and polymerization are heavily dependent on the temperature of the components and the ambient conditions. High temperatures accelerate hardening time,

low temperatures prolong it. Viscosity increases as polymerization progresses, i.e. as the processing phase progresses. The cement mixture is drawn into a syringe as soon as it is mixed, as it is then still low-viscosity and easy to inject. To prevent a vascular escape of cement, it should be applied when soft. **OSTEOPAL®V** can be injected into the vertebral body using a cement syringe authorized for percutaneous vertebroplasty which allows constant and controlled injection. For instructions on using the cement syringe, users should refer to the manufacturer's instructions. Intravertebral application must be supported by fluoroscopy (laterolateral) with real-time display. In the event of paravertebral cement escape, cement injection must be stopped immediately and can be continued after increasing the viscosity of the cement. If vertebral filling is insufficient, further contralateral access is possible. After augmentation, a mandarin should be inserted in the injection needle, to prevent any cement residue remaining in the soft tissues after the injection needle is removed. The processing times and their temperature dependency are shown in the table. These values apply when using CementoSet (Optimed GmbH, 76275 Ettlingen, Germany) with an injection cannula of diameter 3.2 mm. For other syringe systems, different processing methods may be used. Cannulas with a diameter less than 1.8 mm should not be used. The patient must be immobilized until the cement has completely hardened.

#### **Storage**

Do not store above 25 °C.

#### **Shelf life/Sterility**

The expiration date is printed on the folded box, the protective aluminium packaging and the inner sachet. **OSTEOPAL®V** must not be used after the expiration date. The contents of unused but opened or damaged packs must be discarded and must not be resterilized. **OSTEOPAL®V** has been sterilized with ethylene oxide gas and must not be resterilized. The polymer powder must not be used if it exhibits yellow.

#### **Disposal**

Single components of the bone cement, cured solid material as well as (uncleaned) packaging material must be disposed by following the regulations of the local authorities. Dispose of the polymer component in an authorized waste facility. The liquid component should be evaporated under a well-ventilated hood or absorbed by an inert material and transferred in a suitable container for disposal.

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Heraeus Medical GmbH