Clinical Abstract

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Literatur

1. The Orthopaedic Industry Annual Report. Orthoworld Inc. 2010.
2. Hart JA. Joint replacement surgery. *Med J Aust* 2004; 180:S27–30.
3. Coventry MB, Beckenbaugh RD, Nolan DR, Ilstrup DM. 2,012 total hip arthroplasties. A study of postoperative course and early complications. *J Bone Joint Surg Am* 1974; 56:273–284.
4. Ereth MH, Weber JG, Abel MD, et al. Cemented versus noncemented total hip arthroplasty-embolism, hemodynamics, and intrapulmonary shunting. *Mayo Clin Proc* 1992; 67:1066–1074.
5. Parvizi J, Holiday AD, Ereth MH, Lewallen DG. The Frank Stinchfield Award. Sudden death during primary hip arthroplasty. *Clin Orthop Relat Res* 1999; 369:39–48.
6. Deitelzweig SB, McKean SC, Amin AN, Brotman DJ, Jaffer AK, Spyropoulos AC. Prevention of venous thromboembolism in the orthopedic surgery patient. *Cleve Clin J Med* 2008; 75 Suppl 3:S27–36.
7. Freedman KB, Brookenthal KR, Fitzgerald RH, Jr., Williams S, Lonner JH. A meta-analysis of thromboembolic prophylaxis following elective total hip arthroplasty. *J Bone Joint Surg Am* 2000; 82-A:929–938.
8. Geerts WH, Pineo GF, Heit JA, et al. Prevention of venous thromboembolism: the Seventh ACCP Conference on Antithrombotic and Thrombolytic Therapy. *Chest* 2004; 126:338S–400S.
9. Cohen AT, Tapson VF, Bergmann JF, et al. Venous thromboembolism risk and prophylaxis in the acute hospital care setting (ENDORSE study): a multinational cross-sectional study. *Lancet* 2008; 371:387–394.
10. Stratton MA, Anderson FA, Bussey HI, et al. Prevention of venous thromboembolism: adherence to the 1995 American College of Chest Physicians consensus guidelines for surgical patients. *Arch Intern Med* 2000; 160:334–340.
11. Kalka C, Spirk D, Siebenrock KA, et al. Lack of extended venous thromboembolism prophylaxis in high-risk patients undergoing major orthopaedic or major cancer surgery. Electronic Assessment of VTE Prophylaxis in High-Risk Surgical Patients at Discharge from Swiss Hospitals (ESSENTIAL). *Thromb Haemost* 2009; 102:56–61.
12. Darouiche RO. Treatment of infections associated with surgical implants. *N Engl J Med* 2004; 350:1422–1429.
13. Trampuz A, Zimmerli W. Prosthetic joint infections update in diagnosis and treatment. *Swiss Medical Weekly* 2005; 135:243–251.
14. Schierholz JM, Morsczech C, Brenner N, et al. [Special aspects of implant-associated infection in orthopedic surgery. From the pathophysiology to custom-tailored prevention strategies]. *Orthopade* 2004; 33:397–404.
15. Nasser S. The incidence of sepsis after total hip replacement arthroplasty. *Semin Arthroplasty* 1994; 5:153–159.
16. Donaldson AJ, Thomson HE, Harper NJ, Kenny NW. Bone cement implantation syndrome. *Br J Anaesth* 2009; 102:12–22.
17. Weiler-Lorentz A. Komplikationen. In: Frietsch T, Weiler-Lorentz A, eds. *Anästhesie in der Traumatologie und Orthopädie*. München: Urban & Fischer/Elsevier, 2009; 92–119.
18. Lafont ND, Kalonji MK, Barre J, Guillaume C, Boogaerts JG. Clinical features and echocardiography of embolism during cemented hip arthroplasty. *Can J Anaesth* 1997; 44:112–117.
19. Peebles DJ, Ellis RH, Stride SD, Simpson BR. Cardiovascular effects of methylmethacrylate cement. *Br Med J* 1972; 1:349–351.
20. Modig J, Busch C, Waernbaum G. Effects of graded infusions of monomethylmethacrylate on coagulation, blood lipids, respiration and circulation. An experimental study in dogs. *Clin Orthop Relat Res* 1975; 113:187–197.
21. Modig J, Busch C, Olerud S, Saldeen T, Waernbaum G. Arterial hypotension and hypoxaemia during total hip replacement: the importance of thromboplastic products, fat embolism and acrylic monomers. *Acta Anaesthesiol Scand* 1975; 19:28–43.
22. Orsini EC, Byrick RJ, Mullen JB, Kay JC, Waddell JP. Cardiopulmonary function and pulmonary microemboli during arthroplasty using cemented or non-cemented components. The role of intramedullary pressure. *J Bone Joint Surg Am* 1987; 69:822–832.
23. Clarius M, Heisel C, Breusch SJ. Pulmonary embolism in cemented total hip arthroplasty. In: Breusch SJ, Malchau H, eds. *The well-cemented total hip arthroplasty*. Berlin: Springer, 2005; 320–331.
24. Breusch SJ, Lukoschek M, Schneider U, Ewerbeck V. „State of the art“ der zementierten Hüftendoprothetik. *Dt Arztebl* 2000; 97:A2030–2033.
25. Clarius M, Heisel C, Breusch SJ (2005) Pulmonary embolism in cemented total hip arthroplasty. In: Breusch SJ, Malchau H, eds. *The well-cemented total hip arthroplasty*. Berlin: Springer, 2005. ISBN 978-3-540-24197-3; Fig. 15.9 a,b auf Seite 326.
26. Krebs J, Ferguson SJ, Nuss K, et al. Plasma levels of endothelin-1 after a pulmonary embolism of bone marrow fat. *Acta Anaesthesiol Scand* 2007; 51:1107–1114.
27. Tryba M, Linde I, Voshage G, Zenz M. [Histamine release and cardiovascular reactions to implantation of bone cement during total hip replacement]. *Anaesthesist* 1991; 40:25–32.
28. Bengtson A, Larsson M, Gammer W, Heideman M. Anaphylatoxin release in association with methylmethacrylate fixation of hip prostheses. *J Bone Joint Surg Am* 1987; 69:46–49.
29. Gammer W, Bengtson A, Heideman M. Inhibition of complement activation by high-dose corticosteroids in total hip arthroplasty. *Clin Orthop Relat Res* 1988; 236:205–209.
30. Lewis RN. Some studies of the complement system during total hip replacement using bone cement under general anaesthesia. *Eur J Anaesthesiol* 1997; 14:35–39.
31. Lamade WR, Friedl W, Schmid B, Meeder PJ. Bone cement implantation syndrome. A prospective randomised trial for use of antihistamine blockade. *Arch Orthop Trauma Surg* 1995; 114:335–339.
32. Mitsuhashi H, Saitoh J, Saitoh K, et al. Methylmethacrylate bone cement does not release histamine in patients undergoing prosthetic replacement of the femoral head. *Br J Anaesth* 1994; 73:779–781.
33. Rinecker H. New clinico-pathophysiological studies on the bone cement implantation syndrome. *Arch Orthop Trauma Surg* 1980; 97:263–274.
34. Frost PM. Cardiac arrest and bone cement. *Br Med J* 1970; 3:524.
35. Heisel C, Clarius M, Schneider U, Breusch SJ. [Thromboembolic complications related to the use of bone cement in hip arthroplasty-pathogenesis and prophylaxis]. *Z Orthop Ihre Grenzgeb* 2001; 139:221–228.
36. Breed AL. Experimental production of vascular hypotension, and bone marrow and fat embolism with methylmethacrylate cement. *Traumatic hypertension of bone*. *Clin Orthop Relat Res* 1974; 102:227–244.
37. Pietak S, Holmes J, Matthews R, Petrusek A, Porter B. Cardiovascular collapse after femoral prosthesis surgery for acute hip fracture. *Can J Anaesth* 1997; 44:198–201.
38. Koessler MJ, Pitto RP. Fat and bone marrow embolism in total hip arthroplasty. *Acta Orthop Belg* 2001; 67:97–109.
39. Herrenbruck T, Erickson EW, Damron TA, Heiner J. Adverse clinical events during cemented long-stem femoral arthroplasty. *Clin Orthop Relat Res* 2002; 395:154–163.
40. Criteria Committee of the New York Heart Association. Nomenclature and criteria for diagnosis of diseases of the heart and great vessels. Boston: Little, Brown & Co, 1994.
41. McCaskie A. Femoral pressurisation. In: Breusch SJ, Malchau H, eds. *The well-cemented total hip arthroplasty*. Berlin: Springer, 2000; 160–163.
42. Patterson BM, Healey JH, Cornell CN, Sharrock NE. Cardiac arrest during hip arthroplasty with a cemented long-stem component. A report of seven cases. *J Bone Joint Surg Am* 1991; 73:271–277.
43. Guest CB, Byrick RJ, Mazer CD, Wigglesworth DF, Mullen JB, Tong JH. Choice of anaesthetic regimen influences haemodynamic response to cemented arthroplasty. *Can J Anaesth* 1995; 42:928–936.
44. Werner L. A response to ‚Sudden deaths during hip hemi-arthroplasty‘, Parry G, *Anaesthesia* 2003; 58: 922-23. *Anaesthesia* 2004; 59:200.
45. Byrick RJ, Bell RS, Kay JC, Waddell JP, Mullen JB. High-volume, high-pressure pulsatile lavage during cemented arthroplasty. *J Bone Joint Surg Am* 1989; 71:1331–1336.

46. Wheelwright EF, Byrick RJ, Wigglesworth DF, et al. Hypotension during cemented arthroplasty. Relationship to cardiac output and fat embolism. *J Bone Joint Surg Br* 1993; 75: 715–723.
47. Huo MH, Brown BS. What's new in hip arthroplasty. *J Bone Joint Surg Am* 2003; 85-A: 1852–1864.
48. Evans RD, Palazzo MG, Ackers JW. Air embolism during total hip replacement: comparison of two surgical techniques. *Br J Anaesth* 1989; 62: 243–247.
49. McCaskie AW, Barnes MR, Lin E, Harper WM, Gregg PJ. Cement pressurisation during hip replacement. *J Bone Joint Surg Br* 1997; 79: 379–384.