
Heraeus

Medical

# PALABOWL®

Bone Cement Mixing System

Rx Only

# PALABOWL®

<b>EN</b>	<b>Instructions for use</b> (Revision status: 2023-09)	<b>4</b>
<b>DE</b>	<b>Gebrauchsanweisung</b> (Revisionsstand: 2023-09)	<b>6</b>
<b>EL</b>	<b>Οδηγίες χρήσης</b> (Ημερομηνία αναθεώρησης: 2023-09)	<b>8</b>
<b>ES</b>	<b>Instrucciones de uso</b> (Última revisión: 2023-09)	<b>10</b>
<b>FR</b>	<b>Notice d'utilisation</b> (Date de révision: 2023-09)	<b>12</b>
<b>IT</b>	<b>Istruzioni per l'uso</b> (Ultima revisione: 2023-09)	<b>14</b>
<b>NL</b>	<b>Gebruiksaanwijzing</b> (Revisiestatus: 2023-09)	<b>16</b>
<b>PT</b>	<b>Instruções de utilização</b> (Data da última actualização: 2023-09)	<b>18</b>

EN

# PALABOWL®

## 1. Medical Device

- 1) Device name: **PALABOWL® Bone Cement Mixing System**
- 2) Model No.: **CL006-01**
- 3) REF: 5050166

## 2. Indications (Intended Use)

- 1) The Bone Cement Mixing System is for mixing of the bone cement powder and liquid monomer.
- 2) This device can mix up to three units (total 120g powder) of bone cement of any viscosity.

## 3. Working principles

This device's rotating axis ensures the mixing paddle reaches all angles inside the bowl and hence creates a uniform bone cement mixture. Air, fume and odors are largely removed from the mixing bowl through a vacuum hose to reduce the air bubbles and porosity in the cement mix, and limits exposure of staff to harsh odor and toxic fume.

## 4. Specification

- 1) Mixing bowl size (ID×H): About 106 mm×72 mm;  
Vacuum hose size (length): About 2.4 m.
- 2) Net weight: Not more than 0.3 kg.
- 3) Air tightness: There is no leak under negative pressure 75 kPa±1k Pa.
- 4) Vacuum hose: The vacuum hose shall not shrink under negative pressure 75 kPa±1 kPa.
- 5) Mixing paddle: The Mixing paddle can rotate freely under negative pressure 75 kPa±1 kPa.
- 6) Sterilization: Sterilized using ethylene oxide before delivery.
- 7) Valid period (Shelf life): Three years.
- 8) The maximum allowable negative pressure value during operating: 73 kPa (550 mmHg).
- 9) To reach the negative pressure value 73 kPa (550 mmHg) needed time: about 25s.

## 5. Components


- 1) This device consists of **Vacuum Mixing bowl, Mixing paddle, Spatula, Curette, Vacuum hose** and **Rotating handle**.
- 2) Associated devices: The **Bone Cement** and **Vacuum extractor** are not included in this device.


## 6. Warning and Precautions

- 1) This device shall only be operated by qualified professionals.
- 2) The bone cement powder and liquid monomer are not included in this device but prepared by the hospital. The bone cement must be prepared according to the instructions for use provided by the bone cement manufacturer.
- 3) The Vacuum extractor is not included in this device, but prepared by the hospital. The negative pressure within the range of 5 kPa to 75 kPa shall be adjustable.
- 4) This device is sterilized using ethylene oxide before delivery. Sterilization validity is three years.
- 5) This device is for single patient use only. Do not reuse and do not re-sterilize.
- 6) Follow current local/national regulations to properly dispose the device after use.
- 7) During transportation of this device, prevent severe shock, vibration, rain/snow, direct sunlight, moist and pressure.
- 8) This device should be stored in a non-corrosive, cool, dry, and ventilated environment.


## 7. Operating instructions

- 1) Inspect the package prior to each use. Do not use the device if package is damaged.
- 2) Read the instructions for use, attached documents and labels.
- 3) Open the package, take out the product and operate the device according to following figures:



1. Place the bowl on a flat surface and open the lid.


2. Connect the shorter end of the vacuum hose to the vacuum pump and the opposite end (longer end) to the lid. If conditions permit, adjust negative pressure to the optimal range between 65 kPa to 73 kPa.


3. Pour the single components into the bowl according to the instructions provided by the cement manufacturer.


4. Close the lid and twist clockwise to lock.


5. Start the vacuum and begin mixing for about 45 seconds or according to the instructions for use of the bone cement.


6. Once mixing is complete, stop the vacuum and remove the vacuum hose.

7. Turn the lid counterclockwise and remove the lid.

8. To scrape the bone cement on the mixing paddle into mixing bowl by using spatula or curette.

9. Remove the thoroughly mixed bone cement with the spatula.

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# Heraeus

## Symbol Descriptions

The following symbols are applicable to this device, and will appear on instructions for use, labels and packages.

	Batch code (Lot number)		
	Manufacturer		
	Sterilized using ethylene oxide		
	Do not re-use		
	Do not use if package is damaged and consult instructions for use		
	Date of manufacture		
	Use-by date		Caution
	Consult instructions for use		Authorized representative in the European Community/European Union
	Catalogue number		Model number
	Federal law restricts this device to sale, distribution, and use by or on the order of a physician		Unique device identifier
	Keep dry		Importer
	Do not re-sterilize		Distributor

# CE 0123

	Medical device
	Keep away from sunlight
	Fragile, handle with care
	This way up

	<b>SUNGO Europe B.V.</b> Fascinatio Boulevard 522, Unit 1.7, 2909VA Capelle aan den IJssel, The Netherlands. Tel: +31(0)10 3034500 E-mail: ec.rep@sungogroup.com
--	--

**USA Agent:**  
Regrek LLC  
19 Holly Cove Ln., Dover, Delaware, 19901, USA

CL006-01-08E(Z3), Rev.: I, 2023-09-04

**Heraeus Medical GmbH**  
Philipp-Reis-Straße 8/13, 61273 Wehrheim, Germany  
**Heraeus Medical LLC**  
770 Township Line Road, Yardley, PA 19067, USA

**Manufacturer:**  
**Guangzhou Clean Medical Products Manufacturing Corp.**  
No.163 Meidu Road, Chengjiao, Conghua District,  
Guangzhou, Guangdong, China. P.C.: 510900  
Tel: +86-20-87654060, 87817888  
E-mail: info@cleanmfg.com  
Web: www.cleanmfg.com