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Doc.-nr.: 61342

Pagina 13 van
18**Referenties**

1. Wahlig H, Dingeldein E, Buchholz HW, Buchholz M, Bachmann F. Pharmacokinetic study of gentamicin-loaded cement in total hip replacements. Comparative effects of varying dosage. *J Bone Joint Surg Br.* 1984;66(2):175-9.
2. Trela-Larsen L, Sayers A, Blom AW, Webb JCJ, Whitehouse MR. The association between cement type and the subsequent risk of revision surgery in primary total hip replacement. *Acta Orthop.* 2018;89(1):40-6.
3. van der Voort P, Valstar ER, Kaptein BL, Fiocco M, van der Heide HJL, Nelissen RGHH. Comparison of femoral component migration between Refobacin bone cement R and Palacos R + G in cemented total hip arthroplasty: A randomised controlled roentgen stereophotogrammetric analysis and clinical study. *Bone Joint J.* 2016;98-B(10):1333 – 41.
4. Jørgensen PB, Lamm M, Søballe K, Stilling M. Equivalent hip stem fixation by Hi-Fatigue G and Palacos R + G bone cement: a randomized radiostereometric controlled trial of 52 patients with 2 years' follow-up. *Acta Orthop.* 2019;90(3):237 – 42.
5. Tabori-Jensen S, Mosegaard SB, Hansen TB, Stilling M. Inferior stabilization of cementless compared with cemented dual-mobility cups in elderly osteoarthritis patients: a randomized controlled radiostereometry study on 60 patients with 2 years' follow-up. *Acta Orthop.* 2020;91(3):246-53.
6. Turnbull GS, Marshall C, Nicholson JA, MacDonald DJ, Clement ND, Breusch SJ. The Olympia anatomic polished cemented stem is associated with a high survivorship, excellent hip-specific functional outcome, and high satisfaction levels: follow-up of 239 consecutive patients beyond 15 years. *Arch Orthop Trauma Surg.* 2021.
7. Siepen W, Zwicky L, Stoffel KK, Ilchmann T, Clauss M. Prospective two-year subsidence analysis of 100 cemented polished straight stems - a short-term clinical and radiological observation. *BMC Musculoskelet Disord.* 2016;17(1):395.
8. Jain S, Magra M, Dube B, Veysi VT, Whitwell GS, Aderinto JB, Emerton ME, Stone MH, Pandit HG. Reverse hybrid total hip arthroplasty. *Bone Joint J.* 2018;100-B(8):1010 – 1017.
9. Lindalen E, Dahl J, Nordsletten L, Snorrason F, Høvik, Röhrli S. Reverse hybrid and cemented hip replacement compared using radiostereometry and dual-energy X-ray absorptiometry: 43 hips followed for 2 years in a prospective trial. *Acta Orthop.* 2012;83(6):592 – 8.
10. Sprowson AP, Jensen C, Chambers S, Parsons NR, Aradhyula NM, Carluke I, Inman D, Reed MR. The use of high-dose dual-impregnated antibiotic-laden cement with hemiarthroplasty for the treatment of a fracture of the hip: The Fractured Hip Infection trial. *Bone Joint J.* 2016;98-B(11):1534 – 41.
11. Tyas B, Marsh M, Oswald T, Refaie R, Molyneux C, Reed M. Antibiotic resistance profiles of deep surgical site infections in hip hemiarthroplasty; comparing low dose single antibiotic versus high dose dual antibiotic impregnated cement. *J Bone Jt Infect.* 2018;3(3):123 – 129.
12. Unger AC, Dirksen B, Renken FG, Wilde E, Willkomm M, Schulz AP. Treatment of femoral neck fracture with a minimal invasive surgical approach for hemiarthroplasty - clinical and radiological results in 180 geriatric patients. *Open Orthop J.* 2014;8:225 – 31.
13. Sullivan NPT, Hughes AW, Halliday RL, Ward AL, Chesser TJS. Early complications following cemented modular hip hemiarthroplasty. *Open Orthop J.* 2015;9:15 – 9.
14. Kiran M, Johnston LR, Sripada S, Mcleod GG, Jariwala AC. Cemented total hip replacement in patients under 55 years. *Acta Orthop.* 2018;89(2):152 – 5.

Titel: SSCP PALACOS LV+G

Doc.-nr.: 61342

Pagina 14 van
18

15. Sanz-Ruiz P, Matas-Diez JA, Sanchez-Somolinos M, Villanueva-Martinez M, Vaquero-Martín J. Is the commercial Antibiotic-loaded bone cement useful in prophylaxis and cost saving after knee and hip joint arthroplasty? The transatlantic paradox. *J Arthroplasty*. 2017;32(4):1095 – 9.
16. Pallaver A, Zwicky L, Bolliger L, Bösebeck H, Manzoni I, Schädelin S, Ochsner PE, Clauss M. Long-term results of revision total hip arthroplasty with a cemented femoral component. *Arch Orthop Trauma Surg*. 2018;138(11):1609 – 16.
17. Born P, Manzoni I, Ilchmann T, Clauss M. Is cemented revision total hip arthroplasty a reasonable treatment option in an elderly population. *Orthop Rev (Pavia)*. 2019;11(3):8263.
18. Ebied AM, Ebied AA, Marei S, Smith E. Enhancing biology and providing structural support for acetabular reconstruction in single-stage revision for infection. *J Orthop Traumatol*. 2019;20(1):23.
19. Ilchmann T, Zimmerli W, Ochsner PE, Kessler B, Zwicky L, Graber P, Clauss M. One-stage revision of infected hip arthroplasty: outcome of 39 consecutive hips. *Int Orthop*. 2016;40(5):913-8.
20. Brüggemann A, Mallmin H, Hailer NP. Do dual-mobility cups cemented into porous tantalum shells reduce the risk of dislocation after revision surgery. *Acta Orthop*. 2018;89(2):156 – 62.
21. Mueller Greber P, Manzoni I, Ochsner PE, Ilchmann T, Zwicky L, Clauss M. Excellent long-term results of the Müller acetabular reinforcement ring in primary cup revision. *Acta Orthop*. 2017;88(6):619 – 26.
22. Tootsi K, Heesen V, Lohrengel M, Enz AE, Illiger S, Mittelmeier W, Lohmann CH. The use of antibiotic-loaded bone cement does not increase antibiotic resistance after primary total joint arthroplasty. *Knee Surg Sports Traumatol Arthrosc*. 2021.
23. Birkeland Ø, Espehaug B, Havelin LI, Furnes O. Bone cement product and failure in total knee arthroplasty. *Acta Orthop*. 2017;88(1):75 – 81.
24. Bendich I, Zhang N, Barry JJ, Ward DT, Whooley MA, Kuo AC. Antibiotic-laden bone cement use and revision risk after primary total knee arthroplasty in U.S. veterans. *J Bone Joint Surg Am*. 2020;102(22):1939 – 47.
25. Hauer G, Hörlesberger N, Klim S, Bernhardt GA, Leitner L, Glehr M, Leithner A, Sadoghi P. Mid-term results show no significant difference in postoperative clinical outcome, pain and range of motion between a well-established total knee arthroplasty design and its successor: a prospective, randomized, controlled trial. *Knee Surg Sports Traumatol Arthrosc*. 2021;29(3):827-831.
26. Fuchs A, Häussermann P, Hömig D, Ochs BG, Klopfer T, Müller CA, Helwig P, Konstantinidis L. 10-year follow-up of the Columbus knee prostheses system in a prospective multicenter study. *Arch Orthop Trauma Surg*. 2021.
27. Rassir R, Nolte PA, Lugt JCT, Nelissen RGHH, Sierevelt IN, Verra WC. No differences in cost-effectiveness and short-term functional outcomes between cemented and uncemented total knee arthroplasty. *BMC Musculoskelet Disord*. 2020;21(1):448.
28. Sanz-Ruiz P, Villanueva-Martínez M, Matas-Diez JA, Vaquero-Martín J. Revision TKA with a condylar constrained prosthesis using metaphyseal and surface cementation: a minimum 6-year follow-up analysis. *BMC Musculoskelet Disord*. 2015;16:39.
29. Baur J, Zwicky L, Hirschmann MT, Ilchmann T, Clauss M. Metal backed fixed-bearing unicompartmental knee arthroplasties using minimal invasive surgery: a promising outcome analysis of 132 cases. *BMC Musculoskelet Disord*. 2015;16:177.

Titel: SSCP PALACOS LV+G

Doc.-nr.: 61342

Pagina 15 van
18

30. Böhler C, Kolbitsch P, Schuh R, Lass R, Kubista B, Giurea A. Midterm results of a new rotating hinge knee implant: A 5-year follow-up. *Biomed Res Int.* 2017;2017:7532745.
31. Giurea A, Neuhaus H-J, Miehke R, Schuh R, Lass R, Kubista B, Windhager R. Early results of a new rotating hinge knee implant. *Biomed Res Int.* 2014;2014:948520.
32. Staats K, Wannmacher T, Weihs V, Koller U, Kubista B, Windhager R. Modern cemented total knee arthroplasty design shows a higher incidence of radiolucent lines compared to its predecessor. *Knee Surg Sports Traumatol Arthrosc.* 2019;27(4):1148 – 55.
33. Petursson G, Fenstad AM, Gøthesen, Haugan K, Dyrhovden GS, Hallan G, Röhrli SM, Aamodt A, Nilsson KG, Furnes O. Similar migration in computer-assisted and conventional total knee arthroplasty. *Acta Orthop.* 2017;88(2):166 – 72.
34. Shah N, Iqbal HJ, Brookes-Fazakerley S, Sinopidis C. Shoulder hemiarthroplasty for the treatment of three and four part fractures of the proximal humerus using Comprehensive® Fracture stem. *Int Orthop.* 2011;35(6):861 – 7.
35. Ernstbrunner L, Hingsammer A, Imam MA, Sutter R, Brand B, Meyer DC, Wieser K. Long-term results of total elbow arthroplasty in patients with hemophilia. *J Shoulder Elbow Surg.* 2018;27(1):126-32.
36. Wang X, Luo F, Huang K, Xie Z. Induced membrane technique for the treatment of bone defects due to post-traumatic osteomyelitis. *Bone Joint Res.* 2016;5(3):101-5.
37. Marais LC, Ferreira N. Bone transport through an induced membrane in the management of tibial bone defects resulting from chronic osteomyelitis. *Strategies Trauma Limb Reconstr.* 2015;10(1):27 – 33.
38. Scholz AO, Gehrmann S, Glombitza M, Kaufmann RA, Bostelmann R, Flohe S, Windolf J. Reconstruction of septic diaphyseal bone defects with the induced membrane technique. *Injury.* 2015;46 Suppl 4:S121 – 4.
39. Helbig L, Bechberger M, Aldeeri R, Ivanova A, Haubruck P, Miska M, Schmidmaier G, Omlor GW. Initial peri- and postoperative antibiotic treatment of infected nonunions: results from 212 consecutive patients after mean follow-up of 34 months. *Ther Clin Risk Manag.* 2018;14:59 – 67.
40. Omlor GW, Lohnherr V, Lange J, Gantz S, Merle C, Fellenberg J, Raiss P, Lehner B. Enchondromas and atypical cartilaginous tumors at the proximal humerus treated with intralesional resection and bone cement filling with or without osteosynthesis: retrospective analysis of 42 cases with 6 years mean follow-up. *World J Surg Oncol.* 2018;16(1):139.
41. Kasha S, Rathore SS, Kumar H. Antibiotic cement spacer and induced membrane bone grafting in open fractures with bone loss: a case series. *Indian J Orthop.* 2019;53(2):237 – 45.
42. Tetsworth K, Woloszyk A, Glatt V. 3D printed titanium cages combined with the Masquelet technique for the reconstruction of segmental femoral defects: Preliminary clinical results and molecular analysis of the biological activity of human-induced membranes. *OTA Int.* 2019;2(1):016.
43. de Jong L, Klem T, Kuijper TM, Roukema GR. Factors affecting the rate of surgical site infection in patients after hemiarthroplasty of the hip following a fracture of the neck of the femur. *Bone Joint J.* 2017;99-B(8):1088-94.
44. Kienapfel H, Hildebrand R, Neumann T, Specht R, Koller M, Celik I, Mueller HH, Griss P, Klose KJ, Georg C. The effect of Palamed G bone cement on early migration of tibial components in total knee arthroplasty. *Inflamm Res.* 2004;53 Suppl 2:S159-63.

Titel: SSCP PALACOS LV+G

Doc.-nr.: 61342

Pagina 16 van
18

45. Yen SH, Chen JH, Lu YD, Wang JW. Perioperative complications of total knee arthroplasty in dialysis patients. *J Arthroplasty*. 2018;33(3):872-7.
46. Mjoberg B. Loosening of the cemented hip prosthesis. The importance of heat injury. *Acta Orthop Scand Suppl*. 1986;221:1-40.
47. Mjoberg B, Franzen H, Selvik G. Early detection of prosthetic-hip loosening. Comparison of low- and high-viscosity bone cement. *Acta Orthop Scand*. 1990;61(3):273-4.
48. Mjoberg B, Rydholm A, Selvik G, Onnerfalt R. Low- versus high-viscosity bone cement. Fixation of hip prostheses analyzed by roentgen stereophotogrammetry. *Acta Orthop Scand*. 1987;58(2):106-8.
49. Chotel F, Nguiabanda L, Brailon P, Kohler R, Berard J, Abelin-Genevois K. Induced membrane technique for reconstruction after bone tumor resection in children: a preliminary study. *Orthop Traumatol Surg Res*. 2012;98(3):301-8.
50. Berberich C, Sanz-Ruiz P. Risk assessment of antibiotic resistance development by antibiotic-loaded bone cements: is it a clinical concern? *EFORT Open Rev*. 2019;4(10):576-84.
51. Martínez-Moreno J, Merino V, Nácher A, Rodrigo JL, Clemente M, Merino-Sanjuán M. Antibiotic-loaded bone cement as prophylaxis in total joint replacement. *Orthop Surg*. 2017;9(4):331-41.
52. Xu YM, Peng HM, Feng B, Weng XS. Progress of antibiotic-loaded bone cement in joint arthroplasty. *Chin Med J (Engl)*. 2020;133(20):2486-94.
53. Wahlig H, Dingeldein E. Antibiotics and bone cements. Experimental and clinical long-term observations. *Acta Orthop Scand*. 1980;51(1):49-56. Salvati EA, Callaghan JJ, Brause BD, Klein RF, Small RD. Reimplantation in infection. Elution of gentamicin from cement and beads. *Clin Orthop Relat Res*. 1986(207):83-93.
54. Cambon-Binder A, Revol M, Hannouche D. Salvage of an osteocutaneous thermonecrosis secondary to tibial reaming by the induced membrane procedure. *Clin Case Rep*. 2017;5(9):1471 – 1476.
55. Schmitt JW, Benden C, Dora C, Werner CML. Is total hip arthroplasty safely performed in lung transplant patients? Current experience from a retrospective study of the Zurich lung transplant cohort. *Patient Saf Surg*. 2016;10:17.
56. Grubhofer F, Imam MA, Wieser K, Achermann Y, Meyer DC, Gerber C. Staged revision with antibiotic spacers for shoulder prosthetic joint infections yields high infection control. *Clin. Orthop. Relat. Res*. 2018;476(1):146 – 152.
57. Wilairatana V, Sinlapavilawan P, Honsawek S, Limpaphayom N. Alteration of inflammatory cytokine production in primary total knee arthroplasty using antibiotic-loaded bone cement. *J Orthop Traumatol*. 2017;18(1):51 – 7.
58. Agni NR, Costa ML, Achten J, Peckham N, Dutton SJ, Png ME, Reed MR, WHITE 8 Investigators. High-dose dual-antibiotic loaded cement for hip hemiarthroplasty in the UK (WHITE 8): a randomised controlled trial. *Lancet*. 2023;402(10397):196-202.
59. Baryeh K, Wang C, Sochart DH. Periprosthetic femoral fractures around the original cemented polished triple-tapered C-stem femoral implant: a consecutive series of 500 primary total hip arthroplasties with an average follow-up of 15 years. *Arch Orthop Trauma Surg*. 2023;143(7):4511-8.
60. Bourbotte-Salmon F, Ferry T, Cardinale M, Servien E, Rongieras F, Fessy M-H, Bertani A, Laurent F, Buffe-Lidove M, Batailler C, Lustig S, Lyon Bone and Joint Infections Study Group. Rotating hinge knee

Titel: SSCP PALACOS LV+G

Doc.-nr.: 61342

Pagina 17 van
18

- arthroplasty for revision prosthetic-knee infection: good functional outcomes but a crucial need for superinfection prevention. *Front Surg.* 2021;8:551814.
61. Crego-Vita D, Aedo-Martín D, Garcia-Cañas R, Espigares-Correa A, Sánchez-Pérez C, Berberich CE. Periprosthetic joint infections in femoral neck fracture patients treated with hemiarthroplasty - should we use antibiotic-loaded bone cement? *World J Orthop.* 2022;13(2):150-9.
 62. EPRD. Detail Report for Heraeus Medical GmbH PALACOS® LV+G (B1.1 Bone Cement). EPRD Deutsche Endoprothesenregister gGmbH; 2022.
 63. Matar HE, Bloch BV, James PJ. High ten-year implant survivorship and low patellofemoral complication rate for S-ROM rotating-hinge implants in revision total knee arthroplasty: a single-centre study. *Bone Jt Open.* 2022;3(3):205-10.
 64. Mukka S, Hailer NP, Möller M, Gordon M, Lazarinis S, Rogmark C, Östlund O, Sköldenberg O, Wolf O, DAICY study group. Study protocol: The DAICY trial-dual versus single-antibiotic impregnated cement in primary hemiarthroplasty for femoral neck fracture-a register-based cluster-randomized crossover-controlled trial. *Acta Orthop.* 2022;93:794-800.
 65. Ott N, Müller C, Jacobs A, Paul C, Wegmann K, Müller LP, Kabir K. Outcome of geriatric proximal humeral fractures: a comparison between reverse shoulder arthroplasty versus open reduction and internal fixation. *OTA Int.* 2022;5(2 Suppl):188.
 66. Reinbacher P, Wittig U, Hauer G, Draschl A, Leithner A, Sadoghi P. Impact of the COVID-19 pandemic on early clinical outcome after total knee arthroplasty: a retrospective comparative analysis. *Arch Orthop Trauma Surg.* 2023;143(6):3319-26.
 67. National Joint Registry. Implant Summary Report for Elbows using PALACOS LV+G Cement [in primary elbow procedures]. National Joint Registry; 2024.
 68. National Joint Registry. Implant Summary Report for Elbows using PALACOS LV+G Cement [in revision elbow procedures]. National Joint Registry; 2024.
 69. National Joint Registry. Implant Summary Report for Hips using PALACOS LV+G Cement [in primary hip procedures]. National Joint Registry; 2024.
 70. National Joint Registry. Implant Summary Report for Hips using PALACOS LV+G Cement [in revision hip procedures]. National Joint Registry; 2024.
 71. National Joint Registry. Implant Summary Report for Knees using PALACOS LV+G Cement [in primary knee procedures]. National Joint Registry; 2024.
 72. National Joint Registry. Implant Summary Report for Knees using PALACOS LV+G Cement [in revision knee procedures]. National Joint Registry; 2024.
 73. National Joint Registry. Implant Summary Report for Shoulders using PALACOS LV+G Cement [in primary shoulder procedures]. National Joint Registry; 2024.
 74. National Joint Registry. Implant Summary Report for Shoulders using PALACOS LV+G Cement [in revision shoulder procedures]. National Joint Registry; 2024.
 75. EPRD. EPRD Data Summary for PALACOS® LV+G (B1.1 Bone Cement). EPRD Deutsche Endoprothesenregister gGmbH; 2023.
 76. Angelomenos V, Mohaddes M, Kärrholm J, Malchau H, Shareghi B, Itayem R. A prospective randomized study of Refobacin Bone Cement R versus Palacos R + G. *Bone Joint J.* 2024;106-B(5):435-41.

Titel: SSCP PALACOS LV+G

Doc.-nr.: 61342

Pagina 18 van
18

77. Duncan ST, Sabatini F. The use of calcium sulfate/hydroxyapatite bone graft substitute to restore acetabular bone loss in revision total hip arthroplasty. *Arthroplast Today*. 2023;23:101217
78. Giovanoulis V, Koutserimpas C, Lepidas N, Vasiliadis AV, Batailler C, Ferry T, Lustig S. Restoring the anatomy of long bones with large septic non-union defects with the masquelet technique. *Maedica (Bucur)*. 2023;18(3):413-9.
79. Heraeus Medical GmbH, Register Report PALACOS® R+G. 2024
80. Holzer LA, Finsterwald MA, Sobhi S, Jones CW, Yates PJ. Application of bone cement directly to the implant in primary total knee arthroplasty. Short-term radiological and clinical follow-up of two different cementing techniques. *Arch Orthop Trauma Surg*. 2024;144(1):333-340.
81. Howie C, McCrosson M, Padgett AM, Sanchez T, McGwin G, Naranje S. The economic and clinical impact of fast- versus slow-setting cement in primary total knee arthroplasty. *Arch Orthop Trauma Surg*. 2024;144(1):15-21.
82. Matar HE, Bloch BV, James PJ. Satisfactory medium- to long-term outcomes of cemented rotating hinge prosthesis in revision total knee arthroplasty. *Bone Jt Open*. 2023;4(10):776-81.
83. O'Donovan P, McAleese T, Harty J. Does lucency equate to revision? A five-year retrospective review of Attune and Triathlon total knee arthroplasty. *Knee Surg Sports Traumatol Arthrosc*. 2023;31(11):4773-4781.
84. Png ME, Costa M, Nickil A, Achten J, Peckham N, Reed MR, WHITE-8 Investigators. Cost-utility analysis of dual-antibiotic cement versus single-antibiotic cement for the treatment of displaced intracapsular hip fractures in older adults. *Bone Joint J*. 2023;105-B(10):1070-1077.
85. Saccomanno MF, Lädemann A, Collin P. Two-stage exchange arthroplasty for periprosthetic reverse shoulder arthroplasty infection provides comparable functional outcomes to primary reverse shoulder arthroplasty. *J Clin Med*. 2024;13(3).
86. Shi J, Yang X, Song M, Zhang X, Xu Y. Clinical effects of early debridement, internal fixation, and Masquelet technique for childhood chronic haematogenous osteomyelitis of long bones. *J Orthop Surg Res*. 2023;18(1):11
87. Turnbull G, Blacklock C, Akhtar A, Dunstan E, Ballantyne JA. Experience of an anatomic femoral stem in a UK orthopaedic centre beyond 20 years of follow-up. *Eur J Orthop Surg Traumatol*. 2024;34(4):2155-62.
88. Turnbull GS, Akhtar MA, Dunstan ERR, Ballantyne JA. Experience of an anatomic femoral stem in a United Kingdom center - excellent survivorship and negligible periprosthetic fracture rates at mean 12 years following primary total hip arthroplasty. *J Arthroplasty*. 2024;39(1):187-192.
89. Wagh A, Tandel J, Ballyapally D, Jagtap K, Bharadwaj B. Can "aseptic" looking tibia non-union be result of an unrecognized subclinical infection. *J Orthop Case Rep*. 2023;13(12):75-9.
90. Yang X, Xu X, Li J, Song M, Sun H, Zhang H, Zhang X, Xu Y, Shi J. Management of infected bone defects of the femoral shaft by Masquelet technique: sequential internal fixation and nail with plate augmentation. *BMC Musculoskelet Disord*. 2024;25(1):552.