


Heraeus

Medical

# OSTEOPAL<sup>®</sup> V

Low-viscosity, radiopaque bone cement

A stylized, dark gray silhouette of a pterosaur is shown in flight, positioned horizontally across the middle of the page. The pterosaur's wings are spread wide, and its long, pointed beak is directed towards the right. The background consists of two overlapping, light gray curved shapes that create a sense of depth and movement. The entire illustration is set against a white background.

Instructions for use

**Rx** Only

# OSTEOPAL® V

|           |                               |          |
|-----------|-------------------------------|----------|
| <b>EN</b> | <b>Instructions for use</b>   | <b>3</b> |
|           | <b>Working times</b>          | <b>9</b> |
|           | (Revision status: 2024-10-25) |          |

## 1 General information

### 1.1 Device description

**OSTEOPAL®V** is a standard-setting, low-viscosity, radiopaque, poly(methyl methacrylate)-based (PMMA) bone cement. It contains the X-ray contrast medium zirconium dioxide. To improve visibility in the surgical field, it has been colored with chlorophyll-copper-complex (E141). The bone cement consists of two components and is prepared immediately before use by mixing the polymer powder (= powder) with the monomer liquid (= liquid). The low viscosity increases with progressing polymerization. A paste forms that is inserted into the vertebral body using an application system and sets within a few minutes.

**OSTEOPAL®V** is intended for single use and is supplied sterile.

### 1.2 Device sizes and packaging

#### Device sizes

**OSTEOPAL®V** is available in the size 20.

#### Package content

- 1 pouch filled with powder
- 1 brown glass ampule filled with liquid
- 1 set of product stickers for documentation
- 1 instructions for use

#### Packaging design

The powder is triple packaged: it is packed first in an inner polyethylene (PE) paper pouch and secondly in an outer PE-paper pouch (peel-off). Both PE-pouches are sterilized using ethylene oxide. The peel-off pouch contains the sterile inner PE-paper pouch with powder and is sterile on the inside and non-sterile on the outside. Afterwards the sterilized pouches are packed in a protective aluminum pouch.

The liquid is single packaged. The ampule blister is sterilized using ethylene oxide and contains a brown glass ampule with sterile-filtered liquid.

### 1.3 Composition

| Composition               |      |
|---------------------------|------|
| <b>Powder:</b>            |      |
| PMMA copolymer            | 54 % |
| zirconium dioxide         | 45 % |
| benzoyl peroxide          | 1 %  |
| <b>Liquid:</b>            |      |
| methyl methacrylate       | 98 % |
| N, N-dimethyl-p-toluidine | 2 %  |

The data is rounded

#### Other constituents:

- Powder: chlorophyll-copper-complex (E141)
- Liquid: chlorophyll-copper-complex (E141), hydroquinone

**OSTEOPAL®V** does not contain a radiation source.

| Device size    | 20     |
|----------------|--------|
| <b>Powder:</b> | 26.3 g |
| <b>Liquid:</b> | 10 ml  |

The data is rounded

The mass ratio is about 27 % liquid to 73 % powder.

### 1.4 Supporting equipment

Mixing equipment is required to mix the two components of **OSTEOPAL®V**. The mixing bowl with spatula distributed by Heraeus Medical GmbH is suitable.

**OSTEOPAL®V** can be introduced into the vertebral body using an application system that has been approved for vertebroplasty or kyphoplasty and enables steady and controlled injection. The user must ensure that the application system is compatible with the ingredients of the **OSTEOPAL®V** composition. For the correct usage and surgical technique, refer to the instructions of the manufacturer. The handling of the system is described in the manufacturer's instructions for use.

Cannulae with an inner diameter of less than 1.8mm (13G) should not be used.

The instructions for use of the supporting equipment must be followed.

**Note:** Heraeus Medical GmbH has not tested the compatibility of **OSTEOPAL®V** with devices of other manufacturers and does not assume any liability for this. The use of mixing equipment of other manufacturers is done in the sole discretion and responsibility of the user.

## 1.5 MRI safety information

**OSTEOPAL®V** is MR safe.

## 2 Directions for use

### 2.1 Intended use

**OSTEOPAL®V** bone cement is intended for the treatment of pathological fractures of the vertebral body due to osteoporosis, cancer, or benign lesions using a vertebroplasty or balloon kyphoplasty procedure.

### 2.2 Intended patient population

Adult population.

### 2.3 Intended user

Healthcare professionals in a clinical context and experienced in the handling of the product. The surgeon and nurse must be thoroughly familiar with the properties and handling characteristics of **OSTEOPAL®V**. As the handling of **OSTEOPAL®V** varies with temperature, humidity and mixing technique, a test mix should be performed to ensure familiarity with its characteristics.

### 2.4 Indications

**OSTEOPAL®V** bone cement is indicated for the treatment of pathological fractures of the vertebral body due to osteoporosis, cancer, or benign lesions using a vertebroplasty or balloon kyphoplasty procedure.

### 2.5 Contraindications

**OSTEOPAL®V** must not be used

- in patients with hemorrhagic diathesis and infections
- in patients with lesions of the vertebral body with epidural extension due to the risk of spinal cord compression
- in patients with suspected or known hypersensitivity to constituents of the bone cement
- spinal stenosis (> 20% by retropulsed fragments)
- compromised of the vertebral fractures due to posterior involvement
- patient clearly improving on medical therapy
- prophylaxis in metastatic or osteoporotic patients with no evidence of acute fracture

- non-pathological, acute traumatic fractures of the vertebra
- vertebral plana (collapse > 90 %)
- compromise of the vertebral body or the walls of the pedicles
- unstable vertebral fractures due to posterior involvement

PMMA bone cement is contraindicated in the presence of active or incompletely treated infection, at the site where the bone cement is to be applied.

## 2.6 Adverse events and residual risks

### Adverse Events

- Pneumonia
- Pneumothorax
- Fracture of a pedicle
- Rib fracture in patients with diffuse osteopenia, especially during thoracic vertebroplasty procedures, due to the significant downward force exerted during needle insertion

Serious adverse events, some with fatal outcome, associated with the use of acrylic bone cements for vertebroplasty or kyphoplasty include myocardial infarction, cardiac arrest, cerebrovascular accident, pulmonary embolism, and cardiac embolism. Although the majority of these adverse events present early within the post-operative period, there have been some reports of diagnoses beyond a year or more after the procedure.

Other reported adverse events for acrylic bone cements intended for vertebroplasty or kyphoplasty include: Leakage of the bone cement beyond the site of its intended application with introduction into the vascular system resulting in embolism of the lung and/or heart or other clinical sequelae.

The most frequent adverse reactions reported with acrylic bone cements are a transitory fall in blood pressure, thrombophlebitis, hemorrhage and hematoma, superficial or deep wound infection and short-term cardiac conduction irregularities. Other reported adverse reactions include heterotopic new bone formation.

Other reported adverse events for acrylic bone cements include pyrexia due to an allergy to the bone cement, hematuria, dysuria, bladder fistula and delayed sciatic nerve entrapment due to extrusion of the bone cement beyond the region of its intended application, and adhesions and stricture of the ileum due to the heat released during polymerization.

Cement leakage between intervertebral discs.

Furthermore, application of **OSTEOPAL®V** can result in temporary neurological deficits.

### Residual risk

Residual risks listed below are procedure related risks which are beyond the control of the manufacturer because they are procedure or user related.

#### Vascular System, Heart, Respiratory System, Blood and Lymphatic System

##### Bone Cement Implantation Syndrome (BCIS):

Insertion of bone cement may produce a high intravertebral pressure that forces bone marrow constituents into the venous vascular system, resulting in fat and marrow emboli.

In case of pulmonary or cardiovascular events it is necessary to monitor blood volume and possibly increase it.

In the case of acute respiratory failure, anesthesiologic measures should be taken. In general, adverse reactions of BCIS might include low blood pressure/hypotension, hypoxia, bradycardia, tachycardia, pulmonary hypertension, thrombosis, embolism, pulmonary embolism, myocardial infarction, cerebrovascular accident, respiratory arrest, and cardiac arrest.

- BCIS grade 1  
moderate hypoxia ( $\text{SpO}_2 < 94\%$ ) or hypotension  
[fall in systolic blood pressure  $> 20\%$ ]
- BCIS grade 2  
severe hypoxia ( $\text{SpO}_2 < 88\%$ ) or hypotension  
[fall in systolic blood pressure  $> 40\%$ ] or unexpected loss of consciousness.
- BCIS grade 3  
cardiovascular collapse, requiring CPR
- Asymptomatic cement leakage
- Cement leakage, which can cause myocardial infarction, cardiac arrest, and cerebrovascular accident
- Symptomatic pulmonary embolism
- Symptomatic cardiac embolism

#### Nervous System

- Permanent neurological deficit

#### Musculoskeletal System

- Loss of range of motion
- Ambulation difficulties
- Re-collapse of treated or adjacent vertebrae/secondary fractures

#### Infection

- Bacterial infection including cellulitis, and/or osteomyelitis

#### Generalized Disorders

- Inflammation
- Swelling/edema
- Fibrosis
- Pain

### 2.7 Warnings

#### Regarding intended users

Caution should be exercised during the mixing of the two components of **OSTEOPAL®V** to prevent excessive exposure to the concentrated monomer vapors, which may produce irritation of the respiratory tract, eyes, and possibly the liver. Personnel wearing contact lenses should not be near or involved in mixing this bone cement. Manufacturers of soft contact lenses recommend removing the lenses in the presence of damaging or irritant vapors. Since soft contact lenses are permeable to liquids and gases, they should not be worn in the operating room if methyl methacrylate is being used.

Eye protection is recommended when opening the ampoule during the sterile preparation steps of the bone cement to protect the eye from any glass fragments or monomer liquid.

The monomer is a powerful lipid solvent and should not come into direct contact with the body. When handling **OSTEOPAL®V** it is essential to wear gloves that provide the necessary protection against penetration of the monomer into the skin. Three-layered PVP gloves (polyethylene, ethylene vinyl alcohol copolymer, and polyethylene) and Viton®/butyl gloves have proved to provide good protection over an extended period. It is recommended that two pairs of gloves be worn over one another, e.g., a polyethylene surgical glove over an inner pair of standard latex surgical gloves. Do not allow the monomer to contact latex or polystyrene-butadiene gloves. Request confirmation from your glove supplier that the respective gloves are suitable for use with this bone cement.

The monomer is not for injection.

Polymerization of the bone cement is an exothermic reaction, which occurs while the bone cement is hardening in situ. The released heat may damage bone or other tissues.

Avoid over-pressurizing the bone cement because this may lead to extrusion of the bone cement beyond the site of its intended application and damage to the surrounding tissue. Cement leakage may cause tissue damage, nerve or circulatory problems, and other serious adverse events. If bone cement is seen outside the vertebral body or in the circulatory system during the procedure, immediately stop the injection.

Unanticipated postoperative events may affect the cement-bone interface and lead to micro motion of bone cement against bone surface. A fibrous tissue layer may develop between the bone cement and the bone. Long-term follow-up is advised for all patients on a regularly scheduled basis.

**Note:** **OSTEOPAL®V** is a single-use device and must never be re-used! Re-use may result in diminished safety, performance, and compliance with relevant specifications.

#### Regarding the intended patient population

There is very little evidence with pregnant and lactating women, children and adolescents. Therefore, it is not recommended to use **OSTEOPAL®V**. If no other option is available, the decision whether to use **OSTEOPAL®V** lies with the attending physician. The benefits for the mother should be weighed against the potential risk to the child before using **OSTEOPAL®V** during pregnancy and lactation.

Monitor patients carefully for any change in blood pressure during and immediately after the application of bone cement. Adverse patient reactions involving the cardiovascular system are in particular linked to the pressurization of bone cement. Hypotensive reactions have occurred shortly after application of bone cement. However, consequences such as cardiac arrest are only reported in very few cases.

Few publications reported gradual deterioration of baseline cardiovascular parameters during consecutive augmentation of several vertebral bodies with PMMA. To avoid this, not more than 3 augmentations in a single operation are recommended.

In addition, paravertebral structures may be damaged by cement escaping. Complications such as spinal cord compression, intercostal neuralgia, cement escaping into the intervertebral space, peri-vertebral filling of veins and arteries (danger of embolism), infections and post-procedural pain are possible. To prevent cement escape and in order to detect unwanted occurrences in good time, application must be performed using imaging techniques (real-time display). In order to detect a vertebral extension, a phlebography must be performed before applying the cement. Immediate surgical intervention may also be necessary to counteract the described complications. Prior to surgery, a careful radiological investigation must be performed to assess possible risks (e.g. existing vertebral-body lesions, vascular supply of the vertebral body or edema).

## **2.8 Precautions**

#### Regarding intended users

Do not use **OSTEOPAL®V** after the expiration date printed on the folding box. This device may not be safe or effective beyond its expiration date.

Follow the handling and mixing instructions to avoid contact dermatitis. Strict adherence to the instructions for mixing the powder and liquid components may reduce the incidence of this complication.

Adequately ventilate the operating room to eliminate as much monomer vapor as possible.

The liquid is highly volatile and flammable. Ignition of monomer fumes caused by use of electrocautery devices in surgical sites near freshly implanted bone cements has been reported.

Do not use the bone cement after the application phase.

Do not use the bone cement if its consistency is inhomogeneous.

#### Regarding the intended patient population

Blood pressure, pulse, and breathing must be monitored carefully during and immediately after introduction of the bone cement. Any significant change in these vital signs must be resolved without delay by taking appropriate action.

Incomplete filling of the vertebral body with bone cement may result in an insufficient reduction in acute pain and reduced long-term stability of the treated vertebral body.

## **3 Handling**

Perform mixing and application of **OSTEOPAL®V** under sterile conditions. **OSTEOPAL®V** is mixed in a mixing bowl (see 1.4 Supporting equipment).

Prior to using **OSTEOPAL®V**, consult the instructions for use of the selected supporting equipment. Make sure that the handling steps of the supporting equipment are followed adequately (see 1.4 Supporting equipment).

If **OSTEOPAL®V** is used pre-chilled, it is recommended that the bone cement components be pre-chilled at 4 °C – 7 °C for at least 24 hours. It should be removed from cooling and filled into the mixing bowl just before mixing.

The mixing, waiting, application, and setting times of **OSTEOPAL®V** are shown in the diagrams at the end of the instructions for use (see 7.1 Working times). The values are valid for application using a 3 ml single-use syringe and a cannula with a diameter of 3.0 mm and a length of 10 cm. The processing may be different for other application systems, syringe diameter, cannula length and cannula diameter.

**Note:** These are stated for guidance only because the application and setting times depend not only on temperature but also on mixing method used and the humidity in the operating room. The storage temperature also influences the application time. In general, higher temperatures during storage and in the operating room, higher humidity, and vigorous mixing of the bone cement lead to shorter application times and vice versa.

Apart from the handling procedure described below, **OSTEOPAL®V** must not be modified in any way. Any admixing of substances, especially aqueous (e.g., antibiotic containing) solutions, has a considerable detrimental effect on the mechanical properties of the bone cement.

**OSTEOPAL®V** contains the following non-sterile/sterile packaging components:

| Packaging component  | Condition                               |
|----------------------|-----------------------------------------|
| <b>Powder:</b>       |                                         |
| aluminum pouch       | non-sterile                             |
| outer PE-paper pouch | outside: non-sterile<br>inside: sterile |
| inner PE-paper pouch | sterile                                 |
| <b>Liquid:</b>       |                                         |
| ampoule blister      | outside: non-sterile<br>inside: sterile |
| ampoule(s)           | sterile                                 |

### 3.1 Amount required

The amount of bone cement required depends on the specific surgical intervention, on the particular anatomical conditions, on the technique being used as well as the size of the particular vertebral body and the defect.

To avoid gradual deterioration of baseline cardiovascular parameters, not more than 3 augmentations in a single operation are recommended.

At least one additional pack of **OSTEOPAL®V** should be available before commencing the operation.

### 3.2 Non-sterile preparation steps

#### Non-sterile user

- Open the folding box, take out the aluminum pouch, the ampoule blister, the instructions for use, and the set of product stickers.
- Before opening the aluminum pouch, move the content down by shaking or tapping to ensure that when the aluminum pouch is torn open at the top, the content is not damaged.
- Open the aluminum pouch and take out the outer PE-paper pouch.
- Visually inspect the outer PE-paper pouch immediately prior to use in order to determine if breaches in sterile barrier system integrity are evident. Do not use if outer PE-paper pouch is damaged.
- Break the orange sterility seal, then open the outer PE-paper pouch as shown in Figure 1, and present the sterile inner PE-paper pouch to the sterile user for sterile removal.
- Open the ampoule blister using its flap and present the sterile ampoule to the sterile user for sterile removal.

#### Open the outer PE-paper pouch




Figure 1

The opening flaps of the pouch top help to detach the foil from the paper. To maximize the area of the opening flaps that can be grasped, the side of the paper/foil should be kept between thumb, index, and middle finger. Use the whole thumb surface to grasp the foil and paper side and peel off each side evenly.

### 3.3 Sterile preparation steps

It is advisable to first add the liquid and then the powder. If this order is reversed, nests of unpolymerized powder are more likely to form due to polymerization commencing immediately at the surface.

#### Sterile user

- Remove the ampoule and open it under sterile conditions: the liquid is ready to fill into the mixing system.
- Remove the inner PE-paper pouch.
- Before opening the inner PE-paper pouch, move the content down by shaking or tapping to ensure that when the inner PE-paper pouch is cut open at the top, no powder is lost.
- Open the inner PE-paper pouch under sterile conditions: the powder is ready to fill into the mixing system.

**Note:** Do not open the ampoule over the mixing device to prevent contamination of the bone cement with glass fragments. To make it easier to open the ampoule, it is provided with a predetermined breaking point at the transition to the head of the ampoule.

The ampoules are provided with a snap-off device (tube) to facilitate the opening procedure. To do so, grasp the fitted snap-off device instead of the ampoule head and break off the ampoule head using the device. When the ampoule head has snapped off, it remains inside the tube.

### 3.4 Mixing and application

For mixing of **OSTEOPAL®V**, the instructions for use of the suitable system must be followed (see 1.4 Supporting equipment).

Mixing the powder and liquid together produces a paste that is used to fill and stabilize vertebral bodies.

Powder and liquid should only be added to the mixing system just before mixing. During the mixing time of 30 seconds, the two components are mixed with one another while stirring evenly. Always mix the entire contents of a pouch with the entire contents of an ampoule of liquid.

If powder or liquid is spilled, a new package must be used.

### Processing

The viscosity increases with progressing polymerization. The bone cement paste should be placed in an application system immediately after mixing because its viscosity is still low and it is easy to aspirate at this point. To prevent vascular leakage of bone cement, the bone cement should be applied when it is a paste-like consistency.

During intravertebral application, systematic real-time fluoroscopy (latero-lateral) is necessary. In the event of paravertebral cement leakage, cement injection must be immediately stopped and can be continued after the cement's viscosity has been increased. If the vertebral filling is not sufficient, additional contralateral access is viable. After augmentation, a mandrel should be inserted into the hypodermic needle so that no traces of bone cement remain in the soft tissue after removal of the hypodermic needle.

The patient must remain immobilized until the bone cement has fully set.

## **4 Storage, expected lifetime, transport, shelf life, sterilization**

### Storage

Do not store above 25 °C (77 °F).

**OSTEOPAL®V** must be stored in dry conditions and must not be exposed to direct sunlight, ionizing radiation, extremes of temperature, or particulate contamination.

### Expected lifetime

There is no general factor influencing the expected lifetime of **OSTEOPAL®V**. The actual lifetime of **OSTEOPAL®V** can be influenced by factors such as the medical situation and lifestyle of the patient.

Removal of the implant is not necessary unless it is clinically required. Implant removal is in the sole discretion of the responsible surgeon based on the medical condition of the patient.

### Transport

Care shall be exercised during transport and handling of the **OSTEOPAL®V** to avoid any damage or alteration to the performance characteristics of the **OSTEOPAL®V** and its packaging as received. Do not remove the **OSTEOPAL®V** from the sterile packaging until immediately before use. Do not use if packaging is damaged.

### Shelf life

The shelf life of **OSTEOPAL®V** is printed on the folding box. Do not use **OSTEOPAL®V** if the date indicated has expired.

### Sterilization

Powder and packaging have been sterilized using ethylene oxide and the liquid by sterile filtration. The product must not be re-sterilized. Non-sterility may cause an infection in the patient. If the powder has turned yellow, do not use **OSTEOPAL®V**.

## **5 Disposal**

Dispose of the polymer component in an authorized waste facility. The liquid component should be evaporated under a well-ventilated hood or absorbed by an inert material and transferred in a suitable container for disposal.

## **6 Safety and performance information**


In the event of adverse events or malfunctions, contact the manufacturer via [hm.vigilance.medical@heraeus.com](mailto:hm.vigilance.medical@heraeus.com).

Follow the national reporting regulations in your country and also report serious events to the responsible authorities, if required.

## 7 Appendix

### 7.1 Working times

#### Working times for mixing without vacuum (not pre-chilled bone cement)



Test conditions: Not pre-chilled mixing bowl, 55% humidity

**Note:** The working times are stated for guidance only because the application and setting times depend not only on temperature but also on mixing method used and the humidity in the operating room.



# Symbols

|                                                                                     |                                                                   |                                                                                     |                                                                                                    |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------|-------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------|
|    | Manufacturer                                                      |    | Keep dry                                                                                           |
|    | Country and date of manufacture (DE=Germany)                      |    | Upper limit of temperature<br>25°C (77° F)                                                         |
|    | Use-by date                                                       |    | Do not re-use                                                                                      |
|    | Batch code                                                        |    | Consult instructions for use                                                                       |
|    | Catalogue number                                                  |    | Caution                                                                                            |
|     | Sterilized using aseptic processing techniques                    |    | Medical Device                                                                                     |
|    | Do not re-sterilize                                               |    | Unique Device Identifier                                                                           |
|    | Do not use if package is damaged and consult instructions for use |    | Causes skin irritation                                                                             |
|     | Sterilized using ethylene oxide                                   |    | Flammable liquid                                                                                   |
|  | Single sterile barrier system                                     |  | Federal law restricts this device to sale, distribution, and use by or on the order of a physician |
|  | Single sterile barrier system with protective packaging outside   |  | MR safe                                                                                            |
|  | Keep away from sunlight                                           |                                                                                     |                                                                                                    |



Heraeus Medical GmbH  
Philipp-Reis-Straße 8/13  
61273 Wehrheim, Germany  
Phone: +49 (0) 6181 35 33 99  
[www.heraeus-medical.com](http://www.heraeus-medical.com)

OSTEOPAL® is a trademark of Heraeus Medical GmbH.  
Made in Germany.