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## Literatur

1. Robertsson O et al. The Swedish Knee Arthroplasty Register 1975–1997: an update with special emphasis on 41,223 knees operated on in 1988–1997. *Acta Orthop Scand.* 2001; 72(5): 503–13.
2. Jansen E et al. Risk factors for infection after knee arthroplasty. A register-based analysis of 43,149 cases. *J Bone Joint Surg Am.* 2009; 91(1): 38–47.
3. Gandhi R et al. Survival and clinical function of cemented and uncemented prostheses in total knee replacement: a meta-analysis. *J Bone Joint Surg Br.* 2009; 91(7): 889–95.
4. Cawley DT et al. Cementing techniques for the tibial component in primary total knee replacement. *Bone Joint J.* 2013; 95-B(3): 295–300.
5. Bourne RB. Prophylactic use of antibiotic bone cement: an emerging standard in the affirmative. *J Arthroplasty.* 2004; 19(4 Suppl 1): 69–72.
6. Australian Orthopaedic Association. National Joint Registry. Annual Report. Adelaide. 2009.
7. Namba RS et al. Outcomes of routine use of antibiotic-loaded cement in primary total knee arthroplasty. *J Arthroplasty.* 2009; 24(6 Suppl): 44–7.
8. Waddell J et al. Orthopaedic practice in total hip arthroplasty and total knee arthroplasty: results from the Global Orthopaedic Registry (GLORY). *Am J Orthop (Belle Mead NJ).* 2010; 39(9 Suppl): 5–13.
9. Hernandez-Castanos DM, Ponce VV, Gil F. Release of ischaemia prior to wound closure in total knee arthroplasty: a better method? *Int Orthop.* 2008; 32(5): 635–8.
10. Horlocker TT et al. Anesthetic, patient, and surgical risk factors for neurologic complications after prolonged total tourniquet time during total knee arthroplasty. *Anesth Analg.* 2006; 102(3): 950–5.
11. Tai TW et al. Tourniquet use in total knee arthroplasty: a meta-analysis. *Knee Surg Sports Traumatol Arthrosc.* 2010; DOI: 10.1007/s00167-010-1342-7.
12. Miskovsky C, Whiteside LA, White SE. The cemented unicompartmental knee arthroplasty. An in vitro comparison of three cement techniques. *Clin Orthop Relat Res.* 1992(284): 215–20.
13. Amirfeyz R, Bannister G. The effect of bone porosity on the shear strength of the bone-cement interface. *Int Orthop.* 2009; 33(3): 843–6.
14. Reading AD, Miles AW, Bannister GC. Keyholes in acetabular preparation. *Hip Int* 2002; 12: 163–5.
15. Röhrig H. Operation der Kniegelenksendoprothese In: Wirtz DC, ed. *AE-Manual der Endoprothetik Knie.* Berlin: Springer, 2011; p. 127–41.
16. Breusch SJ et al. Lavage technique in total hip arthroplasty: jet lavage produces better cement penetration than syringe lavage in the proximal femur. *J Arthroplasty.* 2000; 15(7): 921–7.
17. Breusch SJ et al. [Significance of jet lavage for in vitro and in vivo cement penetration]. *Z Orthop Ihre Grenzgeb.* 2001; 139(1): 52–63.
18. Breusch SJ, Malchau H. *Bone preparation: femur. The well-cemented total hip arthroplasty.* Berlin: Springer, 2005; p. 125–40.
19. Seeger JB et al. The effect of bone lavage on femoral cement penetration and interface temperature during Oxford unicompartmental knee arthroplasty with cement. *J Bone Joint Surg Am.* 2013; 95(1): 48–53.
20. Schlegel UJ et al. Pulsed lavage improves fixation strength of cemented tibial components. *Int Orthop.* 2010; (im Druck).
21. Clarius M et al. Pulsed lavage reduces the incidence of radiolucent lines under the tibial tray of Oxford unicompartmental knee arthroplasty: pulsed lavage versus syringe lavage. *Int Orthop.* 2009; 33(6): 1585–90.
22. Breusch SJ et al. Cemented hip prosthesis implantation decreasing the rate of fat embolism with pulsed pressure lavage. *Orthopäde.* 2000; 29(6): 578–86.
23. Norton MR, Eyres KS. Irrigation and suction technique to ensure reliable cement penetration for total knee arthroplasty. *J Arthroplasty.* 2000; 15(4): 468–74.
24. Phillips AM, Goddard NJ, Tomlinson JE. Current techniques in total knee replacement: results of a national survey. *Ann R Coll Surg Engl.* 1996; 78(6): 515–20.
25. British Orthopaedic Association and British Association for Surgery of the Knee. *Knee Replacement: A Guide to Good Practice.* London: British Orthopaedic Association, 2011.
26. Kuehn KD, Ege W, Gopp U. Acrylic bone cements: mechanical and physical properties. *Orthop Clin North Am.* 2005; 36(1): 29–39, v-vi.
27. Smeds S, Goertzen D, Ivarsson I. Influence of temperature and vacuum mixing on bone cement properties. *Clin Orthop Relat Res.* 1997(334): 326–34.
28. Dunbar MJ. Antibiotic bone cements: their use in routine primary total joint arthroplasty is justified. *Orthopedics.* 2009; 32(9).
29. Ensing GT et al. Copal bone cement is more effective in preventing biofilm formation than Palacos R-G. *Clin Orthop Relat Res.* 2008; 466(6): 1492–8.
30. Neut D et al. The effect of mixing on gentamicin release from polymethylmethacrylate bone cements. *Acta Orthop Scand.* 2003; 74(6): 670–6.
31. Lewis G. Properties of antibiotic-loaded acrylic bone cements for use in cemented arthroplasties: a state-of-the-art review. *J Biomed Mater Res B Appl Biomater.* 2009; 89B(2): 558–74.
32. Frommelt L, Kühn KD. Antibiotic-loaded cement. In: Breusch SJ, Malchau M, eds. *The well-cemented total hip arthroplasty.* Heidelberg: Springer, 2005; p. 86–92.
33. Jiranek W. Antibiotic-loaded cement in total hip replacement: current indications, efficacy, and complications. *Orthopedics.* 2005; 28(8 Suppl): s873–7.
34. Dunne NJ, Orr JF. Influence of mixing techniques on the physical properties of acrylic bone cement. *Biomaterials.* 2001; 22(13): 1819–26.
35. Wang JS et al. Porosity of bone cement reduced by mixing and collecting under vacuum. *Acta Orthop Scand.* 1993; 64(2): 143–6.
36. Wixson RL, Lautenschlager EP, Novak MA. Vacuum mixing of acrylic bone cement. *J Arthroplasty.* 1987; 2(2): 141–9.
37. Geiger MH et al. The clinical significance of vacuum mixing bone cement. *Clin Orthop Relat Res.* 2001(382): 258–66.
38. Lidgren L, Bodelind B, Moller J. Bone cement improved by vacuum mixing and chilling. *Acta Orthop Scand.* 1987; 58(1): 27–32.
39. Eitenmuller J, Wolbert R, Eisen E. [The effect of circulation on the polymerizing temperature of Palacos. An experimental study (author's transl)]. *Arch Orthop Trauma Surg.* 1981; 98(1): 61–7.
40. Toksvig-Larsen S, Franzen H, Ryd L. Cement interface temperature in hip arthroplasty. *Acta Orthop Scand.* 1991; 62(2): 102–5.
41. Biehl G, Harms J, Hanser U. [Experimental studies on heat development in bone during polymerization of bone cement. Intraoperative measurement of temperature in normal blood circulation and in bloodlessness]. *Arch Orthop Unfallchir.* 1974; 78(1): 62–9.

42. Scuderi GR, Clarke H. Optimizing cementing technique. In: Bellemans J, Ries MD, Victor JMK, eds. Total knee arthroplasty. Berlin Heidelberg, : Springer, 2005; p. 223–7.
43. Vaninbroux M et al. Cementing the femoral component in total knee arthroplasty: which technique is the best? *Knee*. 2009; 16(4): 265–8.
44. Walker PS et al. Control of cement penetration in total knee arthroplasty. *Clin Orthop Relat Res*. 1984(185): 155–64.
45. Kopec M et al. Effect of hand packing versus cement gun pressurization on cement mantle in total knee arthroplasty. *Can J Surg*. 2009; 52(6): 490–4.
46. Vanlommel J et al. Cementing the tibial component in total knee arthroplasty: which technique is the best? *J Arthroplasty*. 2011; 26(3): 492–6.
47. Bert JM, McShane M. Is it necessary to cement the tibial stem in cemented total knee arthroplasty? *Clin Orthop Relat Res*. 1998(356): 73–8.
48. Skwara A et al. Primary stability of tibial components in TKA: in vitro comparison of two cementing techniques. *Knee Surg Sports Traumatol Arthrosc*. 2009; 17(10): 1199–205.
49. Uhlenbrock, A et. al. Influence of time in-situ and implant type on fixation strength of cemented tibial trays – A post mortem retrieval analysis. *Clinical Biomechanics*. 2012; (27): 929–935