
Operator and Service Manual

palaJet® AC Power Adapter


AC Power Supply for Disposable Electronically Pulsed Lavage Suction Apparatus/ Disposable Pulse Lavage System (AC Power Adapter)

palaJet® AC Power Adapter

EN	Operator and Service Manual	3
CS	Příručka pro obsluhu a servis	9
DE	Bedienungs- und Wartungsanleitung	16
EL	Εγχειρίδιο χειριστή και σέρβις	23
ES	Manual del operador y de servicio	30
FI	Käyttö- ja huolto-opas	37
FR	Manuel d'utilisation et d'entretien	44
IT	Manuale operativo e di assistenza	51
NL	Bedienings- en onderhoudshandleiding	58
NO	Bruks- og vedlikeholdsanvisning	65
PL	Instrukcja obsługi i serwisowania	72
SK	Návod na obsluhu a servis	79
SV	Användar- och servicehandbok	86

palaJet® AC Power Adapter

1. Device

Name: AC Power Supply for Disposable Electronically Pulsed Lavage Suction Apparatus/Disposable Pulse Lavage System (AC Power Adapter)

W-205

REF 5203576 (European Plug)
5206240 (British Plug)

2. Specifications

Mains Line Voltage Range:	100 – 240 V ~ (AC), 50/60 Hz, 120 VA, Auto switching
Line Current:	0.50 A RMS @ 120 ~ (AC) (Max.)
Output Power:	48 W
Output Voltage:	DC 7.5 – 12 V
Size:	Height x Wide x Depth: About 21 cm x 9 cm x 6 cm
Net weight:	About 0.8 kg
Valid period (Shelf life):	Three years

3. IEC 60601-1 Classification

Type of protection against electric shock:	Class I
Classification according to the degree of protection against ingress of water:	IPX 0
Mode of operation:	Continuous operation

CAUTION: THIS DEVICE IS NOT SUITABLE FOR USE IN THE PRESENCE OF FLAMMABLE ANESTHETIC OR GASES.

4. Emissions and Immunity

The **AC Power Supply** complies with EMC criteria set forth in IEC 60601-1-2 and EN 60601-1-2.

5. Notice

HIGH VOLTAGE WILL BE PRESENT ON THE POWER CONTROL BOARD.

THERE ARE NO SERVICEABLE PARTS WITHIN THE POWER SUPPLY ENCLOSURE.

Should service be required please contact your local Clean Medical representative for service help and/or warranty information.

6. Intended Use

The **AC Power Supply** is intended to be used for supplying power to the AC powered palaJet® AC Disposable Pulse Lavage System (Disposable Electronically Pulsed Lavage Suction Apparatus) (#: W-205).

7. Warning

- 1) The AC Power Supply should be operated only by qualified professionals.
- 2) Please refer to the user manual of AC powered palaJet® AC Disposable Pulse Lavage System (Disposable Electronically Pulsed Lavage Suction Apparatus) (#: W-205).
- 3) This device (AC Power Adapter) should not be used after the expiry date shown on the package.


8. Initial Inspection

Upon receipt, unpack the **AC Power Supply** and inspect the unit for any obvious damage that may have occurred during shipment. We recommend that this inspection be performed by a professional familiar with electronic medical devices. If the unit is damaged, notify the carrier and your representative immediately.

9. Equipotential


- 1) The equipotential cable is not included in this device, but prepared by medical institution.
The equipotential cable specification: diameter 2mm²
- 2) Connect the socket of the equipotential cable to the equipotential column located at the lower right-hand corner on the back of the AC Power Supply by plugging them together (notice there should be resistance to ensure tight connection).
- 3) Connect the other end of the equipotential cable to the equipotential busbar in the operating room or any other common equipotential connection.

One example of connectivity in the operating room is as follows:

10. Functional Check

If the initial inspection results are satisfactory, connect the power plug of the unit to a properly polarized and grounded power source with voltage and frequency characteristics compatible with the specifications listed in Clause-2. Using a digital volt meter set for DC measurement, place the volt meter ground probe on the middle terminal (2: DC GND) of the power supply connector. Place the other volt meter probe on the left terminal lead (1) and observe the voltage is approximately $12V \pm 1V$ DC. Next move the volt meter probe to the right terminal lead (3) and observe the voltage is also approximately $12V \pm 1V$ DC. If both results are within the limits described the unit can be placed into operational service.


Functional check

- 1: Left terminal lead
- 2: DC GND
- 3: Right terminal lead

11. Check the Device Prior Each Use

- 1) To visual check the device should not have any obvious damage.
- 2) Prior each clinical use of the device, please check whether the device is working properly. Connect the Pulsed lavage kit and the AC Power supply (AC Power Adapter) together, and start the trial operation without connecting the Irrigation bag to check whether it can work normally.
Please refer to the user manual of AC powered palaJet® AC Disposable Pulse Lavage System (Disposable Electronically Pulsed Lavage Suction Apparatus) (#: W-205).
- 3) Only the device that meets the requirements can be put into use.
If it does not meet the requirements, replace it with a new unit.

12. Operation

WARNING: To avoid the risk of electric shock, the unit must only be connected to a supply main with protective earth.

- 1) Hang the unit on an IV pole using the unit's built in hanger.
Note: THE UNIT MUST ALWAYS BE HUNG ON AN IV POLE AND SHOULD NEVER BE ALLOWED TO REST ON A TABLE OR FLOOR. THE TABLE AND/OR FLOOR MAY BECOME WET WHICH COULD DAMAGE THE UNIT.
- 2) Connect the power plug of the unit to a properly polarized and grounded power source with voltage and frequency characteristics compatible with the specifications listed in Clause-2.
- 3) Apply one end of the equipotential cable to the equipotential column on the back of the AC adapter. Connect the other end of the equipotential cable to the equipotential busbar in the operating room or any other common equipotential connection point. Please see Clause 9 above.
- 4) The unit is now ready for use.
- 5) Please refer to the user manual of AC powered palaJet® AC Disposable Pulse Lavage System (Disposable Electronically Pulsed Lavage Suction Apparatus) (#: W-205).

13. Maintenance, Cleaning and Disinfection

- 1) Always clean and disinfect the exterior of the unit after every case, per hospital protocol. The unit must be cleaned (wiped) with a cloth that has been dampened (not dripping) with a mild disinfectant detergent.
- 2) This device does not contain user repairable parts and does not require repair.

14. Storage Environment

The following are environmental conditions for storage:
Ambient temperature range: -4 °F to 104°F (-20 °C to +40 °C).
Relative humidity range: 10% RH to 90% RH (Non-Condensing).

15. Labels

palaJet® AC Power Adapter
AC Power Supply for Disposable Electronically Pulsed Lavage Suction Apparatus

W-205 REF 5203576 (European Plug)
 5206240 (British Plug)

CE 0123 **EMC**  

Guangzhou Clean Medical Products Manufacturing Corp.
 No.163 Meidu Road, Chengjiao, Conghua District, 510900
 Guangzhou, Guangdong, China. Tel: +86-20-87654060, 87817888

POWER INPUT: 120VA OUTPUT: DC 7.5-12V
100-240V ~, 50/60Hz POWER: 48W

DANGER:
EXPLOSION HAZARD, DO NOT USE IN THE PRESENCE OF FLAMMABLE ANESTHETICS OR GASES.

 **CAUTION:**
RISK OF ELECTRIC SHOCK. DO NOT REMOVE COVER. REFER SERVICING TO QUALIFIED SERVICE PERSONNEL.

MEDICAL EQUIPMENT
WITH RESPECT TO ELECTRIC SHOCK, FIRE, AND MECHANICAL HAZARDS ONLY, IN ACCORDANCE WITH IEC 60601-1.

  **yyyy-mm**

LOT Lyyymmdd-x **SN** yymmdd-x.xxx

 **CAUTION**
 Use with **palaJet® AC**
 Pulse Lavage Handpiece only.

 **CAUTION**
 Do not sterilize!
 Wipe with disinfectant only!

 **CAUTION**
 GROUND RELIABILITY CAN ONLY BE ACHIEVED WHEN POWER SUPPLY CORD IS CONNECTED TO A RECEPTACLE MARKED "HOSPITAL GRADE"

palaJet® AC Power Adapter
 AC Power Supply for Disposable Electronically Pulsed Lavage Suction Apparatus

W-205

REF 5203576 (European Plug)
 5206240 (British Plug)

Guangzhou Clean Medical Products Manufacturing Corp.
 No.163 Meidu Road, Chengjiao, Conghua District,
 510900 Guangzhou, Guangdong, China.
 Tel: +86-20-87654060, 87817888

POWER INPUT: 120VA
 100-240V~, 50/60 Hz
 OUTPUT: DC 7.5-12V
 POWER: 48W

DANGER:
 EXPLOSION HAZARD, DO NOT USE IN THE PRESENCE OF FLAMMABLE ANESTHETICS OR GASES.

CAUTION:
 RISK OF ELECTRIC SHOCK, DO NOT REMOVE COVER. REFER SERVICING TO QUALIFIED SERVICE PERSONNEL.

MEDICAL EQUIPMENT
 WITH RESPECT TO ELECTRIC SHOCK, FIRE, AND MECHANICAL HAZARDS ONLY, IN ACCORDANCE WITH IEC 60601-1.

Barcode:  Lot number
 Date of manufacture
 Serial number

CAUTION
 Use with **palaJet® AC**
 Pulse Lavage Handpiece only.

CAUTION
 Do not sterilize!
 Wipe with disinfectant only!


CAUTION
 GROUND RELIABILITY CAN ONLY BE ACHIEVED WHEN POWER SUPPLY CORD IS CONNECTED TO A RECEPTACLE MARKED "HOSPITAL GRADE"

16. Symbol Descriptions


 Manufacturer	 Fragile, handle with care
 Date of manufacture	 Keep dry
 Catalogue number	 This way up
 Batch code (Lot number)	 Keep away from sunlight
 Serial number	 Temperature limit
 Use-by date	 Electric hazard
 Medical device	 Equipotentiality
 Consult instructions for use	 Direct current
 Non-sterile	 Alternating current
 CE marking of conformity	 Caution
 Declaration that a product meets the EMC regulations per marking requirements.	 Authorized Representative in the European Community/European Union.
 Model number	 Importer
 Unique device identifier	 Distributor
QTY QTY Quantity	 Type BF applied part
 UK Responsible Person	
 Disposal of waste electrical and electronic equipment. This product contains electrical or electronic materials. The presence of these materials may, if not disposed properly, have potential adverse effects on the environment. Presence of this label on the product means it must not be disposed of in normal household waste and must be disposed of separately. Please follow local/national current regulations to properly recycle or dispose the failed and expired AC power supply.	

17. Cable Plug

Europe/
Plug Type E & F

United Kingdom & Ireland/
Plug Type G

18. Manufacturer

Name: **Guangzhou Clean Medical Products Manufacturing Corp.**

Address: No.163 Meidu Road, Chengjiao, Conghua District, Guangzhou, Guangdong, China.

P.C.: 510900

Tel: +86-20-87654060, 87817888

E-mail: info@cleanmfg.com

Web: www.cleanmfg.com

19. EU Authorized Representative

Name: **SUNGO Europe B.V.**

Address: Fascinatio Boulevard 522, Unit 1.7, 2909VA Capelle aan den IJssel, The Netherlands.

Tel: +31(0)10 3034500

E-mail: ec.rep@sungogroup.com

20. Importer and Distributor in EU

Name: **Heraeus Medical GmbH**

Address: Philipp-Reis-Straße 8/13, 61273 Wehrheim, Germany.

21. Importer and Distributor in UK

Name: **Heraeus Medical UK Ltd.**

Address: Oxford House, 12-20 Oxford Street, Newbury, Berkshire, RG14 1JB, England, United Kingdom.

22. UK Responsible Person

Name: **SUNGO Certification Company Limited**

Address: 3rd Floor, 70 Gracechurch Street, London, EC3V 0HR, England, United Kingdom.

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⚡ CAUTION:
RISK OF ELECTRIC SHOCK. DO NOT REMOVE COVER. REFER SERVICING TO QUALIFIED SERVICE PERSONNEL.

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  yyyy-mm

LOT Lyyymmdd-x SN yymmdd-x.xxx

⚠ CAUTION
 Use with **palaJet® AC**
 Pulse Lavage Handpiece only.

⚠ CAUTION
 Do not sterilize!
 Wipe with disinfectant only!

⚠ CAUTION
 GROUND RELIABILITY CAN ONLY BE ACHIEVED WHEN POWER SUPPLY CORD IS CONNECTED TO A RECEPTACLE MARKED "HOSPITAL GRADE"

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Guangzhou Clean Medical Products Manufacturing Corp.
 No.163 Meidu Road, Chengjiao, Conghua District,
 510900 Guangzhou, Guangdong, China.
 Tel.: +86-20-87654060, 87817888

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W-205 REF 5203576 (European Plug)
5206240 (British Plug)

CE 0123 EMC

Guangzhou Clean Medical Products Manufacturing Corp.
No.163 Meidu Road, Chengjiao, Conghua District, 510900
Guangzhou, Guangdong, China. Tel: +86-20-87654060, 87817888

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CAUTION:
RISK OF ELECTRIC SHOCK. DO NOT REMOVE COVER. REFER SERVICING TO QUALIFIED SERVICE PERSONNEL.

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yyyy-mm

LOT Lyyymmmdd-x SN yymmmdd-x.xxx

CAUTION
Use with **palaJet® AC**
Pulse Lavage Handpiece only.

CAUTION
Do not sterilize!
Wipe with disinfectant only!

CAUTION
GROUND RELIABILITY CAN ONLY BE ACHIEVED WHEN POWER SUPPLY CORD IS CONNECTED TO A RECEPTACLE MARKED "HOSPITAL GRADE"

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Guangzhou Clean Medical Products Manufacturing Corp.
No.163 Meidu Road, Chengjiao, Conghua District,
510900 Guangzhou, Guangdong, Kiina.
Puh.: +86-20-87654060, 87817888

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

18. Valmistaja

Nimi: **Guangzhou Clean Medical Products Manufacturing Corp.**

Osoite: No.163 Meidu Road, Chengjiao, Conghua District, Guangzhou, Guangdong, Kiina

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Nimi: **SUNGO Certification Company Limited**

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Namn: **Guangzhou Clean Medical Products Manufacturing Corp.**
Namn: **SUNGO Europe B.V.**
Namn: **Heraeus Medical GmbH**
Namn: **Heraeus Medical UK Ltd.**
Namn: **SUNGO Certification Company Limited**