
Heraeus

COPAL[®] knee moulds

Moulds for temporary knee spacer

COPAL® knee moulds

EN	Instructions for use (Revision status: 2020-02)	4
BG	Инструкции за употреба (Дата на актуализиране на текста: 2020-02)	11
CS	Návod k použití (Poslední revise: 2020-02)	18
DA	Brugsanvisning (Aktualiseret: 2020-02)	25
DE	Gebrauchsanweisung (Revisionsstand: 2020-02)	32
EL	Οδηγίες χρήσης (Ημερομηνία αναθεώρησης: 2020-02)	39
ES	Instrucciones de uso (Última revisión: 2020-02)	46
ET	Kasutusjuhend (Viimati muudetud: 2020-02)	53
FI	Käyttöohje (Revisiostatus: 2020-02)	60
FR	Notice d'utilisation (Date de révision: 2020-02)	67
HR	Upute za uporabu (Stanje revizije: 2020-02)	74
HU	Alkalmazási utasítás (Aktualizálás időpontja: 2020-02)	81
IT	Istruzioni per l'uso (Ultima revisione: 2020-02)	88
KO	사용 설명서 (2020년 02월 수정)	95
LT	Naudojimo instrukcija (Tikrinta: 2020-02)	102
NL	Gebruiksaanwijzing (Revisiestatus: 2020-02)	109
NO	Bruksinformasjon (Revisjonsstatus: 2020-02)	116

PL	Instrukcja użytkowania (Ostatnia aktualizacja: 2020-02)	123
PT	Instruções de utilização (Data da última actualização: 2020-02)	130
RO	Instrucţiuni de utilizare (Stadiu revizie: 2020-02)	137
RU	Инструкция по применению (Ревизионный статус: 2020-02)	144
SK	Návod na použitie (Posledná revízia textu: 2020-02)	151
SL	Navodilo za uporabo (Datum revizije: 2020-02)	158
SV	Bruksanvisning (Revision status: 2020-02)	165
TR	Kullanma Talimatı (Bilgilerin son güncellenme tarihi: 2020-02)	172
ZH	使用说明书 (修订状态: 2020 年 02 月)	第179页

COPAL® knee moulds

Product description

COPAL® knee moulds are sterile single-use moulds used for the preparation of spacers that are intended as temporary knee replacements as part of two-stage septic joint endoprosthesis revision.

COPAL® knee moulds comprise a tibial component and a femoral component, which together form a bearing and move against one another. They can be used in both the right and the left knee joint.

The spacer function provides that after removal of the endoprosthesis the existing joint space is retained and contraction of the musculature and the surrounding tissues is prevented.


COPAL® knee moulds are intended for single use and must not be re-used or resterilised.

Package contents

COPAL® knee moulds are available in 3 different sizes.

Size	Femur (medial-lateral)	Tibia (medial-lateral)
S	60 mm	65 mm
M	70 mm	75 mm
L	80 mm	85 mm

The height of the tibia can be adjusted in all 3 sizes from 12 mm to 37 mm in 5-mm increments.

One set of **COPAL® knee moulds** includes:

- 1 mould to prepare the femoral component (a)
- 1 mould to prepare the tibial component (b)
- 1 optional shaft for the tibial component, included in (b)
- 1 plug to force any remaining cement out of the filling opening on the femoral component (c)
- 1 trial spacer of 12 mm (d)
- 2 extensions for the trial spacer of 10 mm (e)
- 1 extension for the trial spacer of 5 mm (f)

COPAL® knee moulds are double packaged. The outer blister is non-sterile on the outside and sterile on the inside. Inside the outer blister is a sterile inner blister that contains the set. Both blisters are sealed with Tyvek.

Composition

All components of the **COPAL® knee moulds** are made of medical-grade plastic and are latex-free and PVC-free.

Indication

COPAL® knee moulds are designed to prepare spacers by filling the moulds with bone cement. The spacers prepared with **COPAL® knee moulds** are indicated for temporary knee replacement as part of two-stage joint endoprosthesis revision in cases of infection and inflammatory reactions following explantation of an existing joint endoprosthesis. They are used in the first step following the explantation of the existing endoprosthesis and are removed prior to the implantation of a new endoprosthesis in a second step.

Contraindications

The spacers prepared with **COPAL® knee moulds** must not be implanted in patients in whom:

- their health precludes two-stage endoprosthesis revision.
- the loss of bone tissue from the proximal tibia and/or the distal femur is so extensive that anchoring the spacer with bone cement is not possible.
- the bone quality has been so negatively affected, by osteoporosis, for example, that loosening or migration of the spacer or an inadequate mechanical weight-bearing capacity of the bone bed can be expected.
- the musculature and ligaments do not have sufficient functional competence or there are neuromuscular or vascular disorders present.
- myasthenia gravis has been diagnosed.
- osteosynthesis material may interfere with the bearing of the tibial and femoral components.
- an additional systemic and/or chronic infection is present.
- the infected knee joint endoprosthesis cannot be removed either fully or in part.
- there is no total knee replacement present and the infection is secondary to trauma, septic arthritis or other diseases that require surgical treatment.
- the skeleton is not yet mature.

Side effects

The following adverse effects are possible when using spacers prepared with **COPAL® knee moulds**:

- Damage to the femur and tibia
- Damage to the surrounding blood vessels and nerves
- Dislocation of the joint or the spacer
- Changes to the surrounding tissue
- Difference in leg lengths

Warnings and precautions

User

Before using **COPAL® knee moulds** the user should be well acquainted with their properties, handling and application. We recommend that users practice the complete procedure before using it for the first time. It is also essential that users are familiar with the bone cement application system used to fill the **COPAL® knee moulds**. Two mixing and application systems for bone cement are required for the

filling process. The femoral component can only be filled using an application system with an application nozzle having a maximum diameter of 15 mm. We recommend filling the **COPAL® knee moulds** with a **PALAMIX®** mixing system. The stability of the spacer is reduced if there are air pockets and cracks in the cement; these can be prevented by slowly filling the moulds and using vacuum mixing for the bone cement. We do not recommend the use of fast-curing bone cements. A bone cement with a processing time of at least 3.5 minutes under surgical temperature conditions should be used.

COPAL® knee moulds may only be used if they are undamaged. The spacers prepared with **COPAL® knee moulds** must not be positioned by striking with a hammer or similar instruments.

The user must ensure that the bone cement used for fixation completely covers the surfaces of the tibial and femoral components intended for fixation. In doing so, prevent the bone cement reaching the bearing surfaces of the tibial and femoral components. After using **COPAL® knee moulds** the presence of all parts noted on the counting chart should be checked.

After surgery

After surgery the patient may only stand up, walk and sit down with the help of walking aids with partial weight bearing and maximum bending of the knee of 90°.

The patient must be informed of this beforehand and warned that there is a high risk of the spacer breaking with full weight-bearing or high levels of activity. If the function of the musculature and the ligaments is limited, the use of an orthosis to stabilise the joint function may be helpful. The decision as to whether and to what extent partial weight bearing can occur is made by the surgeon and should be determined by the anatomical situation, the health status of the patient and the stability of the spacer after insertion.

Handling

Preparation

Prior to the procedure, an X-ray is recommended to determine the necessary size of the **COPAL® knee moulds** set using the X-ray templates that are available from Heraeus Medical GmbH. The sizes of the **COPAL® knee moulds** are also available in most digital surgery planning systems.

The user must adjust pre-operatively to the anticipated situation (e.g. defect type and expected stability of the bone bed) and have appropriate instruments at hand. The **COPAL® knee moulds** set does not include any instruments. The user must also be familiar with the bone cement and the mixing and application system that will be used to fill the **COPAL® knee moulds**. The instructions for use for all the products used must be noted beforehand.

Before using the **COPAL® knee moulds** set, it should be checked for any damage and completeness. The blister must only be opened under sterile conditions. For this purpose, the inner blister is available sterile after the removal of the outer blister Tyvek. If the sterile barrier is damaged, the **COPAL® knee moulds** set must not be used and must be disposed of.

The femoral component is made up of two pre-assembled halves that are held together by sliding closures. There is a filling port for the bone cement in the centre of the back of the mould. Integrated into the port is an eject mechanism that enables the prepared spacer to be released from the mould after the bone cement has cured by screwing it out of the mould in a clockwise direction.

The tibial component is made up of a height-adjustable mould in which the height of the tibial spacer can be set in 5-mm increments (from 12 mm to 37 mm). The trial spacer is used to determine the appropriate height of the tibia. The tibial component can optionally be fitted with a shaft.

The following quantities of bone cement are needed to fill the **COPAL® knee moulds**: (Measured using the **PALAMIX®** vacuum mixing system and Heraeus bone cement – please note that the volume of other bone cements may deviate from these figures and there may be different quantities of residual cement in other mixing systems.)

Filling volume for the femoral component

Size	Recommended number of packages of bone cement (1x40)*
S	1
M	2
L	2

Filling volume for the tibial component

Size	Height	Recommended number of packages of bone cement (1x40)*
S	Up to 32 mm	2
S	From 37 mm	3
M	Up to 27 mm	2
M	From 32 mm	3
L	Up to 22 mm	2
L	From 27 mm	3

* 1x40 **PALACOS®/COPAL®** bone cement yields ~58 g bone cement. Of this, ~10 g remains in the **PALAMIX®** vacuum mixing system. For the cement quantity produced, please refer to the instructions for use of the bone cement manufacturer.

Fixation with bone cement is used primarily to stabilise the spacer. Deep penetration of the bone cement into the bone structure is not desirable and can make subsequent explantation of the spacer more difficult. Excess bone cement must be removed.

After the bone cement used for hardening is cured, the joint space must be cleaned of any bone cement particles. During the cleaning, the spacers should not come into contact with aqueous solution nor should the pulse lavage touch the spacers because any substances present in the bone cement may be flushed out as a result.

Explantation

The spacers prepared with **COPAL® knee moulds** may remain in the patient for a maximum of 180 days and must then be explanted. After the explantation, the knee joint must be provided with an appropriate knee joint endoprosthesis or another suitable treatment (e.g. arthrodesis of the knee joint) must be carried out.

Storage

COPAL® knee moulds must be stored unopened in a dry and clean location in the original packaging.

Shelf-life

The expiry date is shown on the outer and inner blisters.

COPAL® knee moulds must not be used after the expiry date.


Sterility

COPAL® knee moulds are sterile. Ethylene oxide gas is used for sterilisation. They are intended for single use and must not be re-used or re-sterilised because successfully preparing spacers cannot be guaranteed if there have been changes to the material. The contents of an opened set or damaged packages must be disposed of.


Disposal

COPAL® knee moulds and the outer packaging can be disposed of in normal hospital waste. Please comply with local environmental legislation and waste disposal regulations.


Preparation of the spacer Femoral component


Check the contents of the **COPAL® knee moulds** for completeness. Prepare the mixing system and the bone cement.


When the cement reaches the application phase*, slowly and evenly fill the mould. The nozzle must completely touch the bottom throughout the entire filling process. The mould is completely filled when small amounts of cement exit the vents.


Remove the application nozzle and attach the plug ('click!'). The plug has a guide and two recesses to lock it in place.

When the cement has completely cured*, break off the sliding closures.


Compress the mould between your hands to release the spacer from the surfaces. A cracking sound indicates that it has detached.


1. Remove the top half.
2. Screw the spacer off in a clockwise direction.

If necessary, mechanically detach the top mould using a suitable tool and repeatedly screwing the screw in and out.


* For the exact processing times, refer to the instructions for use of the bone cement manufacturer.

Break the base plate off the spacer.


If necessary, remove any burrs or unevenness.


Insert the spacer without fixing it. Determine the distance between the femoral component and the tibia using the modular trial spacer. The height of the bone cement used for fixation must be noted here.

Tibia component

Prepare the mixing system and the bone cement. Adjust the distance for the height of the tibial spacer by pushing up the base plate according to the scale.


If necessary, correct the height: Push the base plate all the way up and remove it. Insert the base plate again with the bearing surface facing upwards.


Attach the lid and all 4 sliding closures and close them.


If required, attach the shaft. ('click'!)


Mix the bone cement in the mixing system and attach the shortest possible application nozzle. As soon as the cement reaches the application phase*, slowly and evenly fill the mould.


Remove any excess bone cement.


Close the shaft with the lid. If the shaft is not needed, the slot provided is covered with the large lid located on the shaft.


When the cement has completely cured*, break off the sliding closures.

Compress the mould between your hands to release the spacer from the surfaces. A cracking sound indicates that it has detached.

Mechanically release the cover plate and remove it. Press the spacer out from below.


If necessary, remove any burrs or unevenness.


After using **COPAL® knee moulds** the presence of all parts noted on the counting chart should be checked.

* For the exact processing times, refer to the instructions for use of the bone cement manufacturer.


Use in patients
Implantation


1. Completely remove the infected endoprosthesis and residual cement.


2. Thoroughly debride the joint space.


3. Apply bone cement onto the distal, non-sliding half of the tibial spacer and the proximal bone surface of the tibia.


4. Position the tibial spacer onto the proximal surface of the tibia and carefully remove any excess bone cement residue.

5. Apply bone cement onto the non-sliding side of the femoral spacer.

6. Position the femoral spacer onto the distal femur and hold until the bone cement is almost cured.

7. Before the bone cement has completely cured, check that both spacers are correctly seated by repeatedly bending the knee joint.

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