
User Manual

palaJet® AC – Disposable Pulse Lavage System

Disposable Electronically Pulsed Lavage Suction Apparatus


palaJet® AC – Disposable Pulse Lavage System

EN	User Manual	3
CS	Uživatelská příručka	13
DE	Gebrauchsanweisung	23
EL	Εγχειρίδιο χρήσης	33
ES	Manual de usuario	43
FI	Käyttöopas	53
FR	Manuel de l'utilisateur	63
IT	Manuale d'uso	73
NL	Gebruikershandleiding	83
NO	Bruksanvisning	93
PL	Instrukcja użytkowania	103
SK	Príručka pre používateľa	113
SV	Användarhandbok	123

User Manual

palaJet® AC – Disposable Pulse Lavage System


Disposable Electronically Pulsed Lavage Suction Apparatus

Guangzhou Clean Medical Products Manufacturing Corp.

AF-205-09E(Z), Rev.: B

palaJet® AC – Disposable Pulse Lavage System

1. Medical Device

Name: palaJet® AC – Disposable Pulse Lavage System
(Disposable Electronically Pulsed Lavage Suction Apparatus)

W-205

REF 5203575

2. Indications (Intended Use)


- 1) This device is used together with normal sterile saline or other applicable solutions and vacuum extractor. Vacuum extractor is not a must.
- 2) This device can be used for open wound cleaning, soft tissue cleaning, and surgical operation location cleaning.
- 3) This device is suitable to be used in operation room, emergency room and treatment room.

3. Components

- 1) This device mainly consists of AC Power supply (B) (Including AC Power wire (B.1)) and Pulsed lavage kit (A), Pulsed lavage kit (A) is the applied part.

- 2) The Pulsed lavage kit (A) include Handpiece (1) (Including DC motor, Compression pump, Trigger (1.2), Nozzle lock (1.1), Enclosure), Nozzle (2) (Such as Fan Spray Nozzle (Short Nozzle) (2a), Femoral Nozzle (Long Nozzle) (2b), Femoral Brush Nozzle (Long Nozzle) (2c), Shower Spray Nozzle (Short Nozzle) (2d), Large Splash Shield (2e) or other nozzles. See the package and labels for details.), Waste pipe (3) (Including Clamp (3.1), Connecting pipe (3.2)), Irrigation pipe (4) (Including Clamp (3.1), Luer fitting (4.1)), Power wire (5).
- 3) Associated devices: Vacuum extractor (7) and Irrigation bag (6) are not included in this device, but prepared by medical institution. The negative pressure of Vacuum extractor (7) within the range of 5 kPa to 40 kPa shall be adjustable. The Irrigation bag (6) should be used with normal sterile saline or other applicable solutions (not suitable for viscous solutions).
- 4) Accessories: Nozzles (2).
- 5) Detachable parts: Nozzle (2).
- 6) Applied parts: Nozzle (2).

Device structure diagram:

4. Working Principle

- 1) This device is powered by AC power supply (AC power adapter).
- 2) This device consists of DC motor, compression pump, handpiece, enclosure, nozzles and AC power supply.
- 3) When squeezing the trigger, power is on. When power is turned on, the DC motor drives the Compression pump in reciprocating motion, engender cavity volume changes, thus it can breathe in and out the solution. There are a low power position, a high power position and a stop tap position. The waste pipe integrated to the handpiece can be connected to a vacuum extractor to extract waste. Different nozzles can be fitted on during operation.


1 Install nozzle
2 Connect vacuum extractor
3 Spike irrigation bag
4 Hang on an IV pole
5 Plug AC power supply to hospital socket
6 Plug output terminal lead
7 Squeeze trigger start use

1. Unlock (Lock button: push up = unlock, push down = lock)

2. Insert Nozzle

Push down to lock nozzle

- 4) Connect the **Waste pipe (3)** to the **Vacuum extractor (7)** (if prepared). Vacuum extractor is not a must.
- 5) Connect the **Luer fitting (4.1)** on the **Irrigation pipe (4)** to the **Irrigation bag (6)**.


- 6) Hang the **AC Power supply (B)** on an IV pole using the unit's built in hanger. To avoid the liquid from the **Irrigation bag (6)** spillage on the **AC Power supply (B)**.

Note:

The **AC Power supply (B)** always be hung on an IV pole and should never be allowed to rest on a table or floor. The table and/or floor may become wet which could damage the unit.


- 7) Plug the **Power wire (B.1)** of **AC power supply (B)** to the socket.
- 8) Plug the **Power wire (5)** of **Pulsed lavage kit (A)** to the **AC Power supply (B)** output terminal lead

Plug the power wire of Pulsed lavage kit to the AC Power supply output terminal lead.

Notice direction!

- 9) Use the **Nozzle (2)** to aim the treatment area. Start use and adjust working condition by squeezing **Trigger (1.2)**. The **Trigger (1.2)** have a low power position with symbol short "—", a high power position with symbol long "—", and a stop tap position with symbol "O".

Trigger control


1: Stop

2: High pulse

3: Low Pulse

- 10) Slide the **Clamps (3.1)** can control the liquid flow of both the **Irrigation pipe (4)** and **Waste pipe (3)**. Please lock the **Clamps (3.1)**, when you stop operating the device.

The Clamps can control the liquid flow of both the Irrigation pipe and Waste pipe

- 11) It can suck the waste solution concurrently while spraying or after, if the **Vacuum extractor (7)** is prepared.
- 12) To press **Nozzle lock (1.1)** upward first, when withdrawing the **Nozzle (2)**.

- 13) Optional: The **Large Splash Shield (2e)** is to be used in conjunction with the **Short nozzle (2a, 2d)**. Its function is to prevent the expansion of the splashing range of the washing liquid, and it can be used or not used according to the application part.
- To install: Pinch the soft head of the **Short nozzle**, insert into the opening of the **Large Splash Shield** and slightly pull in opposite directions to fixate the **Large Splash Shield** on the nozzle head.
- To uninstall: Pinch the soft head of the **Short nozzle** and remove the **Large Splash Shield**.

6. Specification

No.	Item	Specification
1)	Power supply (Ratings)	100–240V ~ (AC), 50/60Hz, 120VA Line Current: 0.50A RMS @ 120 ~ (AC) (Max.)
2)	Working voltage	DC 7.5–12V
3)	Power rating	48W
4)	Handpiece dimension	About 185mm x 55mm x 40mm (Long x Wide x High)
5)	Net weight	About 0.6 kg
6)	Flow rate (high power position)	Not less than 0.9L/minute.
7)	Maximum irrigation length (high power position)	Not less than 2.5meters.
8)	Noise (Decibel)	Less than 75 dB(A)
9)	Solution form (The shape of the spray solution)	Fog-like and thread-like.
10)	Sterilization	This device is sterilized using ethylene oxide before delivery.
11)	Valid period (Shelf life)	Two years.

7. IEC Classification

- 1) Class I equipment.
- 2) Type BF applied part.
- 3) IPX0.
- 4) Not category AP/APG equipment.

8. Contraindications

Not found.

9. Adverse Effect

Not found.

10. Warning and Precautions

- 1) This device should be operated only by qualified professionals.
- 2) This device (**Pulsed lavage kit**) should not be used after the expiry date shown on the package.
- 3) This device (**Pulsed lavage kit**) is for single patient use only. Do not reuse the device.
- 4) Prior to each use, inspect the package. Do not use the device if the package is damaged.
- 5) This device is not evaluated as category AP or APG equipment. This device is not suitable for use in the presence of flammable anesthetic mixtures with air or oxygen or nitrous oxide.
- 6) No modification of this equipment is allowed.

- 7) Except the Nozzle, there is no replaceable part and accessory within the unit.
- 8) Follow current local/national regulations to properly dispose the device after use.

11. Disposal of the Power Supply

- 1) The AC power supply is a non-disposable part. Clean in accordance to the hospital's regulations before each use.
- 2) Follow local/national current regulations to properly recycle or dispose the power supply after expiration date.
- 3) Do NOT cut the power wire at any time before removing the AC power supply to prevent short-circuit, which may cause overheating of the AC power supply and result in a potential fire risk.

12. Maintenance, Cleaning and Sterilization (**Pulsed lavage kit**)

- 1) This device is maintenance free and does not require routine maintenance.
- 2) This device does not contain user repairable parts and does not require repair.
- 3) This device is for single patient use only. There is no cleaning or sterilization requirement for the user. Do not resterilize.

13. Symbol Descriptions

The following symbols will appear on user manual, labels and packages. Some of the symbols represent standards and compliances associated with the device and its use.

 Sterilized using ethylene oxide	 Do not re-use	
 Consult instructions for use	 Do not re-sterilize	
 Do not use if package is damaged and consult instructions for use	 Caution	
 CE marking of conformity	 DISPOSAL: Do not dispose this product as unsorted municipal waste. Collection of such waste separately for special treatment is necessary.	
 Authorized representative in the European Community/European Union		
 Batch code (Lot number)	 Catalogue number	 Model number
 Date of manufacture	 Serial number	
 This symbol on the trigger side means stop tap position.	 Use-by date	
 This symbol on the trigger side means low power position.	 Double sterile barrier system	
 This symbol on the trigger side means high power position.	 Temperature limit	
 Unique device identifier	 Medical device	 Manufacturer
 This way up	 Fragile, handle with care	 Importer
 Keep dry	 Keep away from sunlight	 Distributor
 Alternating current	 Electric hazard	 Type BF applied part
 Quantity	 Direct current	 UK Responsible Person

14. Electromagnetic Compatibility Descriptions

Conform to the requirements of IEC 60601-1-2:2014 and EN 60601-1-2:2015.

Guidance and manufacture's declaration – electromagnetic emissions – for all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacture's declaration – electromagnetic emission		
The W-205 are intended for use in the electromagnetic environment specified below. The customer of the user of the W-205 should assure that they are used in such an environment.		
Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 14-1	Group 1, Class A	The W-205 are not suitable for interconnection with other equipment.
Harmonic emissions IEC 61000-3-2	Comply	The product is suitable for use in all facilities, including household facilities and directly connected to the residential public low-voltage power supply network.
Voltage fluctuations/flicker emissions IEC 61000-3-3	Comply	

Guidance and manufacture's declaration – electromagnetic immunity – for all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacture's declaration – electromagnetic immunity			
The W-205 are intended for use in the electromagnetic environment specified below. The customer or the user of W-205 should assure that they are used in such an environment.			
Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	Contact: ± 8 kV Air: ± 2 kV, ± 4 kV, ± 6 kV, ± 8 kV, ± 15 kV	Contact: ± 8 kV Air: ± 15 kV	The floor shall be of wood, concrete, or tile, and if the ground is covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	Power supply lines: ± 2 kV input/output lines: ± 1 kV 100 kHz repetition frequency	Power supply lines: ± 2 kV	The power supply should have a typical commercial or hospital environment.
Surge IEC 61000-4-5	line(s) to line(s): ± 0.5 kV, ± 1 kV. line(s) to ground: ± 0.5 kV, ± 1 kV, ± 2 kV.	Line to line: ± 1 kV	The power supply should have a typical commercial or hospital environment.
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	0% U_i , 0.5 cycle At 0° , 45° , 90° , 135° , 180° , 225° , 270° and 315° 0% U_i , 1 cycle And 70% UT 25/30 cycles Single phase: at 0 0% U_i , 250 cycle	0% for 0.5 cycle 0% for 1 cycle 70% for 25 cycles 0% for 250 cycles	The power supply should have a typical commercial or hospital environment.
Power frequency magnetic field IEC 61000-4-8	50 Hz: 30 A/m 60 Hz: 30 A/m	50 Hz: 30 A/m 60 Hz: 30 A/m	The power frequency magnetic field should have the characteristics of power frequency magnetic field in typical commercial or hospital environment.
NOTE: UT is the a.c. mains voltage prior to application of the test level.			

**Guidance and manufacture's declaration – electromagnetic immunity –
for ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING**

Guidance and manufacture's declaration – electromagnetic immunity			
The W-205 are intended for use in the electromagnetic environment specified below. The customer or the user of W-205 should assure that they are used in such an environment.			
Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	150 KHz to 80 MHz 3 Vrms; ISM bands between 150 KHz to 80 MHz 6 Vrms	3 Vrms 6 Vrms (in ISM bands) 80% Am at 1 kHz	Portable and mobile RF communications equipment should be used no closer to any part of the W-205 including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance
Radiated RF IEC 61000-4-3	80 MHz to 2.700 MHz 3 V/m (rms) 385 MHz 27 V/m (rms) 450 MHz 28 V/m (rms) 710 MHz, 745 MHz, 780 MHz 9 V/m (rms) 810 MHz, 870 MHz, 930 MHz 28 V/m (rms) 1,720 MHz, 1,845 MHz, 1,970 MHz 28 V/m (rms) 2,450 MHz 28 V/m (rms) 5,240 MHz, 5,500 MHz, 5,785 MHz 9 V/m (rms)	3 V/m, 80% Am at 1 kHz 27 V/m PM at 18 Hz 28 V/m FM $\pm$ 5 kHz deviation at 1 kHz sine 9 V/m PM at 217 Hz 28 V/m PM at 18 Hz 28 V/m PM at 217 Hz 9 V/m PM at 217 Hz	$d = [3.5/V_1] \sqrt{P}$ $d = [3.5/E_1] \sqrt{P}$ 80 MHz to 800 MHz $d = [7/E_1] \sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the W-205 are used exceeds the applicable RF compliance level above, the W-205 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the W-205 .			
b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

15. Environmental Requirements

1) Storage conditions:

This device should be stored in non-corrosive gases, cool, dry, and ventilated place.
Temperature: -10 °C ~ +40 °C;
Humidity: 20 % RH ~ 90 % RH;
Air pressure: 50 kPa ~ 106 kPa.

2) Operating conditions:

Temperature: 5 °C ~ 40 °C;
Humidity: ≤ 80 % RH;
Air pressure: 86 kPa ~ 106 kPa.

16. Manufacturer

Name: **Guangzhou Clean Medical Products**

Manufacturing Corp.

Address: No.163 Meidu Road, Chengjiao, Conghua
District, Guangzhou, Guangdong, China.
P.C.: 510900
Tel: +86-20-87654060, 87817888
E-mail: info@cleanmfg.com
Web: www.cleanmfg.com

17. EU Authorized Representative

Name: **SUNGO Europe B.V.**

Address: Fascinatio Boulevard 522, Unit 1.7, 2909VA
Capelle aan den IJssel, The Netherlands.
Tel: +31(0)10 3034500
E-mail: ec.rep@sungogroup.com

18. Importer and Distributor in EU

Name: **Heraeus Medical GmbH**

Address: Philipp-Reis-Straße 8/13, 61273 Wehrheim,
Germany.

19. Importer and Distributor in UK

Name: **Heraeus Medical UK Ltd.**

Address: Oxford House, 12-20 Oxford Street, Newbury,
Berkshire, RG14 1JB, England, United Kingdom.

20. UK Responsible Person

Name: **SUNGO Certification Company Limited**

Address: 3rd Floor, 70 Gracechurch Street, London,
EC3V 0HR, England, United Kingdom.
Tel: +44-(0) 7868 234015
E-mail: ec.rep@sungoglobal.com

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Namn: **Guangzhou Clean Medical Products Manufacturing Corp.**
Namn: **SUNGO Europe B.V.**
Namn: **Heraeus Medical GmbH**
Namn: **Heraeus Medical UK Ltd.**
Namn: **SUNGO Certification Company Limited**