

palajet®

AVAILABLE ATTACHMENTS

Fan Spray Nozzle/Knee Nozzle*

- Attachment for knee surgery and other wound cleaning
- Dimensions: 97 mm x 32.5 mm



Shower Spray Nozzle

- Nozzle with reduced pressure for spine surgery
- Dimensions: 97 mm x 32.5 mm




Large Splash Shield

- Additional splash shield for the knee and shower spray nozzles for intensive wound cleaning
- Dimensions: 120 mm x 58 mm

Femoral Nozzle/Hip Nozzle*

- Attachment for thorough cleaning of the femoral bone bed
- Dimensions: 252 mm x 8.4 mm



Femoral Brush Nozzle

- Attachment for hip revisions
- Dimensions: 252 mm x 8.4 mm



palaJet®

REMOVING & IMPROVING

All endoprosthesis interventions pose various challenges for the surgical team. When cleaning the bone bed, the individual anatomy of the joints and the level of contamination require special tools to achieve an optimal cleaning result and to ensure that the prosthesis has a long-term survival rate.

The palaJet® pulse lavage system offers different attachments for the specific requirements of different endoprosthesis surgeries, which can be easily attached to the palaJet® handpiece using a click mechanism, thus facilitating thorough cleaning.



Product	Description	Content	REF
palaJet® Disposable Pulse Lavage System	Pulse Lavage System for single use, Knee and Hip Nozzles included	10	5156831
palaJet® Femoral Brush Nozzle	Attachment for hip revisions (with brush head)	5	5167807
palaJet® Shower Spray Nozzle	Attachment for spinal surgery	5	5167808
palaJet® Large Splash Shield	Splash Shield for Knee Nozzle and Shower Spray Nozzles	5	5167809
palaJet® Fan Spray Nozzle/Knee Nozzle	Attachment for knee surgery and other wound cleaning	5	5202875
palaJet® Femoral Nozzle/Hip Nozzle	Attachment for thorough cleaning of the femoral bone bed	5	5202878

Simply order from Heraeus.

www.heraeus-medical.com



12626 EN