



**YOU CONTROL THE SOLUTION.
WE MAKE IT EASY.**

**EASY TO CUSTOMIZE.
EASY TO DELIVER.
OPTIMIZED FOR PERFORMANCE.**



EASY TO CUSTOMIZE

Synthecure® is a synthetic calcium sulfate that offers superior handling, rapid curing, and complete resorption within 30–60 days. This biocompatible material is indicated for direct placement at infection sites in bone voids and defects. Delivered as a complete kit, Synthecure® makes preparation simple and helps streamline the procedure.

SYNTHECURE® SYNTHETIC CALCIUM SULFATE

Designed to simplify bead preparation, support precise placement, and deliver consistent performance.



FEATURES:

- 100% fully synthetic calcium sulfate
- Implant-grade calcium sulfate hemihydrate powder
- Sterile saline mixing solution included
- Mixing bowl included in every kit
- Dedicated bead mats available in 3 mm, 4.8 mm, and 6 mm sizes
- Single-sided bead mats eliminate flipping and refilling



Everything you need for
efficient bead preparation in
one intuitive kit.

EASY TO DELIVER. OPTIMIZED FOR PERFORMANCE.

DESIGNED FOR PRECISION PLACEMENT

The Synthecure® system helps maintain control throughout delivery with thoughtful design features that simplify placement and reduce disruption in the OR.


EASY TO DELIVER:


- Integrated mixing bowl converts to a precision delivery funnel
- Targeted bead placement
- Reduced bead scatter
- Intact beads support controlled placement
- Collection tray helps minimize bead damage

OPTIMIZED FOR PERFORMANCE:

- Over 2x less acidic than other synthetic calcium sulfate products
- Biocompatible and biodegradable
- Excellent batch-to-batch consistency
- Intended to fill bony voids and gaps of the skeletal system

A proven synthetic calcium sulfate designed to support healing in infected and non-infected treatment sites.





AVAILABLE BEAD SIZES

3 mm | 4.8 mm | 6 mm



PRODUCT OFFERING

PASTE VOLUME	BEAD VOLUME	SETTING TIME	ORDER CODE
3cc	7cc	~4 min	50-003
5cc	12cc	~4 min	50-005
10cc	25cc	~4 min	50-010
20cc	50cc	~4 min	50-020

For more information, visit www.heraeus-medical-usa.com

Heraeus Medical, LLC
770 Township Line Rd
Yardley, PA 19067



14185
MKT-0007